# EasySep™ Red Blood Cell Lysis Buffer, 10X Concentrate

Red blood cell lysis buffer

**REF** #100-0749 20 mL

Document #10000017274 | Version 01

Scientists Helping Scientists™ | [WWW.STEMCELL.COM](http://WWW.STEMCELL.COM)

TOLL FREE PHONE 1 800 667 0322 • PHONE +1 604 877 0713

[INFO@STEMCELL.COM](mailto:INFO@STEMCELL.COM) • [TECHSUPPORT@STEMCELL.COM](mailto:TECHSUPPORT@STEMCELL.COM)

FOR GLOBAL CONTACT DETAILS VISIT OUR WEBSITE

## Intended Use

EasySep™ Red Blood Cell (RBC) Lysis Buffer, 10X Concentrate is intended to be used with EasySep™ IVD cell isolation kits for the lysis of red blood cells prior to the enrichment of targeted cells for downstream in vitro diagnostic assays.

## Product Description

EasySep™ RBC Lysis Buffer, 10X Concentrate is developed to allow osmotic lysis of red blood cells while providing a physiologically relevant pH optimized to maintain cell survival. It also has osmolality similar to the natural environment of the cells to maintain balanced osmotic pressure. It is typically used manually with EasySep™ cell isolation kits.

## Properties

**Storage:** Store at 15 - 25°C.

**Shelf Life:** EasySep™ RBC Lysis Buffer, 10X Concentrate is stable until the expiry date (use-by date) on the label.

NOTE: EasySep™ RBC Lysis Buffer (1X) can be stored at 2 - 8°C for up to 3 months. Do not exceed the expiry date (use-by date) of the original 10X buffer.

## Warnings and Precautions

1. FOR IN VITRO DIAGNOSTIC USE.
2. Refer to the Product Information Sheet (PIS) of the EasySep™ cell isolation kit for the list of materials required and detailed instructions on the use of this product with the kit for the enrichment of targeted cells.
3. This product should be handled by trained personnel observing good laboratory practices. Once this product is added to human cells, treat the suspension as potentially biohazardous. Handling and disposal of wastes should observe all local, state, or national regulations.
4. Insufficient lysis or excess lysis buffer can interfere with the EasySep™ cell isolation kit's performance.
5. Ensure that bottle is capped properly during storage. Proper handling and storage of this product is essential. Shelf life and performance may be adversely affected when stored outside of the recommended storage conditions.
6. Do not use the product if packaging container's integrity has been compromised upon arrival or if there are any signs of contamination, cloudy appearance, discoloration, drying, leaks or any other evidence of deterioration. Unused product may be disposed of according to standard laboratory procedures for non-hazardous liquids.
7. Do not use the product beyond the use-by date indicated on its label.
8. Do not refrigerate or freeze EasySep™ RBC Lysis Buffer, 10X Concentrate.
9. Important information regarding the safe handling, transport, and disposal of this product is contained in the Safety Data Sheet (SDS).

NOTE: SDS of this product is available at [www.stemcell.com](http://www.stemcell.com).

STEMCELL Technologies Canada Inc. | 1618 Station Street, Vancouver, BC | V6A 1B6 | Canada | [www.stemcell.com](http://www.stemcell.com)

### For Technical Assistance

Tel: +1.604.877.0713

European toll-free number: 00800 7836 2355

North American toll-free number: +1.800.667.0322

email: [techsupport@stemcell.com](mailto:techsupport@stemcell.com)

MDSS GmbH

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30175 Hannover • Germany

Document #10000017274

Version 01

Issue Date 2024-02

## Directions for Use

### PREPARATION OF EASYSEP™ RBC LYSIS BUFFER (1X)

This buffer is supplied as a 10X concentrate. EasySep™ RBC Lysis Buffer (1X) should be prepared at least 1 hour before use.

1. Prepare EasySep™ RBC Lysis Buffer (1X) by adding 1 part EasySep™ RBC Lysis Buffer, 10X Concentrate to 9 parts distilled or Type 1 water.
2. Mix gently and completely before use.

NOTE: If not used immediately, store EasySep™ RBC Lysis Buffer (1X) at 2 - 8°C for up to 3 months. Do not freeze.

\*Type I water refers to ultrapure water suitable for use in analytical procedures. It is defined by the American Society for Testing and Materials (ASTM) as having a resistivity of > 18 MΩ-cm, a conductivity of < 0.056 μS/cm and < 50 ppb of total organic carbons (TOC).

### LYSIS OF RBCs

Specimen: Bone Marrow

1. Prepare the sample within the volume range of 0.5 - 4.5 mL and add sample to required tube. If the volume is > 4.5 mL, split the sample between two 14 mL round-bottom tubes (e.g. Catalog #38008).
2. Add an equal volume of EasySep™ RBC Lysis Buffer (1X) to the sample. For example, if the sample volume is 2 mL, add 2 mL of EasySep™ RBC Lysis Buffer (1X) to the sample. Mix Well.
3. Proceed immediately to cell labeling and separation. Add reagents according to the relevant EasySep™ cell isolation kit PIS.

## Technical Assistance

This PIS is also available in various national languages in electronic format on STEMCELL's website and can be provided in paper copy upon request.

For technical support, contact [techsupport@stemcell.com](mailto:techsupport@stemcell.com) or call toll-free either +1.604.877.0713 (Canada), +1.800.667.0322 (North America), or 00800 7836 2355 (Europe). For more information, visit [www.stemcell.com](http://www.stemcell.com).

Any serious incident that has occurred in relation to the device shall be reported to the manufacturer and the competent authority of the member state in which the user is established.


### For Technical Assistance

Tel: +1.604.877.0713

European toll-free number: 00800 7836 2355

North American Toll-Free Number: 1.800.667.0322

email: [techsupport@stemcell.com](mailto:techsupport@stemcell.com)

MDSS GmbH

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30175 Hannover • Germany

|  |   |   |
|--|---|---|
| <p>Catalog or reference number</p>                                   | <p>Batch code</p>                         | <p>Use-by date:<br/>YYYY-MM-DD</p>                            |
| <p>Caution, consult accompanying documents</p>                       | <p>In Vitro Diagnostic Medical Device</p> | <p>For storage within temperature limits</p>                  |
| <p>Consult Instruction for Use or electronic Instruction for Use</p> | <p>Do not use if packaging is damaged</p> | <p>Unique device identifier</p>                               |
| <p>CE Mark</p>   | <p>Manufacturer</p>                       | <p>Authorized EC representative in the European Community</p> |

EASYSEP™ RED BLOOD CELL LYSIS BUFFER, 10X CONCENTRATE IS MANUFACTURED IN COMPLIANCE WITH 21 CFR PART 820 AND ISO 13485. FOR ADDITIONAL INFORMATION ON QUALITY AT STEMCELL, REFER TO [WWW.STEMCELL.COM/COMPLIANCE](http://WWW.STEMCELL.COM/COMPLIANCE).


Copyright © 2024 by STEMCELL Technologies Inc. All rights reserved including graphics and images. STEMCELL Technologies & Design, STEMCELL Shield Design, Scientists Helping Scientists, and EasySep are trademarks of STEMCELL Technologies Canada Inc. All other trademarks are the property of their respective holders. While STEMCELL has made all reasonable efforts to ensure that the information provided by STEMCELL and its suppliers is correct, it makes no warranties or representations as to the accuracy or completeness of such information.

STEMCELL Technologies Canada Inc. | 1618 Station Street, Vancouver, BC | V6A 1B6 | Canada | [www.stemcell.com](http://www.stemcell.com)

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STEMCELL Technologies Canada Inc. | 1618 Station Street, Vancouver, BC | V6A 1B6 | Canada | www.stemcell.com