

Mammary Epithelial Cells

Standardized Media and Reagents




Table of Contents

3	Mammary Epithelial Cell Types and Assays
4	Mammary Epithelial Cell Product Overview
5	Dissociate Mammary Gland Tissues
6	Detect Viable Mammary Precursor Cells
7	Isolate Human Mammary Epithelial Cells
8	Isolate Mouse Mammary Epithelial Cells
9	Assay Human Mammary Epithelial Cells
10	Culture Human Mammary Epithelial Cells
11	Culture Human Mammospheres or Tumorspheres
13	Culture and Assay Mouse Mammary Epithelial Cells
14	Support Products for Epithelial Cell Research
15	References

Serving Mammary Epithelial Cell Research Scientists

STEMCELL Technologies, a privately-owned biotechnology company, is a leader in specialty cell culture media, cell separation products and ancillary reagents for life science research. As scientists helping scientists, STEMCELL delivers over 1500 products to more than 70 countries worldwide. To learn more about how STEMCELL helps make research work, visit www.stemcell.com.

Front cover: Human mammary bipotent epithelial cell colony cultured in EpiCult™-B (Human) and visualized with Wright-Giemsa. Photo courtesy of Dr. C. Eaves and P. Eirew.

Mammary Epithelial Cell Research

Cell Types and Assays

Emerging evidence suggests that the mammary epithelium in both humans and mice may comprise a hierarchy of cells, spanning from mammary stem cells to differentiated luminal and myoepithelial cells (Figure 1).


In the gold standard *in vivo* assay for mammary stem cells, multi-lineage outgrowths are generated in primary and secondary recipient mice from limiting numbers of cells. Self-renewing cells that generate such multi-lineage outgrowths are termed mammary repopulating units (MRUs).

Mammary progenitor cells can be detected by their ability to generate colonies *in vitro*, either in liquid culture on tissue culture plastic or within 3-dimensional gels such as BD Matrigel™ (page 9). These cells are termed mammary colony-forming cells (Ma-CFCs).

Primitive cells within the mammary epithelium can also be detected with an alternate *in vitro* assay: the mammosphere assay (page 11). Some mammosphere-initiating cells are perceived to be highly primitive precursor cells since clonal mammospheres can form secondary mammospheres and generate multiple lineages of daughter cells.

Support your work with these assays and cell types by leveraging premium, standardized mammary cell products from STEMCELL Technologies (www.stemcell.com).

FIGURE 1. Mammary Epithelial Cell Hierarchy



Why Choose Products From STEMCELL Technologies?

COMPREHENSIVE. Culture, characterize and differentiate primary human and mouse mammary epithelial cells with products based on strong research expertise and specialized technical support.

CONSISTENT. Maximize the reproducibility of your experimental results with fully defined and rigorously tested culture media.

CONVENIENT. Save significant time and effort and minimize experiment troubleshooting with our easy-to-use products and pre-tested protocols, optimized for specific applications.

Matrigel™ is a trademark of BD Biosciences.

Mammary Epithelial Cell Research

Product Overview

Mammary Tissue Dissociation

Dissociation Enzymes
(page 5)

Viable Precursor Cell Detection

**ALDEFLUOR™ Enzymatic Assay
for Detection of ALDH^{br} cells**
(page 6)

Cell Isolation or Enrichment

EasySep™ Cell Separation
(page 7)

Functional In Vitro Assays

**EpiCult™-B for Ma-CFC and
3D Morphogenesis Assays**
(Human page 9; Mouse page 13)

Standardized Culture Media

EpiCult™-C for Short-Term Culture
(Human; page 10)

MammoCult™ for Mammospheres
(Human; page 12)

EpiCult™-B for Short-Term Culture
(Mouse; page 13)



WEBINAR

Dr. John Stingl: Assays for Mammary
Epithelial Stem and Progenitor Cells
www.stemcell.com/StinglAssaysWebinar



WALLCHART

Assays for Human Mammary
Stem and Progenitor Cells
www.stemcell.com/MammaryPoster

Dissociation Enzymes

Dissociate Mammary Gland Tissues

Although different techniques are required for mouse and human mammary tissue dissociation, the samples from both species can be effectively dissociated into a single cell suspension containing viable mammary stem or progenitor cells.

Enzymatic dissociation of mouse mammary organoids usually occurs after 6-8 hours of incubation with a collagenase-hyaluronidase cocktail, which can be extended overnight using our new, gentler cocktail. The tissues are further dissociated into a single cell suspension using Trypsin-EDTA, DNase I and Dispase, and the resulting cell suspension is filtered through a 40 μ m cell strainer to remove clumps. Human breast samples are similarly dissociated with collagenase/hyaluronidase enzymes (standard strength) overnight. The human mammary organoids, obtained through differential centrifugation, can then be dissociated into a single cell suspension using Trypsin/EDTA, Dispase and DNase enzymes.

NEW: Overnight Tissue Dissociation

With the new Gentle Collagenase/Hyaluronidase cocktail that dissociates mouse mammary gland tissues in 12 hours (rather than the usual 6-8), experiments can be completed the next day rather than late at night.

Why Optimize Tissue Dissociation?

Following an optimized dissociation protocol and using quality dissociation enzymes will yield higher cell numbers as well as increased cell viability. Download a guide to the dissociation of solid mammary tissue at www.stemcell.com/tissuedissocguide.

To Dissociate Mammary Gland Tissue

PRODUCT	QUANTITY	CATALOG #
Collagenase/Hyaluronidase (10X)	10 mL	07912
Gentle Collagenase/Hyaluronidase (10X)	10 mL	07919
Trypsin-EDTA	500 mL	07901
DNase I (1 mg/mL)	1 mL	07900
Dispase (5 mg/mL)	100 mL	07913
40 μ m Cell Strainer	50/case	27305

FIGURE 2. Human and mouse mammary tissue samples for dissociation



Freshly extracted mouse mammary glands ready for enzymatic dissociation to a single cell suspension



A human reduction mammoplasty sample

ALDEFLUOR™

Detect Viable Mammary Precursor Cells

ALDEFLUOR™ detects both primitive mammary cells and breast cancer stem cells, based on their expression of aldehyde dehydrogenase (ALDH). Although this fluorescent, non-toxic assay was first developed for hematopoietic stem cells, optimized protocols for assaying some normal and malignant mammary cells with ALDEFLUOR™ are available from techsupport@stemcell.com.

Why Choose ALDEFLUOR™?

- The only non-immunological assay for detecting primitive mammary cells.
- Non-destructive, allowing for the downstream analysis of ALDH^{br} cells.
- Compatible with immunophenotyping and counterstaining of cells with other fluorescent antibodies.

To Detect Viable Precursor Cells

PRODUCT: ALDEFLUOR™ Kit
CATALOG #: 01700 40 Tests/Kit


RECOMMENDED FOR:

Identification, enumeration and isolation of viable normal and cancer stem/progenitor cells based on their ALDH activity.

CONTAINS:

- ALDEFLUOR™ (Dry reagent, 50 µg)
- DEAB (1.5 mM Diethylaminobenzaldehyde in 95% ethanol, 1 mL)
- HCl (2N Hydrochloric Acid, 1.5 mL)
- DMSO (Dimethylsulphoxide, 1.5 mL)
- ALDEFLUOR™ Assay Buffer (4 x 25 mL)



FIGURE 3. SKBR3 breast cancer cells visualized with ALDEFLUOR™



FACS Profile Results With ALDEFLUOR™ Kit

DEAB Control: 1.34% ALDH^{br} cells

ALDEFLUOR™: 9.55% ALDH^{br} cells




Starting with primary human mammary epithelial cells, ALDH^{br} cell content of the selected cells typically ranges around 10%.

Technical Bulletin: ALDEFLUOR™ and Cancer Stem Cells

Obtain a technical bulletin detailing publications, technical data and protocols on using the ALDEFLUOR™ assay with the MCF7 and SKBR3 breast cancer cell lines at www.stemcell.com.

EasySep™ Isolate Human Mammary Epithelial Cells

FIGURE 4. Schematic of EasySep™ immunomagnetic cell isolation.




Investigating Circulating Tumor Cells (CTCs)?

Try the RosetteSep™ Human Circulating Epithelial Tumor Cell Enrichment Cocktail (Catalog #15127).


To Isolate Human CD10⁺ Cells

The CD10⁺ cell fraction isolated by EasySep™ Human CD10 Positive Selection Kit (Catalog #18358) contains bipotent and myoepithelial-restricted progenitors, as well as non-clonogenic myoepithelial cells.

Start: 45.6% CD10⁺ Cells



Selected: 95.0% CD10⁺ Cells




The CD10⁺ cell content of the enriched fraction typically ranges from 84 - 98%. The purity of the CD10⁺ cells can be assessed with FITC-labeled anti-dextran (which recognizes the dextran on the magnetic particle).


To Isolate Human EpCAM⁺ Cells

The EpCAM⁺ cell fraction isolated by EasySep™ Human EpCAM Positive Selection Kit (Catalog #18356) contains bipotent and luminal-restricted progenitors, as well as luminal cells.

Start: 49.5% 5E11⁺ Cells



Selected: 94.6% 5E11⁺ Cells




Starting with cultured mammary tissue, the typical purity of the selected fraction is 90 - 96%. Purity has been assessed by staining with the monoclonal antibody 5E11-FITC (Catalog #10110) which has a distribution identical to EpCAM.


To Isolate Human MUC1⁺ Cells

The MUC1⁺ cell fraction isolated by EasySep™ Human MUC1 Positive Selection Kit (Catalog #18359) contains luminal-restricted progenitors as well as non-clonogenic luminal epithelial cells.

Start: 19.5% MUC1⁺ Cells



Selected: 94.9% MUC1⁺ Cells




Starting with cultured mammary tissue, the typical purity of the enriched fraction is 95 - 99%. The purity of MUC1⁺ cells can be assessed with FITC-labeled anti-dextran (which recognizes the dextran on the magnetic particle).

EasySep™ Isolate Mouse Mammary Epithelial Cells

To Enrich For Mouse Mammary Stem Cells

EasySep™ Mouse Mammary Stem Cell Enrichment Kit (Catalog #19757) is optimized for the isolation of CD31⁺CD45⁺TER119⁺CD24⁺CD49f⁺ population; the sub-population of CD24⁺CD49f⁺ is enriched for mammary epithelial stem cells and is isolated via fluorescence-activated cell sorting.




Flow cytometry profile of freshly dissociated mouse mammary epithelial cells distributed according to their expression of CD49f and CD24. The mammary repopulating unit (MRU)-rich fraction is found in the CD49f⁺CD24⁺ subset. The mammary colony-forming cells (Ma-CFCs) are found in the CD49f⁺CD24⁺ subset.²


To Enrich For Mouse Epithelial Cells

EasySep™ Mouse Epithelial Cell Enrichment Kit (Catalog #19758) removes contaminating hematopoietic and endothelial cells from cell preparations. Cells expressing CD45, Ter119, BP-1 or CD31 cell surface antigens are depleted, leaving behind an unlabeled epithelial cell-enriched population.

Before Depletion:
10.37% ± 0.50% unwanted cells



After Depletion:
1.12% ± 0.22% unwanted cells



Purity has been assessed by staining with goat anti-rat IgG FITC, which recognizes the antibodies used to deplete unwanted cells expressing CD45, Ter119, BP-1 and/or CD31.

EasySep™ Mammary Epithelial Cell Isolation Kit Listing

EASYSEP™ SELECTION KIT	SPECIES	COMPONENTS	CAPACITY	CATALOG #
CD10 Positive Selection	Human	EasySep™ Human CD10 Positive Selection Cocktail EasySep™ Magnetic Particles	1 x10 ⁹ cells	18358
EpCAM Positive Selection	Human	EasySep™ Human EpCAM Positive Selection Cocktail EasySep™ Magnetic Particles	1 x10 ⁹ cells	18356
MUC1 Positive Selection	Human	EasySep™ Human MUC1 Positive Selection Cocktail EasySep™ Magnetic Particles	1 x10 ⁹ cells	18359
Mammary Stem Cell Enrichment	Mouse	EasySep™ Mouse Mammary Epithelial Cell Enrichment Cocktail EasySep™ Biotin Selection Cocktail EasySep™ Magnetic Particles Anti-CD24-PE, Anti-CD49f-FITC	1 x10 ⁹ cells	19757
Epithelial Cell Enrichment	Mouse	EasySep™ Mouse Epithelial Cell Enrichment Cocktail EasySep™ Biotin Selection Cocktail EasySep™ Magnetic Particles	1 x10 ⁹ cells	19758


EpiCult™-B (Human)




Assay For Human Mammary Epithelial Cells

EpiCult™-B is a versatile medium that can assay for the presence of human mammary colony-forming cells (Ma-CFCs) in either adherent monolayer or in 3D Matrigel™-based cultures. Quantifying Ma-CFC frequency provides a functional measure of the number of mammary progenitor cells in the sample, while the 3D morphogenesis assay enables more representative studies of the in vivo mammary duct structure.


The Ma-CFC Assay



Mixed Epithelial Colony



Luminal Epithelial Colony



Myoepithelial Colony

Why Choose EpiCult™-B (Human)?

- Fully defined culture medium.
- Promotes clonogenic growth of unipotent and bipotent progenitors in the Ma-CFC assay.
- Promotes mammary progenitor cell differentiation in the 3D morphogenesis assay.

To Assay For Mammary Colony-Forming Cells

PRODUCT:	EpiCult™-B (Human)
CATALOG #:	05601 500 mL

RECOMMENDED FOR:

Culture of human mammary luminal and myoepithelial cells. Promotes the clonogenic growth of bipotent, luminal-restricted and myoepithelial-restricted human mammary epithelial progenitors when used in conjunction with an irradiated feeder layer (NIH 3T3). Also used for the enzymatic dissociation of human mammary tissue when supplemented with Collagenase/Hyaluronidase (Catalog #07912).

CONTAINS:

- EpiCult™-B Basal Medium (Human)
- EpiCult™-B Supplement (Human)

REQUIRES:

Supplement with 0.48 µg/mL freshly dissolved Hydrocortisone (Catalog #07904) before use.

The 3D Morphogenesis Assay



3D Mammary Structure

EpiCult™-C (Human) Culture Human Mammary Epithelial Cells

EpiCult™-C is a serum-free liquid culture medium optimized for the short term culture of primary human mammary epithelial cells as adherent monolayers (Figure 5). It provides balanced, robust growth of both myoepithelial and luminal cells.

Why Choose EpiCult™-C?

- Fully defined culture medium.
- Promotes robust growth of luminal and myoepithelial cells.
- Uniquely optimized for the short-term culture of primary mammary epithelial cells.

To Culture Human Mammary Epithelial Cells

PRODUCT: EpiCult™-C
CATALOG #: 05630 500 mL

RECOMMENDED FOR:

Routine monolayer culture of primary human mammary epithelial cells and a variety of breast cancer cell lines.


CONTAINS:

- EpiCult™-C Basal Medium (Human)
- EpiCult™-C Proliferation Supplements (Human)

REQUIRES:

Supplement with 0.48 µg/mL freshly dissolved Hydrocortisone (Catalog #07904) before use.

FIGURE 5. EpiCult™-C Day 6 culture



Mammary human epithelial cells cultured in EpiCult™-C for 6 days and visualized with Wright-Giemsa



VIDEO


Introducing EpiCult™-C Medium
www.stemcell.com/EpiCultCVideo

MammoCult™




Culture Human Mammospheres or Tumorspheres

Mammosphere culture involves seeding mammary epithelial cells at low density in an environment that prevents their adherence to a substratum and enables their proliferation in suspension as spherical clusters. Primitive mammary cells have been detected over multiple generations of mammospheres, and cells from dissociated mammospheres demonstrate multi-lineage differentiation potential in the Ma-CFC assay. The mammosphere culture system is thus capable of propagating mammary stem and progenitor cells in vitro.



Mammospheres can be classified into different groups based on morphology



Hollow Mammosphere

This hollow mammosphere was cultured in suspension in MammoCult™ for 7 days.



Solid Mammospheres

These solid mammospheres were cultured in suspension in MammoCult™ for 7 days.



Large and Small Mammospheres

These large ($>70\ \mu\text{m}$) and small ($<70\ \mu\text{m}$) mammospheres were cultured in suspension in MammoCult™ for 7 days.

Protocol: Tumorsphere Culture

Cells from breast cancer cell lines may also be grown in nonadherent liquid suspension cultures. Obtain a technical bulletin on tumorsphere culture with detailed protocols at www.stemcell.com.


MammoCult™ **Culture Human Mammospheres or Tumorspheres**

MammoCult™ is the only commercially available medium for the culture of mammospheres and tumorspheres (Figures 6 and 7). It supports robust cultures that produce significantly larger numbers of mammospheres and tumorspheres than published formulations. Use this defined culture medium to maintain primary mammary epithelial progenitors or breast cancer cell lines.

Why Choose MammoCult™?


- Fully defined culture medium.
- Generates a large number of robust mammospheres and tumorspheres.
- Supports short-term maintenance of primary mammary epithelial progenitors.
- Supports long-term maintenance of breast cancer cell lines MCF7, SKBR3, MDA -MB-231, AU 565, SUM149 and BT474.

FIGURE 6. MammoCult™ mammosphere culture



Mammospheres derived from primary mammary tissue and cultured in liquid suspension in MammoCult™ for 8 days.

FIGURE 7. MammoCult™ tumorsphere culture



Tumorspheres derived from the MCF-7 breast cancer cell line and cultured in liquid suspension in MammoCult™ for 8 days.

To Culture Mammospheres and Tumorspheres from Human Mammary Cells

PRODUCT:	MammoCult™
CATALOG #:	05620 500 mL

RECOMMENDED FOR:

Culture of primary human mammary stem and progenitor cells and breast cancer cell lines as mammospheres and tumorspheres.

CONTAINS:

- MammoCult™ Basal Medium (Human)
- MammoCult™ Proliferation Supplement (Human)

REQUIRES:

Supplement with 0.48 µg/mL freshly dissolved Hydrocortisone (Catalog #07904) and 4 µg/mL Heparin (Catalog #07980) before use.

MammoCult™ Scientific Poster

MammoCult™ outperforms the published formulation (Dontu et al. 2003), by maintaining a larger number of mammospheres and preserving both unipotent and bipotent progenitor cells longer.

Download the scientific poster under Technical Resources at www.stemcell.com.

EpiCult™-B (Mouse) Culture and Assay Of Mouse Mammary Epithelial Cells

EpiCult™-B (Mouse) is a robust culture medium for mouse mammary epithelial cells (Figure 8), and is ideal for assaying the presence of mouse mammary colony-forming cells (Ma-CFCs) in either adherent monolayer or in 3D Matrigel™-based cultures. As this medium is formulated specifically for mouse cells, it delays the epithelial-to-mesenchymal transition (EMT) commonly observed in the mouse mammary epithelial in vitro cultures (Figure 9).

Why Use EpiCult™-B (Mouse)?

- Defined culture medium for primary mammary epithelial cells (short-term).
- Promotes clonogenic growth of mammary epithelial cells in the Ma-CFC assay.
- Promotes mammary progenitor cell growth in the 3D morphogenesis assay.
- Delays and decreases EMT occurrence in cell cultures, compared to the traditional medium.

To Assay For Mammary Colony-Forming Cells

PRODUCT:	EpiCult™-B (Mouse)
CATALOG #:	05610 500 mL

RECOMMENDED FOR:

Colony-forming cell assays with mouse mammary epithelial cells. Promotes the clonogenic growth of mouse mammary progenitors when used in conjunction with an irradiated feeder layer (NIH 3T3). It is also used for the enzymatic dissociation of mouse mammary tissue when supplemented with Collagenase/Hyaluronidase (Catalog #07912).


CONTAINS:

- EpiCult™-B Basal Medium (Mouse)
- EpiCult™-B Proliferation Supplements (Mouse)

REQUIRES:


Supplement with 10 ng/mL rh Epidermal Growth Factor (Catalog #02633), 10 ng/mL rh Basic Fibroblast Growth Factor (Catalog #02634) and 4 µg/mL Heparin (Catalog #07980) before use.

FIGURE 8. EpiCult™-B mouse mammary epithelial cell colonies



(A) EpiCult™-B (Mouse) generates twice as many colonies compared to commonly used formulation. (B) and (C) are examples of mouse mammary epithelial colonies grown in EpiCult™-B (Mouse).

FIGURE 9. EpiCult™-B (Mouse) cell cultures display decreased EMT



Mammary epithelial cells display a clear delay and decrease of epithelial-to-mesenchymal transition (EMT) when cultured in (A) EpiCult™-B (Mouse) versus (B) traditional formulation (photo courtesy of Dr. John Stingl). EMT is defined by emergence of a fibroblastic-like morphology.

Antibodies and Miscellaneous Support Products

Antibodies

ANTIGEN	CLONE	EXPRESSION	ISOTYPE	FORMAT*	APPLICATIONS**	QUANTITY	CATALOG#
CD10 (CALLA)	FR4D11	Basal-positioned cells	Mouse IgG ₁	Purified	CS, IH	100 µg	01431
CD24	32D12	Luminal cells	Mouse IgG ₁ , κ Mouse IgG ₁ , κ Mouse IgG ₁ , κ	Purified Biotin FITC	CS CS CS	100 µg 100 tests 100 tests	01434 10231 10424
CD44	IM7	Human mammary tumor stem cells	Rat IgG _{2b} , κ Rat IgG _{2b} , κ Rat IgG _{2b} , κ	Purified Biotin FITC	CS CS CS	50 µg 500 µg 50 µg	01465 10232 10432
CD49f (α6 Integrin)	GoH3	Basal-positioned cells, human and mouse progenitors 1,3,4 and mouse stem cells	Rat IgG _{2a}	FITC	FC	100 tests	10111
CD90 (Thy1)	5E10	Myoepithelial progeny of human mammary epithelial cells	Mouse IgG ₁ , κ Mouse IgG ₁ , κ	Purified FITC	FC, CS FC, CS	100 µg 100 tests	01437 10427
EpCAM	VU-1D9	Bipotent and luminal-restricted progenitors and luminal cells	Mouse IgG ₁ Mouse IgG ₁	Purified FITC	CS, IH FC, IF	100 µg 100 tests	01420 10109
	323/A3	Same as VU-1D9 but 323/A3 has a slightly broader distribution in basal compartment	Mouse IgG ₁	Purified	CS, IH	100 µg	01421
Epithelial Cell Surface	5E11	Luminal cells. Variable staining of malignant breast epithelium	Mouse IgG ₁ Mouse IgG ₁	Purified FITC	CS, IH FC	100 µg 100 tests	01422 10110
Keratin 18	DC-10	Luminal cells	Mouse IgG ₁	Purified	IH	100 µg	01425
Keratin 19	A53-B/A2	Subpopulation of luminal cells	Mouse IgG _{2a} , κ	Purified	IH	100 µg	01432
Ly-6A/E (SCA1)	E13-161.7	Mouse mammary stem cells (weak) and some progenitors	Rat IgG _{2a} , κ Rat IgG _{2a} , κ	FITC PE	FC, CS FC, CS	500 µg 200 µg	10716 10816
MUC1 Glycoprotein	214D4	Luminal-restricted progenitors and luminal cells	Mouse IgG ₁	Purified	CS, IH	200 µg	01423
Progesterone Receptor	SP2	Subpopulation of luminal cells	Rabbit IgG	Supernatant	IH	500 µg	01466

Format:** FITC - Fluorescein Isothiocyanate; PE - Phycoerythrin *Applications:** CS - Cell Separation; FC - Flow Cytometry; IH - Immunocytochemistry IF - Immunofluorescence

Support Products

PRODUCT	QUANTITY	CATALOG #
Collagen Solution (Bovine)	35 mL	04902
Hanks' Balanced Salt Solution Modified	500 mL	37150
Tissue Dissociation Flask (250 mL)	1 flask	27300
Tissue Culture Dishes (35 mm)	20 dishes/pack	27114
Ultra-Low Adherent Plates (6-well plates)	5 plates/pack	27145

References

1. Stingl J, et al., Purification and Unique Properties of Mammary Epithelial Stem Cells. *Nature* 439:993-997, 2006
2. Stingl J. Detection and Analysis of Mammary Gland Stem Cells. *J Pathology* 217:229-241, 2009
3. Stingl J, et al., Characterization of Bipotent Mammary Epithelial Progenitor Cells in Normal Adult Human Breast Tissue. *Breast Cancer Research* 2:93-109, 2001
4. Viladsen R, et al., Evidence for a Stem Cell Hierarchy in the Adult Human Breast. *J. Cell Biol* 1:87-101, 2007
5. DiGiovanna MP. Activation State-Specific Monoclonal Antibodies Detects Tyrosine Phosphorylated p185neu/erbB-2 in a Subset of Human Breast Tumors Over-expressing this Receptor. *Cancer Res* 9:1846-55, 1995
6. Shackelton M, et al., Generation of a Functional Mammary Gland from a Single Stem Cell. *Nature* 439:84-88, 2006
5. Shafee N, et al., Cancer Stem Cells Contribute to Cisplatin Resistance in Brca1/p53-Mediated Mouse Mammary Tumors. *Cancer Res* 2008; 68: 3243
6. Eirew P, et al., A Method for Quantifying Normal Human Mammary Epithelial Stem Cells with in vivo Regenerative Ability. *Nature Med* 12:1384-1389, 2008
7. Raouf A, et al., Transcriptome Analysis of the Normal Human Mammary Cell Commitment and Differentiation Process. *Cell Stem Cell* 1:10-118, 2008
8. Stingl J, et al., Purification and Unique Properties of Mammary Epithelial Stem Cells. *Nature* 439:993-997, 2006
9. Stingl J, et al., Methods in Mammary Gland Biology and Breast Cancer (MM Ip and BB Asc. Eds), Kluwer Academic Publishers, New York, pg 177-193, 2000
10. Cooper NL, et al., Correlation of CD44 Expression with Proliferative Activity of Normal Human Breast Epithelial Cells in Culture. *Breast Cancer Res Treat* 2:143-153, 1998
11. Bardy P, et al., Isolation and Analysis of Different Subpopulation of Normal Human Breast Epithelial Cells after their Infection with a Retroviral Vector Encoding a Cell Surface Marker. *Breast Cancer Res Treat* 2:153-165, 1997
12. Claudio JA and Emerman JT. The Effects of Cyclosporin A, Tamoxifen and Medroxyprogesterone acetate on the Enhancement of Adriamycin Cytotoxicity in Primary Cultures of Human Breast Epithelial Cells. *Breast Cancer Res Treat* 2:111-112, 1996
13. Emerman JT, et al., Selective Growth of Freshly Isolated Human Breast Epithelial Cells Cultured at Low Concentrations in the Presence or Absence of Bone Marrow Cells. *Breast Cancer Res Treat* 2:147-159, 1996

Publications Citing EpiCult™-B and MammoCult™

1. Deng S, et al., Distinct Expression Levels and Patterns of Stem Cell Marker, Aldehyde Dehydrogenase Isoform 1 (ALDH1), in Human Epithelial Cancers. *PLoS One* 2010; 5(4): e10277
2. Shimono Y, et al., Down-regulation of miRNA-200c Links Breast Cancer Stem Cells with Normal Stem Cells. *Cell* 2009; 138(3): 592-603
3. Yamagi D, et al., Development of mammary luminal progenitor cells is controlled by the transcription factor STAT5A. *Genes Dev* 23(20): 2382-2387, 2009
4. Diehn M, et al., Association of Reactive Oxygen Species Levels and Radio-resistance in Cancer Stem Cells. *Nature* 458(7239):780-3, 2009

Mammary Training Course

Enhance your understanding of mammary stem and progenitor cells in this annual two day training course, offered in collaboration with the Jackson Laboratory. Practical laboratory sessions include in vitro and in vivo functional assessment of mammary stem and progenitor cells. Contact education@stemcell.com for more information.

Scientists Helping Scientists™ | WWW.STEMCELL.COM

TOLL-FREE PHONE 1 800 667 0322 • PHONE 1 604 877 0713

TOLL-FREE FAX 1 800 567 2899 • FAX 1 604 877 0704

ORDERS@STEMCELL.COM • INFO@STEMCELL.COM

FOR FULL CONTACT DETAILS WORLDWIDE VISIT OUR WEBSITE

FOR RESEARCH USE ONLY. NOT FOR THERAPEUTIC OR DIAGNOSTIC USE. DOCUMENT #29188 VERSION 1.2.2 NOV 2012

