

After EasySep™ Separation

Performance Analysis

Presenter

LU LU

SCIENTIFIC SUPPORT SPECIALIST, IMMUNOLOGY



Learning Objectives

In this session, you will learn:

How to assess the separation performance after EasySep™

- Purity analysis based on flow cytometry
 - General considerations for fluorochrome choices
 - General considerations for staining controls
 - General considerations for gating
 - Specific considerations for flow cytometry analysis of EasySep™ separated cells
- Recovery analysis

Purity Analysis Based on Flow Cytometry

Property of
STEMCELL Technologies

Purity Analysis Based on Flow Cytometry

EasySep™ Human Memory B Cell Isolation Kit




Notes and Tips

ASSESSING PURITY

For purity assessment of CD19+CD27+ cells by flow cytometry, use the following fluorochrome-conjugated antibody clones:




- Anti-Human CD19 Antibody, Clone HB19 (Catalog #60005), and
- Anti-Human CD45 Antibody, Clone H30 (Catalog #60018) or Clone D11 (Catalog #60123), and
- Anti-Mouse CD27 Antibody, Clone LG.3A10 (Catalog #60160), which cross-reacts with human CD27

NOTE: Brilliant Violet™ antibody conjugates should be carefully titrated on EasySep™ Release-isolated cells prior to analysis by flow cytometry or fluorescence microscopy. For purity assessment with Brilliant Violet™ antibody conjugates, use of BD Horizon Brilliant™ Stain Buffer is recommended to reduce non-specific interactions. For more information, refer to the manufacturer's instructions or contact us at techsupport@stemcell.com.



Starting with PBMCs, the memory B cell content (CD19+CD27+) of the isolated fraction is typically $97 \pm 2\%$ (mean \pm SD using the purple EasySep™ Magnet). Using the optional protocol, the naïve B cell content (CD19+CD27-) of the isolated fraction is typically $93 \pm 5\%$ (mean \pm SD).




PRODUCTS ARE FOR RESEARCH USE ONLY AND NOT INTENDED FOR HUMAN OR ANIMAL DIAGNOSTIC OR THERAPEUTIC USES UNLESS OTHERWISE STATED. FOR ADDITIONAL INFORMATION ON QUALITY AT STEMCELL, REFER TO WWW.STEMCELL.COM/COMPLIANCE.
Copyright © 2022 by STEMCELL Technologies Inc. All rights reserved including graphics and images. STEMCELL Technologies & Design, STEMCELL Shield Design, Scientists Helping Scientists, EasySep, EasySep, EasySep, and SepMate are trademarks of STEMCELL Technologies Inc. Brilliant Violet is a trademark of Becton Group LLC. Lymphoprep is a trademark of Stemware Bering AG. All other trademarks are the property of their respective holders. While STEMCELL has made all reasonable efforts to ensure that the information provided by STEMCELL and its suppliers is correct, it makes no warranties or representations as to the accuracy or completeness of such information.



Flow Cytometry - Considerations for Fluorochromes

- **Spectrum**

- Compatibility with your flow cytometer
- Overlap between the fluorochromes



Solid peaks: Cells are stained with mouse anti-human CD3 FITC, clone UCHT1.

Left: Signal in FITC channel

Right: Signal in PE channel

Black outlined peaks: Unstained cells


- **Brightness**

- Bright fluorochrome for a target marker with a low expression level or target cells with a low frequency


Flow Cytometry - Staining Controls

- **Unstained control**

Unstained control shows the background autofluorescence.




(A) Unstained PBMCs;




(B) PBMCs stained with mouse anti-human CD3 (clone UCHT1) APC

- **Isotype control**

Isotype controls show non-specific binding and allow optimization of staining.



(C) Mouse splenocytes labeled with an Armenian hamster IgG Alexa Fluor® 488 isotype control antibody and an anti-mouse MHC class II antibody, APC.




(D) Mouse splenocytes labeled with Anti-Mouse CD11c Antibody, Clone N418, Alexa Fluor® 488 and an anti-mouse MHC class II antibody, APC.

Flow Cytometry - Considerations for staining controls

- **Single stained control**

Single stained controls show spectral overlap and allow fluorescence spillover to be corrected.



PBMCs stained with mouse anti-human CD3 (clone UCHT1) BV605. (A) Uncompensated; (B) Compensation applied.




Flow Cytometry - Considerations for Staining Controls

- **Fluorescence Minus One (FMO) Control**

FMO controls show the fluorescence spillover from all other fluorophores to the one being analyzed and help you determine where to set your gates.

Dot plots of multicolor flow cytometry showing the fluorescence spread into the PE channel (FMO control compared to a no stain control). Black dotted line represents the FMO gating boundary compared to the no stain boundary in blue.



The sample for FMO-PE control is stained with anti-CD45-FITC, anti-CD56-APC antibodies, and PI; the sample for the All Stain group is stained with anti-CD3-PE, anti-CD45-FITC, anti-CD56-APC antibodies, and PI.



Flow Cytometry Analysis - Considerations for Gating


1. Check FSC-SSC plot and exclude debris

2. Exclude doublets or multiplets




Flow Cytometry Analysis - Considerations for Gating


3. Exclude non-viable cells



4. Include a leukocyte marker






Flow Cytometry Analysis - Gating strategy



Flow Cytometry Analysis - Percentage of the Target Cells

The percentage of target cells among live single cells is used to measure the cell percentage. Measuring the percentage of target cells among live single leukocytes provides more accurate data when working with immune cells.

- For target cell types that are defined by one or two markers, you can use the percentage of cells gated based on the markers among live and single leukocytes.
- For target cell types that are defined by more than two markers, the purity is the product of multiplying all the percentages of each marker.



X

0.026%

Flow Cytometry Analysis for EasySep™ Separated Cells

- Check the Notes and Tips in the Product Information Sheets and use recommended antibody clones



Directions for Use – Fully Automated RoboSep™ Protocol


See page 1 for Sample Preparation and Recommended Medium. Refer to Table 1 for detailed instructions regarding the RoboSep™ procedure.

Table 1. RoboSep™ Human CD3 Positive Selection Kit II Protocol

STEP	INSTRUCTIONS	Reagents (Catalog #s and volumes)
1	Prepare samples in the RoboSep™ and incubation within the isolation range.	1 x 100 µL cells
2	Add samples to required tubes.	14 mL (17 x 85 mm) polypropylene round-bottom tube (e.g., Corning®)
3	Mount pipettes.	Human CD3 Positive Selection II 11001
4	Verify Reagent/Program™.	30 seconds
5	Mount Reagent/Program™.	30 seconds
6	Start the program.	Follow on-screen prompts
7	Unblock the pipettes when the run is complete. Remove the tubes containing the isolated cells and reagent to desired location. Set the RoboSep™ to "Off" and remove the tubes from the system.	Press the green "Stop" button
8	Isolated cells are ready for use.	

Notes and Tips

The EasySep™ Human CD3 Positive Selection Cocktail uses an anti-CD3 antibody clone that to our knowledge fully or partially blocks all anti-CD3 antibody clones used to assess purity by flow cytometry. For purity assessment of CD3+ cells by flow cytometry, use fluorochrome-conjugated:
 • Anti-Human CD3 Antibody, Clone UCHL1 (Catalog #60011; partially blocked)
 • Anti-Human CD3 Antibody, Clone OKT4 (Catalog #60016; partially blocked)
 • Anti-Human CD3 Antibody, Clone RPA-T8 (Catalog #60022)
 • Use a fluorochrome-conjugated secondary antibody, such as Goat Anti-Mouse IgG (H+L) Antibody, Polyclonal (Catalog #60138)



Starting with a single cell suspension of human PBMCs, the CD3+ cell content of the isolated fraction is typically 80% ± 10% (mean ± SD) using the pure EasySep™ Aligned™. In the above examples, the portion of the start and final isolated fractions are 50.5% and 80.0%, respectively.

STEMCELL TECHNOLOGIES, INC. 10000 DEER CREEK DRIVE, SUITE 100, DENVER, CO 80231, USA. PRODUCTS ARE FOR RESEARCH USE ONLY AND NOT INTENDED FOR HUMAN OR ANIMAL DIAGNOSTIC OR THERAPEUTIC PURPOSES. © 2013 STEMCELL TECHNOLOGIES, INC.

Notes and Tips

ASSESSING PURITY

The EasySep™ Human CD3 Positive Selection Cocktail uses an anti-CD3 antibody clone that to our knowledge fully or partially blocks all anti-CD3 antibody clones used to assess purity by flow cytometry. For purity assessment of CD3+ cells by flow cytometry, use fluorochrome-conjugated:

- Anti-Human CD3 Antibody, Clone UCHL1 (Catalog #60011; partially blocked)

One of the following methods can also be used:

- Use alternative markers such as fluorochrome-conjugated Anti-Human CD4 Antibody, Clone OKT4 (Catalog #60016) and Anti-Human CD8a Antibody, Clone RPA-T8 (Catalog #60022).
- Use a fluorochrome-conjugated secondary antibody, such as Goat Anti-Mouse IgG (H+L) Antibody, Polyclonal (Catalog #60138).

- If cells are isolated with **EasySep™ Release** products: Brilliant Violet™-conjugated antibodies should be carefully titrated on the cells before analysis. The BD Horizon Brilliant™ Stain Buffer is recommended to reduce non-specific interactions.

Recovery Analysis

Property of
STEMCELL Technologies

Recovery Calculation

- To accurately calculate the recovery of a cell isolation procedure, you will need the following four pieces of data. Accurate determination of each of these components is critical:

#1 Total # of cells that will be separated

Started with 5×10^7 cells

#2 Starting % of desired cells

The frequency of desired cells in the starting sample is 30%

#3 Total # of cells in enriched fraction

Enriched 7.5×10^6 cells

#4 % Purity of enriched fraction

The purity of desired cells in the enriched fraction is 95%

- Cell recovery is calculated using the following formula:

$$\% \text{ Recovery} = [(\#3) \times (\#4)] / [(\#1) \times (\#2)] \times 100\%$$

$$\% \text{ Recovery} = [(7.5 \times 10^6) \times 95\%] / [(5 \times 10^7) \times 30\%] = 47.5\%$$

Summary

- There are a number of factors to consider and incorporate into your flow cytometry staining and gating strategy, including:
 - Using appropriate antibodies and dyes for staining
 - Setting necessary controls
 - Using an optimal staining protocol and flow cytometer
 - Applying appropriate gating strategy
- For cells isolated by EasySep™ positive selection products, it is important to analyze them with the antibody clone recommended in the Product Information Sheet.
- For cells isolated with EasySep™ Release products: Brilliant Violet™ conjugated antibodies should be carefully titrated on the cells before analysis. The BD Horizon Brilliant™ Stain Buffer is recommended to reduce non-specific interactions.
- You can determine the cell recovery of the isolation procedure by assessing the cell counts and frequencies of desired cells in the sample before and after isolation.

Questions?

Scientists Helping Scientists™ | TECHSUPPORT@STEMCELL.COM

Thank you!

Scientists Helping Scientists™ | WWW.STEMCELL.COM

Please complete the quiz and the final assessment.