## Applications of Hepatic Organoids

Lecture 4

#### **Presenter**

Kiran Bhullar Scientist, Scientific Support, Epithelial Cell Biology





## **Learning Objectives**

After this session, you should be able to:

Understand the key research applications for human hepatic organoids and possible experimental endpoints

(EMC)

oropertyots



### **Outline**

- **1.** Key Downstream Applications
- **2.** Application Highlights
- Hepatotoxicity
- Disease Modelling
- Gene Editing





#### Section 1 | Overview of Downstream Applications





- Liver development and function in healthy conditions
- Study mechanisms of disease (disease modelling)

#### **Personalised Medicine**

- Patient-specific drug screening
- Gene therapy (gene editing)

#### **Biobanking**

- Regenerative medicine (transplants)
- Drug screening (patient-derived organoids)
- Toxicology studies (drug-induced liver Injury/DILI)



Figure adapted from Prior et al. Gut 2019

### Section 2 | Application Highlights



#### **Application 1: Hepatotoxicity**



## **Drug-Induced Liver Injury (DILI)**





STEMCELL<sup>™</sup>

Weaver et al. Nature Reviews Drug Discovery 2016



S

## Hepatic Organoids as a Model for Hepatotoxicity



## Hepatic Organoids as a Model for Hepatotoxicity

Gastroenterology 2021;160:831-846

#### BASIC AND TRANSLATIONAL—LIVER

## High-Fidelity Drug-Induced Liver Injury Screen Using Human Pluripotent Stem Cell–Derived Organoids

**Tadahiro Shinozawa**,<sup>1,\*</sup> **Masaki Kimura**,<sup>1,\*</sup> **Yuqi Cai**,<sup>1,\*</sup> Norikazu Saiki,<sup>2</sup> Yosuke Yoneyama,<sup>2</sup> Rie Ouchi,<sup>1</sup> Hiroyuki Koike,<sup>1</sup> Mari Maezawa,<sup>2</sup> Ran-Ran Zhang,<sup>1</sup> Andrew Dunn,<sup>1</sup> Autumn Ferguson,<sup>1</sup> Shodai Togo,<sup>1</sup> Kyle Lewis,<sup>1</sup> Wendy L. Thompson,<sup>1</sup> Akihiro Asai,<sup>1</sup> and Takanori Takebe<sup>1,2,3,4</sup>





# Cytotoxicity Screening Using Hepatic Organoids Cultured in HepatiCult<sup>™</sup> (Human)

A Drug screening using proliferative organoids maintained in HepatiCult™ OGM



B Drug screening using mature organoids differentiated in HepatiCult™ ODM



Schematic of drug screening and cell viability assay protocol

Scientific Poster HepatiCult™ for Human Hepatic Organoids: A Validated Culture System for Drug Toxicity Screening

#### Protocol

How to Perform Cytotoxicity Screening Using Hepatic Organoids Cultured in HepatiCult™ (Human)







Summary of drug toxicity effects in human hepatic proliferative (OGM) and mature (ODM) organoids, HepG2 cells, and PHH.

Scientific Poster HepatiCult™ for Human Hepatic Organoids: A Validated Culture System for Drug Toxicity Screening



## **Applications Using hPSC-Derived Hepatic Cells**

#### Hepatotoxicity testing using 2D and 3D hepatic models

**Objective:** Assess effects of compounds with known hepatotoxicity on 2D HLCs and differentiated HLC-derived organoids

10



## **Applications Using hPSC-Derived Hepatic Cells**

#### hPSC-derived hepatic cells are sensitive to compound-induced hepatotoxicity



HLCs = hepatocyte-like cells; ODM = Organoid Differentiation Medium; PHHs = primary human hepatocytes

## Hepatic Organoids as a Model for Hepatotoxicity

- Liver toxicity is a major cause of clinical drug attrition and acute liver failure
- All drugs (modality- and indication-agnostic) are tested for hepatotoxicity
- Compatibility of hepatic organoids (generated using HepatiCult<sup>™</sup>) with hepatotoxicity testing:
  - Hepatic organoids show amenability/compatibility with hepatotoxicity testing
  - Seven drugs were tested for IC50 values and compared to published IC50 values for HepG2 cell lines and primary human hepatocytes (PHH)
  - Differentiated (ODM) organoids showed higher sensitivity to tested drugs, suggesting an effect on more mature cell populations



#### **Application 2: Disease Modelling**





Image adapted from Nature Oct 2015 illustration and Inspero webinar

## Hepatic Organoids for Disease Modelling (NASH) (NASH)

## HEPATOLOGY



ORIGINAL ARTICLE

#### Functional characterization of organoids derived from irreversibly damaged NASH patient liver

Sarah McCarron, Brooke Bathon, Donna M. Conlon, Deepti Abbey, Daniel J. Rader, Katerina Gawronski, Christopher D. Brown, Kim M. Olthoff, Abraham Shaked, Tobias D. Raabe

First published: 26 April 2021 | https://doi.org/10.1002/hep.31857

Researchers at the University of Pennsylvania have generated organoids from NASH patient liver tissue

- Used protocol described in Huch et al., 2015 (HepatiCult<sup>™</sup> is based on this technology)
- NASH patient derived liver organoids are able to replicate in vivo pathology better than mouse models



# Hepatic Organoids for Disease Modelling (NASH)



#### NASH Patient-Derived Liver Organoids Exhibited:

- Decreased albumin
  production
- Increased free fatty acid induced accumulation
- Increased CYP metabolism
- Increased apoptosis

#### Significance:

- NASH organoids recapitulated key diseases characteristics
- Avenues for drug/therapy development

Figures adapted/modified from McCarron et al. Hepatology 2021



## **Hepatic Organoids for Cancer Modelling**

#### Published: 13 November 2017

## Human primary liver cancer-derived organoid cultures for disease modeling and drug screening

Laura Broutier, Gianmarco Mastrogiovanni, Monique MA Verstegen, Hayley E Francies, Lena Morrill Gavarró, Charles R Bradshaw, George E Allen, Robert Arnes-Benito, Olga Sidorova, Marcia P Gaspersz, Nikitas Georgakopoulos, Bon-Kyoung Koo, Sabine Dietmann, Susan E Davies, Raaj K Praseedom, Ruby Lieshout, Jan N M IJzermans, Stephen J Wigmore, Kourosh Saeb-Parsy, Mathew J Garnett, Luc JW van der Laan & Meritxell Huch 🖂

Nature Medicine 23, 1424–1435 (2017) Cite this article



00

hepatocellular carcinoma (HCC) hepatocellular-cholangiocarcinoma (CHC) cholangiocarcinoma samples (CC)







- Liver tumoroids allow the identification of patient-specific drug sensitivities and recapitulate multiple cancer subtypes
- Primary liver cancer (PLC) derived organoids present opportunities for drug testing, therapeutic development and advances in personalized medicine approaches



## Hepatic Organoids for Modelling Disease Progression

www.nature.com/aps

Acta Pharmacologica Sinica

#### ARTICLE

Hepatocyte-derived VEGFA accelerates the progression of non-alcoholic fatty liver disease to hepatocellular carcinoma via activating hepatic stellate cells

Hao Shen<sup>1</sup>, Han Yu<sup>1</sup>, Qian-yu Li<sup>2</sup>, Ya-ting Wei<sup>1,3</sup>, Jing Fu<sup>1</sup>, Hui Dong<sup>1</sup>, Dan Cao<sup>1</sup>, Lin-na Guo<sup>1</sup>, Lei Chen<sup>1</sup>, Yuan Yang<sup>4</sup>, Ying Xu<sup>1</sup>, Meng-chao Wu<sup>1</sup>, Hong-yang Wang<sup>1,3</sup> and Yao Chen<sup>1</sup>

# Normal liver tissue HBV-fibrosis NAFLD-fibrosis NAFLD-fibrosis Normal liver tis

NAFLD-fibrosis well 200 **HBV-fibrosis** Normal liver tissue >200 (µm) 50-200



Shen et al. Acta Pharmacologica 2022



#### Application 3: Gene Editing



## Hepatic Organoids for Gene Editing Applications

Article

#### **Cell Stem Cell**

#### Probing the Tumor Suppressor Function of BAP1 in CRISPR-Engineered Human Liver Organoids

**Graphical Abstract** 

generation of BAP1 mutant

functional

analyses

tumor

Authors

Benedetta Artegiani, Lisa van Voorthuijsen, Rik G.H. Lindeboom, ..., Jacco van Rheenen, Michiel Vermeulen, Hans Clevers

Correspondence h.clevers@hubrecht.eu

#### In Brief

Artegiani et al. show that BAP1 mutation in human liver organoids coincides with loss of multiple epithelial characteristics through impairment of chromatin accessibility and gene expression, and this is critical for the acquisition of malignant features in a human model of cholangiocarcinoma.



tion of chromatin acces

essential role in developing malignant features

and gene expression

Model liver oncogenic process in specific human tissues using organoids and CRISPR-Cas9 gene editing approach.

loss of multiple

enithelial characteristics

Artegiani B et al. Cell Stem Cell 2019





## Hepatic Organoids for Gene Editing Applications

#### Protocol Published: 27 November 2020

#### Establishment of human fetal hepatocyte organoids and CRISPR–Cas9-based gene knockin and knockout in organoid cultures from human liver

Delilah Hendriks 🖾, Benedetta Artegiani 🖾, Huili Hu, Susana Chuva de Sousa Lopes & Hans Clevers 🖾

Nature Protocols 16, 182-217 (2021) Cite this article

#### Gene editing workflow





Characterization of WT and knockin lines





Experimental Workflow for Hepatic Organoid Genome Editing



29



- Hepatic organoids can be used for a wide variety of applications, including hepatotoxicity/cytotoxicity, disease modelling, and gene editing applications
- Hepatic organoids can be derived from patients with inherited liver disease:
  - Opportunity to establish drug-screening platforms and toxicity assays for personalized medicine
  - Understand the signaling pathways involved in the underlying disease
- Genome editing technologies coupled with the hepatic organoid culture platform:
  - Reintroduction of genetically corrected stem cells into affected patients
  - Provides experimental platform for mechanistic studies of cancer gene function in a human context

roperty





