

# ALDEFLUOR™ Kit

For the identification, evaluation, and isolation of stem and progenitor cells expressing high levels of ALDH

Catalog #01700 1 Kit



Scientists Helping Scientists™ | [WWW.STEMCELL.COM](http://WWW.STEMCELL.COM)

TOLL FREE PHONE 1 800 667 0322 • PHONE +1 604 877 0713

[INFO@STEMCELL.COM](mailto:INFO@STEMCELL.COM) • [TECHSUPPORT@STEMCELL.COM](mailto:TECHSUPPORT@STEMCELL.COM)

FOR GLOBAL CONTACT DETAILS VISIT OUR WEBSITE

## Product Description

ALDEFLUOR™ is a reagent kit that is used to identify human cells that express high levels of the enzyme aldehyde dehydrogenase (ALDH). The activated ALDEFLUOR™ Reagent, BODIPY-aminoacetaldehyde (BAAA), is a fluorescent non-toxic substrate for ALDH, which freely diffuses into intact and viable cells. In the presence of ALDH, BAAA is converted into BODIPY-aminoacate (BAA), which is retained inside the cells. The amount of fluorescent reaction product is proportional to the ALDH activity in the cells and is measured using a flow cytometer. Viable ALDH-bright (ALDH<sup>br</sup>) cells can, in principle, be isolated using a cell sorter. Active efflux of the reaction product is inhibited by an efflux inhibitor in the ALDEFLUOR™ Assay Buffer. A specific inhibitor of ALDH, diethylaminobenzaldehyde (DEAB), is used to control for background fluorescence.

ALDEFLUOR™ is optimized for the detection of ALDH<sup>br</sup> hematopoietic cells in human blood and bone marrow, but it can also be used with non-hematopoietic cells. For a full list of ALDEFLUOR™ products, visit [www.stemcell.com](http://www.stemcell.com).

## Product Information

COMPONENT NAME	COMPONENT #	SIZE	STORAGE	SHELF LIFE
Dry ALDEFLUOR™ Reagent <sup>††</sup>	01703	50 µg	Store at 2 - 8°C.	Stable until expiry date (EXP) on label.
ALDEFLUOR™ Diethylaminobenzaldehyde (DEAB) Reagent, 1.5 mM in 95% ethanol <sup>**</sup>	01705	1 mL	Store at 2 - 8°C.	Stable until expiry date (EXP) on label.
Hydrochloric Acid (HCl, 2 N) <sup>**†</sup>	01704	1.5 mL	Store at 2 - 8°C.	Stable until expiry date (EXP) on label.
Dimethylsulphoxide (DMSO) <sup>**†</sup>	01706	1.5 mL	Store at 2 - 8°C.	Stable until expiry date (EXP) on label.
ALDEFLUOR™ Assay Buffer <sup>**</sup>	01702	2 x 55 mL	Store at 2 - 8°C. Do not freeze.	Stable until expiry date (EXP) on label.

\*ALDEFLUOR™ Reagent is not cytotoxic. The combination of Dry ALDEFLUOR™ Reagent, DMSO, and HCl shows no cytotoxic or phototoxic effects at concentrations 100-fold above those used in this assay.

\*\*Please refer to the Safety Data Sheet for hazard information.

<sup>†</sup>This component is sold as part of the ALDEFLUOR™ Kit (Catalog #01700) and is not available for individual sale.

## Directions for Use

Please read the entire protocol before proceeding.

### A) ALDEFLUOR™ ACTIVATION

Dry ALDEFLUOR™ Reagent is provided in a stable, inactive form (BODIPY-aminoacetaldehyde-diethyl acetate, BAAA-DA). For use, Dry ALDEFLUOR™ Reagent is dissolved in DMSO, converted to the fluorescent-activated ALDEFLUOR™ Reagent (BAAA) by treatment with 2 N HCl, and diluted with ALDEFLUOR™ Assay Buffer:

1. Assemble all necessary supplies and allow kit reagents to come to room temperature (15 - 25°C) before use.
2. Add 25 µL of DMSO to the vial of Dry ALDEFLUOR™ Reagent, mix well, and let it stand for 1 minute at room temperature.  
NOTE: Dry ALDEFLUOR™ Reagent is an orange-red powder that changes to a bright yellow-green color upon addition of DMSO.
3. Add 25 µL of 2 N HCl and mix well. Incubate this mixture for 15 minutes at room temperature.  
NOTE: Adding 2 N HCl before DMSO will render the product inactive.

4. Add 360 µL of ALDEFLUOR™ Assay Buffer to the vial and mix.  
NOTE: Upon addition of the ALDEFLUOR™ Assay Buffer, the solution may appear slightly cloudy. This does not affect the assay performance.
5. Keep the activated ALDEFLUOR™ Reagent at 2 - 8°C during use.
6. Aliquot the remaining activated ALDEFLUOR™ Reagent and store at -20°C.

## B) CELL SAMPLE PREPARATION

1. Prepare fresh or previously frozen cell samples according to standard procedures for the cell type.
2. If using blood samples where the red blood cell (RBC) to leukocyte ratio (RBC:WBC) of the specimen is > 2:1, lyse the erythrocytes with Ammonium Chloride Solution (Catalog #07800).
3. After RBC lysis, centrifuge the sample for 5 minutes at 250 x g. Remove the supernatant and suspend cells in 1 mL of ALDEFLUOR™ Assay Buffer.
4. Perform a cell count.
5. If using hematopoietic cells (e.g. peripheral blood, apheresis product, bone marrow, or cord blood), adjust the sample to a concentration of  $1 \times 10^6$  cells/mL with the ALDEFLUOR™ Assay Buffer.  
NOTE: For other cell types, different cell concentrations may be more appropriate. For optimization of ALDEFLUOR™ staining conditions for non-hematopoietic cells, refer to the Technical Bulletin: ALDEFLUOR™ Assay Optimization (Document #29902), available at [www.stemcell.com](http://www.stemcell.com).

## C) ALDEFLUOR™ ASSAY

1. Label one “test” and one “control” tube for each sample to be tested. Place 1 mL of the adjusted cell suspension (section B) into each “test” sample tube.
2. Add 5 µL of ALDEFLUOR™ DEAB Reagent to the “control” tube. Recap control tube and DEAB vial immediately.  
NOTE: ALDEFLUOR™ DEAB is provided in 95% ethanol. Recap immediately to prevent evaporation.
3. Add 5 µL of the activated ALDEFLUOR™ Reagent per milliliter of sample to the first sample “test” tube. Mix and immediately transfer 0.5 mL of the mixture to the DEAB “control” tube.  
NOTE: The ALDH enzymatic reaction begins immediately upon addition of the activated substrate to the cell suspension. It is imperative that an aliquot of the ALDEFLUOR™-reacted cells be added to the DEAB “control” tube without delay.
4. Add control and substrate solutions as described in steps 2 and 3 above for each sample to be tested.
5. Incubate “test” and “control” samples for 30 to 60 minutes (max) at 37°C. Do not exceed 60 minutes.  
NOTE: Optimal incubation times may vary between cell types. For suggestions on optimization of ALDEFLUOR™ staining conditions for non-hematopoietic cells, cultured cells, and cell lines, refer to the Technical Bulletin: ALDEFLUOR™ Assay Optimization (Document #29902), available at [www.stemcell.com](http://www.stemcell.com).
6. Following incubation, centrifuge all tubes for 5 minutes at 250 x g and remove supernatant. Resuspend cell pellets in 0.5 mL of ALDEFLUOR™ Assay Buffer and store the cells on ice or at 2 - 8°C.  
NOTE: If immunophenotyping is to be performed, add and incubate the antibodies after step 6. To prevent efflux of the ALDEFLUOR™ product, it is important that the antibody incubation is performed in ALDEFLUOR™ Assay Buffer. Whenever possible, keep the cells chilled (2 - 8°C or on ice) to slow down the product efflux.
7. OPTIONAL: Perform a viability cell count. If the sample contains fewer than 90% viable cells, it is recommended to stain cells with a DNA dye such as propidium iodide or 7-actinoaminomycin-D in order to stain dead and late apoptotic cells.

## D) FLOW CYTOMETER SET-UP AND DATA ACQUISITION

Refer to the Technical Bulletin: The Basic FACS on ALDEFLUOR™ (Document #28000), available at [www.stemcell.com](http://www.stemcell.com).

## Notes and Tips

- Fresh or previously frozen samples can be analyzed for ALDH<sup>br</sup> cells. However, the ALDEFLUOR™ kit will only detect ALDH activity in cells that are viable and have intact cell membranes.
- Removal of erythrocytes from the sample is required. Erythrocytes may be removed by lysis using reagents that do not contain detergents or fixatives (e.g. Ammonium Chloride Solution, Catalog #07800). They may also be removed by density centrifugation.
- When frozen aliquots of the activated ALDEFLUOR™ Reagent are thawed, a small precipitate (pellet) may be observed. Before use, mix the thawed reagent to resuspend the precipitate. This precipitate does not affect assay performance.
- The cell lines A549 (lung carcinoma), SKBR3 (breast cancer), and K562 (CML) express ALDH activity and can be used as positive controls for the ALDEFLUOR™ assay. In addition, commercially available bone marrow mononuclear cells can also be used as positive controls.

- Identification of rare ALDH<sup>br</sup> cells in heterogeneous cell samples can be improved by removing mature hematopoietic cells on the basis of lineage antigen expression and enriching for ALDH<sup>br</sup> cells by using RosetteSep™ Human Progenitor Cell Basic Pre-Enrichment Kit, RosetteSep™ Human Hematopoietic Progenitor Cell Enrichment Cocktail, or EasySep™ Human CD34 Positive Selection Kit II prior to performing the ALDEFLUOR™ assay.

PRODUCTS ARE FOR RESEARCH USE ONLY AND NOT INTENDED FOR HUMAN OR ANIMAL DIAGNOSTIC OR THERAPEUTIC USES UNLESS OTHERWISE STATED.

Copyright © 2023 by STEMCELL Technologies Inc. All rights reserved including graphics and images. STEMCELL Technologies & Design, STEMCELL Shield Design, Scientists Helping Scientists, ALDEFLUOR, EasySep, and RosetteSep are trademarks of STEMCELL Technologies Inc. All other trademarks are the property of their respective holders. While STEMCELL has made all reasonable efforts to ensure that the information provided by STEMCELL and its suppliers is correct, it makes no warranties or representations as to the accuracy or completeness of such information.