

StemSpan™ Myeloid Expansion Supplement (100X)

Serum-free culture supplement for expansion of human granulocytes

Catalog #02693

1 mL



Scientists Helping Scientists™ | www.stemcell.com

TOLL FREE PHONE 1 800 667 0322 • PHONE +1 604 877 0713

[INFO@STEMCELL.COM](mailto:info@stemcell.com) • [TECHSUPPORT@STEMCELL.COM](mailto:techsupport@stemcell.com)

FOR GLOBAL CONTACT DETAILS VISIT OUR WEBSITE

Product Description

StemSpan™ Myeloid Expansion Supplement (100X) contains a combination of recombinant human cytokines formulated to selectively promote the expansion and differentiation of human myeloid progenitor cells, generating large numbers of CD15+ granulocytes from CD34+ cells isolated from human cord blood (CB), bone marrow (BM), or mobilized peripheral blood (mPB) samples.

StemSpan™ Myeloid Expansion Supplement (100X) is intended for use in combination with any of the following StemSpan™ media:

- StemSpan™ SFEM (Catalog #09600)
- StemSpan™ SFEM II (Catalog #09605)
- StemSpan™-XF (Catalog #100-0073)
- StemSpan™-AOF (Catalog #100-0130)

Advantages:

- Promotes the production of thousands of CD15+ granulocytes in liquid cultures.
- Supplied as a 100X concentrate. After thawing and mixing, the tube contents can be added directly to any hematopoietic cell expansion medium of choice.

Properties

Storage: Store at -20°C.

Shelf Life: Stable until expiry date (EXP) on label.

Contains:

- Recombinant human stem cell factor (SCF)
- Recombinant human thrombopoietin (TPO)
- Recombinant human granulocyte colony-stimulating factor (G-CSF)
- Recombinant human granulocyte/macrophage colony-stimulating factor (GM-CSF)

Directions for Use

PREPARATION OF COMPLETE MEDIUM

1. Thaw StemSpan™ Myeloid Expansion Supplement (100X) at room temperature (15 - 25°C) until just thawed. Mix thoroughly.
NOTE: If necessary, centrifuge for 30 seconds to recover liquid from cap.
NOTE: Once thawed, store supplement at 2 - 8°C for up to 1 month. Alternatively, aliquot and store at -20°C. After thawing aliquots, do not re-freeze.
2. Add StemSpan™ Myeloid Expansion Supplement (100X) to culture medium at a 1 in 100 dilution (e.g. add 1 mL of Supplement to 99 mL of culture medium). Mix thoroughly.

RECOMMENDED PROTOCOL FOR GRANULOCYTE EXPANSION WITH StemSpan™ MEDIA

For optimal performance, use one of the StemSpan™ media suggested below.

1. Prepare StemSpan™ medium of choice by following instructions in the applicable Product Information Sheet, available at www.stemcell.com, or contact us to request a copy.
 - StemSpan™ SFEM II (Document #10000000430)
 - StemSpan™-AOF (Document #10000016677)

NOTE: If using StemSpan™ SFEM II, ensure the medium is completely thawed and mixed thoroughly before use.

2. Prepare complete medium as described in the Preparation section, using a StemSpan™ medium as the culture medium.
3. Thaw cryopreserved CD34+ cells, or use an EasySep™ kit to isolate CD34+ cells from fresh whole CB, BM, mPB, or frozen mononuclear cells (MNCs), as indicated in Table 1.

Alternatively, source frozen isolated CD34+ cells from BM (e.g. Human Bone Marrow CD34+ Cells, Frozen, Catalog #70002*), CB (e.g. Human Cord Blood CD34+ Cells, Frozen, Catalog #70008*), or mPB (e.g. Mobilized Human Peripheral Blood CD34+ Cells, Frozen, Catalog #70060*).

Table 1. Recommended Cell Isolation Kits for Various Cell Sources

CELL SOURCE	RECOMMENDED CELL ISOLATION KIT
Fresh whole CB	EasySep™ Human Cord Blood CD34 Positive Selection Kit II (Catalog #17896)
Fresh BM (e.g. Human Whole Bone Marrow, Fresh, Catalog #70502*)	EasySep™ Human CD34 Positive Selection Kit II (Catalog #17856)
Fresh peripheral blood mobilized with granulocyte colony-stimulating factor (G-CSF), plerixafor, or a combination of both (e.g. Human Mobilized Peripheral Blood Leukopak, G-CSF, Fresh, Catalog #200-0602*)	EasySep™ Human CD34 Positive Selection Kit II (Catalog #100-1569)
Frozen MNCs from BM (e.g. Human Bone Marrow Mononuclear Cells, Frozen, Catalog #70001*), CB (e.g. Human Cord Blood Mononuclear Cells, Frozen, Catalog #70007*), or mPB (e.g. G-CSF Mobilized Human Peripheral Blood Mononuclear Cells, Frozen, Catalog #70049*)	EasySep™ Human CD34 Positive Selection Kit II (Catalog #17856)

* Some primary cell products are available only in select regions. Contact us at techsupport@stemcell.com for further information.

4. **Day 0:** Plate CD34+ cells in complete medium at 1×10^4 to 4×10^4 cells/mL. Refer to Table 2 for recommended plating concentrations; optimal cell concentrations and cultureware are dependent on experimental objectives and cell quality.

Table 2. Recommended CD34+ Cell Concentrations for Various Cultureware

CULTUREWARE**	VOLUME OF MEDIUM PER WELL	NUMBER OF CELLS PER WELL
6-well plate	2 mL	2×10^4
24-well plate	1 mL	1×10^4
96-well plate	100 μ L	1×10^3

** Both tissue culture-treated and non-tissue culture-treated are suitable.

5. Incubate at 37°C and 5% CO₂.
6. **Day 3 or 4:** Add an equal volume of fresh complete medium (i.e. if using a 24-well plate, add 1 mL complete medium per well).
7. **Day 7:** Harvest and replat in fresh complete medium at $< 1 \times 10^5$ cells/mL. Alternatively, add fresh complete medium to maintain cell concentration at $< 1 \times 10^5$ cells/mL.
8. **Day 10/11:** Add an equal volume of fresh complete medium (i.e. if using a 24-well plate, add 1 mL complete medium per well) to maintain cell concentration at $< 1 \times 10^5$ cells/mL.
9. **Day 14:** Harvest cells for evaluation or downstream applications. Count total viable cells using Trypan Blue (e.g. Catalog #07050) and a hemocytometer (e.g. Catalog #100-1181), or an automated cell counting method, and measure granulocyte progenitor cell surface marker expression by flow cytometry.

NOTE: Cultures can be continued beyond 14 days with periodic dilution every 3 - 4 days to maintain a cell concentration of 1×10^5 cells/mL.

Related Products

For related products, including specialized culture and storage media, supplements, antibodies, cytokines, and small molecules, visit www.stemcell.com/HSPCworkflow or contact us at techsupport@stemcell.com. For available fresh and cryopreserved peripheral blood, cord blood, and bone marrow products in your region, visit www.stemcell.com/primarycells.

PRODUCTS ARE FOR RESEARCH USE ONLY AND NOT INTENDED FOR HUMAN OR ANIMAL DIAGNOSTIC OR THERAPEUTIC USES UNLESS OTHERWISE STATED. FOR ADDITIONAL INFORMATION ON QUALITY AT STEMCELL, REFER TO WWW.STEMCELL.COM/COMPLIANCE.

Copyright © 2025 by STEMCELL Technologies Inc. All rights reserved including graphics and images. STEMCELL Technologies & Design, STEMCELL Shield Design, Scientists Helping Scientists, EasySep, and StemSpan are trademarks of STEMCELL Technologies Inc. All other trademarks are the property of their respective holders. While STEMCELL has made all reasonable efforts to ensure that the information provided by STEMCELL and its suppliers is correct, it makes no warranties or representations as to the accuracy or completeness of such information.