

STEMdiff™ Cardiomyocyte Freezing Medium

Medium for cryopreserving hPSC-derived cardiomyocytes

Catalog #05030

50 mL



Scientists Helping Scientists™ | WWW.STEMCELL.COM

TOLL FREE PHONE 1 800 667 0322 • PHONE +1 604 877 0713

INFO@STEMCELL.COM • TECHSUPPORT@STEMCELL.COM

FOR GLOBAL CONTACT DETAILS VISIT OUR WEBSITE

Product Description

STEMdiff™ Cardiomyocyte Freezing Medium can be used to cryopreserve cardiomyocytes derived from human pluripotent stem cells (hPSCs) and maintained in STEMdiff™ Cardiomyocyte Maintenance Medium (Catalog #05020). STEMdiff™ Cardiomyocyte Dissociation Kit (Catalog #05025) is required for dissociating cardiomyocytes prior to freezing. After thawing using STEMdiff™ Cardiomyocyte Support Medium (Catalog #05027), the cardiomyocytes can be used in various downstream applications and analyses. This medium supports hPSC-derived cardiomyocytes generated using STEMdiff™ Ventricular Cardiomyocyte Differentiation Kit (Catalog #05010) or STEMdiff™ Atrial Cardiomyocyte Differentiation Kit (Catalog #100-0215).

NOTE: To reduce stress on hPSC-derived cardiomyocytes during the dissociation and thawing process, use of STEMdiff™ Cardiomyocyte Plating Supplement (Catalog #100-1120) in combination with STEMdiff™ Cardiomyocyte Support Medium is recommended. The addition of the STEMdiff™ Cardiomyocyte Plating Supplement increases the viability and plating efficiency of hPSC-derived cardiomyocytes.

Product Information

PRODUCT NAME	CATALOG #	SIZE	STORAGE	SHELF LIFE
STEMdiff™ Cardiomyocyte Freezing Medium*	05030	50 mL	Store at 2 - 8°C.	Stable for 12 months from date of manufacture (MFG) on label.

*Please refer to the Safety Data Sheet (SDS) for hazard information. This product contains components dissolved in dimethyl sulfoxide (DMSO). DMSO is a strong solvent and skin penetrant, and can transport many substances through the skin. DMSO can also penetrate some protective glove materials including latex and silicone. Extra caution should be utilized when handling this product.

Materials Required But Not Included

PRODUCT NAME	CATALOG #
Corning® Matrigel® hESC-Qualified Matrix	Corning 354277
STEMdiff™ Cardiomyocyte Dissociation Kit <ul style="list-style-type: none">STEMdiff™ Cardiomyocyte Dissociation MediumSTEMdiff™ Cardiomyocyte Support Medium (Catalog #05027)	05025
STEMdiff™ Cardiomyocyte Plating Supplement (100X)	100-1120
STEMdiff™ Cardiomyocyte Maintenance Kit <ul style="list-style-type: none">STEMdiff™ Cardiomyocyte Maintenance Basal MediumSTEMdiff™ Cardiomyocyte Maintenance Supplement (50X)	05020
Trypan Blue	07050

Preparation of Media

For storage, stability, and preparation instructions for STEMdiff™ Cardiomyocyte Dissociation Medium and STEMdiff™ Cardiomyocyte Plating Supplement, refer to the Product Information Sheet (PIS; Document #10000003446).

For storage, stability, and preparation instructions for STEMdiff™ Cardiomyocyte Maintenance Medium, refer to the corresponding PIS (Document #10000009775).

Directions for Use

Please read the entire protocol before proceeding. Use sterile technique when performing the following protocols.

hPSC-derived cardiomyocytes can be harvested and cryopreserved as early as Day 15 of differentiation/maintenance.

A. CRYOPRESERVING hPSC-DERIVED CARDIOMYOCYTES

1. Dissociate/harvest hPSC-derived cardiomyocytes using STEMdiff™ Cardiomyocyte Dissociation Kit. For more information, refer to the applicable PIS (Document #10000003446).
2. Centrifuge cell suspension at 300 x g for 5 minutes at room temperature (15 - 25°C).
3. Gently pour off the supernatant, taking care not to disrupt the cell pellet.
4. Gently resuspend the cells in STEMdiff™ Cardiomyocyte Freezing Medium to a final concentration of 5×10^5 cells/mL.
5. Transfer 1 mL of cell suspension into a labeled cryovial.
6. Freeze cell suspension using a standard slow rate-controlled cooling protocol that reduces temperatures at approximately -1°C/minute, followed by long-term storage at -135°C (liquid nitrogen) or colder. Long-term storage at -80°C is not recommended.

B. THAWING hPSC-DERIVED CARDIOMYOCYTES

Cryopreserved hPSC-derived cardiomyocytes should be thawed and plated onto Corning® Matrigel®-coated cultureware. For coating plates with Corning® Matrigel®, refer to the Technical Manual for mTeSR™1, mTeSR™ Plus, or TeSR™-E8™, available at www.stemcell.com, or contact us to request a copy.

The following instructions are for a 12-well tissue culture plate. For other cultureware, adjust accordingly.

1. Coat a 12-well tissue culture plate with Corning® Matrigel® hESC-Qualified Matrix and bring to room temperature for at least 1 hour prior to use.
2. Prepare STEMdiff™ Cardiomyocyte Plating Medium (see Preparation of Media) and warm to room temperature.
3. Thaw hPSC-derived cardiomyocytes in a 37°C water bath by gently shaking the cryovial continuously until only a small frozen cell pellet remains.
4. Add 5 - 7 mL of STEMdiff™ Cardiomyocyte Plating Medium to a 15 mL conical tube (e.g. Catalog #38009).
5. Using a 2 mL pipette, gently transfer the contents of the cryovial to the tube from step 4.
6. Centrifuge the cells at 300 x g for 5 minutes at room temperature.
7. Aspirate the supernatant and gently add 1 - 2 mL of STEMdiff™ Cardiomyocyte Plating Medium to resuspend cells.
8. Perform a cell count using an automated cell counter (e.g. NucleoCounter® NC-250™) or with Trypan Blue and a Hausser Scientific™ Bright-Line Hemocytometer (Catalog #100-1181).
9. Aspirate Corning® Matrigel® from the tissue culture plate prepared in step 1. Add 2 mL of STEMdiff™ Cardiomyocyte Plating Medium per well.
10. Add cells at a density appropriate for downstream assays or other applications. Incubate at 37°C for 24 hours.
11. Warm STEMdiff™ Cardiomyocyte Maintenance Medium to room temperature.
12. Remove STEMdiff™ Cardiomyocyte Plating Medium and add 2 mL of STEMdiff™ Cardiomyocyte Maintenance Medium per well. Incubate at 37°C.
13. Every 2 days, perform a full-medium change with 2 mL of STEMdiff™ Cardiomyocyte Maintenance Medium per well.

PRODUCTS ARE FOR RESEARCH USE ONLY AND NOT INTENDED FOR HUMAN OR ANIMAL DIAGNOSTIC OR THERAPEUTIC USES UNLESS OTHERWISE STATED. FOR ADDITIONAL INFORMATION ON QUALITY AT STEMCELL, REFER TO WWW.STEMCELL.COM/COMPLIANCE.

Copyright © 2024 by STEMCELL Technologies Inc. All rights reserved including graphics and images. STEMCELL Technologies & Design, STEMCELL Shield Design, Scientists Helping Scientists, and STEMdiff are trademarks of STEMCELL Technologies Canada Inc. mTeSR, TeSR, and E8 are trademarks of WARF. Corning and Matrigel are registered trademarks of Corning Incorporated. All other trademarks are the property of their respective holders. While STEMCELL has made all reasonable efforts to ensure that the information provided by STEMCELL and its suppliers is correct, it makes no warranties or representations as to the accuracy or completeness of such information.