

NeuroCult™ SM1 Neuronal Supplement

Supplement (50X) for the serum-free culture of neurons



Scientists Helping Scientists™ | WWW.STEMCELL.COM

TOLL FREE PHONE 1 800 667 0322 • PHONE +1 604 877 0713

INFO@STEMCELL.COM • TECHSUPPORT@STEMCELL.COM

FOR GLOBAL CONTACT DETAILS VISIT OUR WEBSITE

Catalog #05711 10 mL
#100-1281 100 mL

Product Description

NeuroCult™ SM1 (STEMCELL Modified-1) Neuronal Supplement is a serum-free supplement for the culture of primary or pluripotent stem cell (PSC)-derived neurons. NeuroCult™ SM1 is used in combination with a basal neural medium such as BrainPhys™ Neuronal Medium (Catalog #05790) or as a serum-replacement supplement for various customizable applications wherever B27 is used.

The NeuroCult™ SM1 formulation was developed based on the published formulation of Brewer's B27 supplement. It has been optimized to reproducibly support increased survival and maturation of functional neurons in both short- and long-term cultures with minimal glial cell contamination. NeuroCult™ SM1 is also a component of several BrainPhys™ kits for primary and PSC-derived neuronal cell culture. For further details and protocols, refer to the Product Information Sheet (PIS) for BrainPhys™, available at www.stemcell.com, or contact us to request a copy.

Ordering Information

PRODUCT NAME	CATALOG #	SIZE	KIT COMPONENTS
NeuroCult™ SM1 Neuronal Supplement	05711	10 mL	Not applicable
NeuroCult™ SM1 Neuronal Supplement	100-1281	100 mL	Not applicable
BrainPhys™ Neuronal Medium and SM1 Kit	05792	1 Kit	<ul style="list-style-type: none">BrainPhys™ Neuronal MediumNeuroCult™ SM1 Neuronal Supplement
BrainPhys™ Neuronal Medium N2-A & SM1 Kit	05793	1 Kit	<ul style="list-style-type: none">BrainPhys™ Neuronal MediumNeuroCult™ SM1 Neuronal SupplementN2 Supplement-A
BrainPhys™ Primary Neuron Kit	05794	1 Kit	<ul style="list-style-type: none">BrainPhys™ Neuronal MediumNeuroCult™ SM1 Neuronal SupplementNeuroCult™ Neuronal Plating Medium

Storage and Stability

PRODUCT NAME	CATALOG #	SIZE	STORAGE	SHELF LIFE
NeuroCult™ SM1 Neuronal Supplement*	05711	10 mL	Store at -20°C.	Stable until expiry date (EXP) on label.
NeuroCult™ SM1 Neuronal Supplement*	100-1281	100 mL	Store at -20°C.	Stable until expiry date (EXP) on label.

*Please refer to the Safety Data Sheet (SDS) for hazard information. Lot-to-lot color variations include light to dark yellow or orange. This will not affect performance.

Directions for Use

NOTE: NeuroCult™ SM1 is provided as a 50X stock solution that may be added to the medium of your choice; it may be used as an alternative to Brewer's B27 supplement.

Thaw one bottle of NeuroCult™ SM1 at room temperature (15 - 25°C) or overnight at 2 - 8°C.

NOTE: If not used immediately, aliquot and store at -20°C. Do not exceed the expiry date (EXP) as indicated on the label. After thawing aliquots, use immediately. Do not re-freeze.

- For culturing primary tissue-derived neurons and for neuronal differentiation of hPSC-derived neural progenitor cells using NeuroCult™ SM1, refer to the PIS for BrainPhys™ Neuronal Medium, available at www.stemcell.com, or contact us to request a copy.
- For detailed procedures on immunolabeling, refer to the Technical Manual: In Vitro Proliferation and Differentiation of Mouse and Rat Neural Stem and Progenitor Cells Using NeuroCult™, available at www.stemcell.com, or contact us to request a copy.

Related Products

For related products, including specialized cell culture and storage media, supplements, antibodies, and small molecules, visit www.stemcell.com, or contact us at techsupport@stemcell.com.

References

Brewer GJ et al. (1993) Optimized survival of hippocampal neurons in B27-supplemented Neurobasal, a new serum-free medium combination. *J Neurosci Res* 35(5): 567–76.

PRODUCTS ARE FOR RESEARCH USE ONLY AND NOT INTENDED FOR HUMAN OR ANIMAL DIAGNOSTIC OR THERAPEUTIC USES UNLESS OTHERWISE STATED. FOR ADDITIONAL INFORMATION ON QUALITY AT STEMCELL, REFER TO WWW.STEMCELL.COM/COMPLIANCE.

Copyright © 2023 by STEMCELL Technologies Inc. All rights reserved including graphics and images. STEMCELL Technologies & Design, STEMCELL Shield Design, Scientists Helping Scientists, and NeuroCult are trademarks of STEMCELL Technologies Canada Inc. BrainPhys is a registered trademark of the Salk Institute for Biological Studies, used under exclusive license. All other trademarks are the property of their respective holders. While STEMCELL has made all reasonable efforts to ensure that the information provided by STEMCELL and its suppliers is correct, it makes no warranties or representations as to the accuracy or completeness of such information.