

# Spleen Dissociation Medium

Medium for dissociation of mouse spleen

Catalog #07915

10 x 4 mL

## Product Description

This product is designed to maximize the viability and recovery of dendritic cells from mouse spleen when combined with EasySep™ cell separation technology. This medium contains collagenase, DNase, and fetal bovine serum (FBS), and has been optimized for maximum viability of isolated spleen dendritic cells. A single 4 mL tube is sufficient for processing up to two mouse spleens.

## Materials Required but Not Included

- 60 mm Tissue Culture-Treated Dishes (Catalog #200-0623)
- Blunt-End Needles, 16 Gauge (Catalog #28110)
- 3 cc Syringes (Catalog #28230)
- 70 µm Reversible Strainers (e.g. Catalog #27216/27260)

## Properties

**Stability and Storage:** Store at -20°C. Stable until expiry date (EXP) on label.

**Contains:**

- Collagenase IV
- DNase
- Fetal bovine serum (FBS)
- Roswell Park Memorial Institute (RPMI) medium

## Directions for Use

Please refer to the appropriate EasySep™ Product Information Sheet (PIS) for recommended medium and cell resuspension concentration.

### DISSOCIATION AT ROOM TEMPERATURE (15 - 25°C)

1. Thaw individual tubes of Spleen Dissociation Medium at room temperature (15 - 25°C) and use immediately. Do not refreeze.
2. In a 60 mm Tissue Culture-Treated Dish, mince 1 - 2 freshly isolated spleens into a homogeneous paste using dissection scissors and forceps. Spleen fragments should be less than 1 mm in size.
3. Pour the contents of a 4 mL tube of Spleen Dissociation Medium into the dish and mix well. Using a 5 mL serological pipette (e.g. Catalog #38003), return all suspended spleen fragments and Spleen Dissociation Medium to the original tube.
4. Incubate the tube at room temperature for 30 minutes.  
NOTE: For best results, place the tube horizontally on a rocking platform with continuous agitation. Alternatively, gentle agitation every 5 minutes during the incubation is acceptable.
5. If performing downstream DNase treatment (see kit-specific PIS for details), skip this step and continue to step 6. Otherwise, add EDTA to a final concentration of 10 mM (e.g. 80 µL of a 0.5 M stock), mix, and incubate the dish at room temperature for 5 minutes.

## Spleen Dissociation Medium

6. Dissociate spleen fragments into a smooth suspension by gently passing several times through a 16 Gauge Blunt-End Needle attached to a 3 cc Syringe.
7. Pour the entire suspension through a primed 70  $\mu$ m Reversible Strainer into a 50 mL conical tube (e.g. Catalog #38010).  
NOTE: To prime, pass 5 mL of recommended medium through the Reversible Strainer.
8. Rinse the empty Spleen Dissociation Medium tube and Reversible Strainer with an additional 10 mL of recommended medium and add to the 50 mL conical tube.
9. Centrifuge the cell suspension at 300 x g for 10 minutes.
10. Discard supernatant and resuspend cells in appropriate amount of recommended medium. The cells are now ready for downstream applications.

### DISSOCIATION AT 37°C

IMPORTANT NOTE: This protocol has been optimized for use with certain EasySep™ cell separation kits. When 37°C spleen digestion is recommended in the kit-specific PIS, follow the protocol below.

1. Thaw individual tubes of Spleen Dissociation Medium at 37°C and use immediately. Do not refreeze.
2. In a 60 mm Tissue Culture-Treated Dish, mince 1 - 2 freshly isolated spleens into a homogeneous paste using dissection scissors and forceps. Spleen fragments should be less than 1 mm in size.
3. Pour the contents of a 4 mL tube of Spleen Dissociation Medium into the dish and mix well.
4. Incubate the dish at 37°C for 30 minutes.
5. If performing downstream DNase treatment (see kit-specific PIS for details), skip this step and continue to step 6. Otherwise, add EDTA to a final concentration of 10 mM (e.g. 80  $\mu$ L of a 0.5 M stock), mix, and incubate the dish at room temperature (15 - 25°C) for 5 minutes.
6. Dissociate spleen fragments into a smooth suspension by gently passing several times through a 16 Gauge Blunt-End Needle attached to a 3 cc Syringe.
7. Pour the entire suspension through a primed 70  $\mu$ m Reversible Strainer into a 50 mL conical tube (e.g. Catalog #38010).  
NOTE: To prime, pass 5 mL of recommended medium through the Reversible Strainer.
8. Rinse the empty dish and Reversible Strainer with an additional 10 mL of recommended medium and add to the 50 mL conical tube.
9. Centrifuge the cell suspension at 300 x g for 10 minutes.
10. Discard supernatant and resuspend cells in appropriate amount of recommended medium (see appropriate EasySep™ PIS). The cells are now ready for downstream applications.

## Related Products

For related products, including specialized cell culture and storage media, matrices, antibodies, cytokines, and small molecules, visit [www.stemcell.com/hPSCworkflow](http://www.stemcell.com/hPSCworkflow) or contact us at [techsupport@stemcell.com](mailto:techsupport@stemcell.com).

PRODUCTS ARE FOR RESEARCH USE ONLY AND NOT INTENDED FOR HUMAN OR ANIMAL DIAGNOSTIC OR THERAPEUTIC USES UNLESS OTHERWISE STATED. FOR ADDITIONAL INFORMATION ON QUALITY AT STEMCELL, REFER TO [WWW.STEMCELL.COM/COMPLIANCE](http://WWW.STEMCELL.COM/COMPLIANCE).

Copyright © 2024 by STEMCELL Technologies Inc. All rights reserved including graphics and images. STEMCELL Technologies & Design, STEMCELL Shield Design, Scientists Helping Scientists, and EasySep are trademarks of STEMCELL Technologies Canada Inc. All other trademarks are the property of their respective holders. While STEMCELL has made all reasonable efforts to ensure that the information provided by STEMCELL and its suppliers is correct, it makes no warranties or representations as to the accuracy or completeness of such information.