

CryoStor® CS10

Animal component-free, defined cryopreservation medium with 10% DMSO

| | | |
|-----------|-------|-------------|
| Catalog # | 07959 | 5 x 10 mL |
| | 07952 | 16 x 10 mL |
| | 07931 | 5 x 16 mL |
| | 07930 | 100 mL |
| | 07955 | 100 mL Bag |
| | 07940 | 1000 mL Bag |



Scientists Helping Scientists™ | WWW.STEMCELL.COM

TOLL FREE PHONE 1 800 667 0322 • PHONE +1 604 877 0713

INFO@STEMCELL.COM • TECHSUPPORT@STEMCELL.COM

FOR GLOBAL CONTACT DETAILS VISIT OUR WEBSITE

Product Description

CryoStor® CS10 is a uniquely formulated, serum-free, animal component-free, and defined cryopreservation medium containing 10% dimethyl sulfoxide (DMSO). Designed to preserve cells in low-temperature environments (-80°C to -196°C), CryoStor® CS10 provides a safe, protective environment for cells and tissues during the freezing, storage, and thawing processes. CryoStor® CS10 is recommended for the cryopreservation of hepatocytes, tissue samples, human peripheral blood, CHO cells, myeloma cell lines, hybridomas, human mesenchymal stem cells, human embryonic stem (ES) and induced pluripotent stem (iPS) cells, and other extremely sensitive cell types.

- Ready-to-use
- Serum-free, protein-free
- Animal component-free
- cGMP manufactured with USP grade/highest-quality components
- FDA master file
- Sterility, endotoxin, and cell-based quality control testing

Product Information

| CATALOG # | SIZE | STORAGE | SHELF LIFE | CONTAINS |
|-----------|-------------|-------------------|---|----------|
| 07959 | 5 x 10 mL | Store at 2 - 8°C. | Stable until expiry date (EXP) on label. Protect from prolonged exposure to light. | 10% DMSO |
| 07952 | 16 x 10 mL | Store at 2 - 8°C. | Stable until expiry date (EXP) on label. Protect from prolonged exposure to light. | 10% DMSO |
| 07931 | 5 x 16 mL | Store at 2 - 8°C. | Stable until expiry date (EXP) on label. Protect from prolonged exposure to light. | 10% DMSO |
| 07930 | 100 mL | Store at 2 - 8°C. | Stable until expiry date (EXP) on label. Protect from prolonged exposure to light. | 10% DMSO |
| 07955 | 100 mL Bag | Store at 2 - 8°C. | Stable for 2 years from date of manufacture (MFG) on label. Protect from prolonged exposure to light. | 10% DMSO |
| 07940 | 1000 mL Bag | Store at 2 - 8°C. | Stable for 2 years from date of manufacture (MFG) on label. Protect from prolonged exposure to light. | 10% DMSO |

Please refer to the Safety Data Sheet (SDS) for hazard information.

Product may be shipped at room temperature (15 - 25°C) and should be refrigerated upon receipt.

Handling / Directions for Use

CRYOPRESERVING CELLS

For cryopreserving human embryonic stem (ES) cells and induced pluripotent stem (iPS) cells, refer to the Technical Manuals for mTeSR™1 or mTeSR™ Plus for further information, available at www.stemcell.com or contact us to request a copy.

1. Wipe down the outside of the CryoStor® CS10 container with 70% ethanol or isopropanol before opening.
2. Obtain a cell suspension using a cell-specific protocol and centrifuge cells to obtain a cell pellet.
3. Carefully remove the supernatant with a pipette, leaving a small amount of medium to ensure the cell pellet is not disturbed. Resuspend the cell pellet by gently flicking the tube.
4. Add cold (2 - 8°C) CryoStor® CS10, mix thoroughly, and transfer the suspension to a cryovial.
5. Incubate cells at 2 - 8°C for 10 minutes.
6. Cryopreserve cells using a standard slow rate-controlled cooling protocol (approximately -1°C/minute) or an isopropanol freezing container and store at liquid nitrogen temperature (-135°C).

NOTE: Long-term storage at -80°C is not recommended.

THAWING CELLS

1. Warm medium of choice in a 37°C water bath.
2. Wipe the outside of the vial of cells with 70% ethanol or isopropanol.
3. In a biosafety cabinet, twist the cap a quarter-turn to relieve internal pressure and then retighten.
4. Quickly thaw cells in a 37°C water bath by gently shaking the vial. Do not submerge the vial. Remove the vial when only a small frozen cell pellet remains. Do not vortex cells.
5. Wipe the outside of the vial with 70% ethanol or isopropanol.
6. Dilute cells 1 in 10 with warmed medium.
7. Centrifuge the cell suspension at 300 x *g* for 10 minutes at room temperature (15 - 25°C).
8. Carefully remove the supernatant with a pipette, leaving a small amount of medium to ensure the cell pellet is not disturbed. Resuspend the cell pellet by gently flicking the tube.
9. Gently add medium to the tube.

CRYOSTOR PRODUCTS MEET USP <71> STERILITY AND USP <85> ENDOTOXIN TESTING STANDARDS, AND ARE MANUFACTURED UNDER CGMP.

Copyright © 2020 by STEMCELL Technologies Inc. All rights reserved including graphics and images. STEMCELL Technologies & Design, STEMCELL Shield Design, and Scientists Helping Scientists are trademarks of STEMCELL Technologies Canada Inc. CryoStor is a registered trademark of BioLife Solutions. mTeSR is a trademark of WARF. All other trademarks are the property of their respective holders. While STEMCELL has made all reasonable efforts to ensure that the information provided by STEMCELL and its suppliers is correct, it makes no warranties or representations as to the accuracy or completeness of such information.