

ArciTect™ Human HPRT Positive Control Kit

Positive control for CRISPR-Cas9 genome editing

Catalog # 76013 1 Kit



Scientists Helping Scientists™ | WWW.STEMCELL.COM

TOLL FREE PHONE 1 800 667 0322 • PHONE +1 604 877 0713

INFO@STEMCELL.COM • TECHSUPPORT@STEMCELL.COM

FOR GLOBAL CONTACT DETAILS VISIT OUR WEBSITE

Product Description

ArciTect™ Human HPRT Positive Control Kit is designed as a positive control for experiments using the ArciTect™ CRISPR-Cas9 genome editing system. The kit comprises ArciTect™ Human HPRT crRNA (2 nmol), ArciTect™ Human HPRT Primer Mix (2 nmol), ArciTect™ tracrRNA (5 nmol), and ArciTect™ Annealing Buffer (5X). ArciTect™ Human HPRT crRNA first requires annealing to ArciTect™ tracrRNA (included in the kit). This must then be combined with an ArciTect™ Cas9 Nuclease (e.g. Catalog #76002) to form a ribonucleoprotein complex. ArciTect™ Human HPRT Primer Mix can be used to amplify genomic DNA isolated from a population of transfected cells, which can subsequently be used in a T7 endonuclease I assay to determine cleavage efficiency.

This kit has been tested and validated for use with the ArciTect™ line of genome editing products. HPRT, or hypoxanthine phosphoribosyltransferase, is a housekeeping gene and a commonly used control. The kit can be used to optimize transfection protocols and act as a positive control that can be used alongside custom ArciTect™ crRNAs (e.g. Catalog #76010).

Product Information

The following components are sold as a kit and are not available for individual sale.

COMPONENT NAME	COMPONENT #	SIZE	STORAGE	SHELF LIFE
ArciTect™ Human HPRT crRNA	76014	2 nmol	Store at -80°C. Alternatively, store at -20°C for up to 6 months.	Stable for 2 years from date of manufacture (MFG) on label.
ArciTect™ Human HPRT Primer Mix	76015	2 nmol	Store at -20°C.	Stable for 2 years from date of manufacture (MFG) on label.
ArciTect™ tracrRNA	76016A	5 nmol	Store at -80°C. Alternatively, store at -20°C for up to 6 months.	Stable for 2 years from date of manufacture (MFG) on label.
ArciTect™ Annealing Buffer (5X)	76019	100 µL	Store at -20°C. Alternatively, store at 2 - 8°C for up to 6 months.	Stable for 2 years from date of manufacture (MFG) on label.

Materials Required But Not Included

PRODUCT NAME	CATALOG #
ArciTect™ Cas9 Nuclease (or variant)	76002 (or 76009)
Genomic DNA purification kit	e.g. Norgen Biotek 24700
PCR tubes	e.g. 38091
Nuclease-Free Water	79001
ArciTect™ High-Fidelity DNA Polymerase Kit <ul style="list-style-type: none">ArciTect™ High-Fidelity DNA PolymeraseArciTect™ High-Fidelity BufferArciTect™ High GC Content BufferdNTP Mix (10 mM)	76026
PCR purification kit	e.g. QIAGEN 28104
Proteinase K Solution	79016
DNA Loading Dye	79018
1 kb DNA Ladder	79017

Directions for Use

A. PREPARATION OF ArciTect™ HUMAN HPRT PRIMER MIX

1. Briefly centrifuge the vial of ArciTect™ Human HPRT Primer Mix before opening.
2. Add 20 μL of nuclease-free water to prepare a 100 μM stock solution. Mix thoroughly.
NOTE: If not used immediately, aliquot and store at -20°C for up to 6 months.
3. Prepare a 10 μM working solution by diluting the 100 μM stock solution 1 in 10. Mix thoroughly.

B. PREPARATION OF ArciTect™ HUMAN HPRT crRNA AND ArciTect™ tracrRNA STOCK SOLUTIONS

1. Prepare stock solutions of ArciTect™ Human HPRT crRNA and ArciTect™ tracrRNA by centrifuging the vials before opening and adding nuclease-free water to give a final concentration of 200 μM , as indicated in Table 1.

Table 1. Preparation of crRNA and tracrRNA Stock Solutions

COMPONENT	AMOUNT OF RNA	VOLUME OF NUCLEASE-FREE WATER
ArciTect™ Human HPRT crRNA	2 mol	10 μL
ArciTect™ tracrRNA	5 nmol	25 μL

2. Mix thoroughly. If not used immediately, aliquot and store at -80°C for up to 6 months. After thawing the aliquots, use immediately. Do not re-freeze.

C. GENOME EDITING OF CELLS WITH HPRT POSITIVE CONTROL

For further information, refer to the Technical Bulletin: Genome Editing of Human Pluripotent Stem Cells (Document #27084) or Genome Editing of Human Primary T Cells (Document #27155), available at www.stemcell.com or contact us to request a copy.

1. Prepare 80 μM gRNA by combining components in a microcentrifuge tube as indicated in Table 2. Mix thoroughly.

Table 2. Preparation of 80 μM * gRNA

COMPONENT	VOLUME
ArciTect™ Human HPRT crRNA	4 μL
ArciTect™ tracrRNA	4 μL
ArciTect™ Annealing Buffer (5X)	2 μL
Total	10 μL

*80 μM is equivalent to 80 pmol/ μL

2. In a thermocycler or heating block, incubate gRNA mixture at 95°C for 5 minutes followed by 60°C for 1 minute. Cool to room temperature (15 - 25°C) and place on ice.
NOTE: If not used immediately, store at -80°C for up to 6 months.
3. Prepare RNP Complex Mix by combining ArciTect™ Cas9 Nuclease (4 $\mu\text{g}/\mu\text{L}$ or 25 μM) with 80 μM gRNA (prepared in steps 1 - 2) in an appropriate transfection buffer.
NOTE: For electroporation reactions, we recommend 1 - 4 μM Cas9 (final concentration in electrolytic buffer), and for chemical transfection reactions we recommend 10 - 100 nM Cas9 (final concentration in plating media). For both electroporation and chemical transfection methods, a 1:2 - 1:8 molar ratio of Cas9:gRNA is recommended. RNP complex formation must be optimized for cell type and transfection method.
4. Incubate RNP Complex Mix at room temperature (15 - 25°C) for 10 - 20 minutes.
5. Transfect cells using your preferred method.
6. Culture cells for 48 - 72 hours after transfection to allow genome editing to occur.

D. PCR AMPLIFICATION OF GENOMIC DNA FROM EDITED CELLS

1. Isolate genomic DNA (gDNA) from edited cells using a genomic DNA purification kit.
2. Prepare Reagent Mix for PCR amplification of target region from 100 ng of gDNA as indicated in Table 3.

NOTE: Indicated reaction volumes are for ArciTect™ High-Fidelity DNA Polymerase Kit. For other DNA polymerases, adjust component concentrations as required.

Table 3. Reagent Mix for PCR Amplification of Target Region

COMPONENT	VOLUME (μL)	FINAL AMOUNT/ CONCENTRATION
ArciTect™ High GC Content Buffer	10	1X
dNTP Mix (10 mM)	1	200 μM each
10 μM ArciTect™ Human HPRT Primer Mix (working solution)	2.5 μL	0.5 μM
DNA template	Variable	100 ng
ArciTect™ High-Fidelity DNA Polymerase	0.5	1 U
Nuclease-free water	Variable	Bring solution to total volume of 50 μL

3. Amplify the target region by PCR, using the conditions indicated in Table 4.

Table 4. PCR Cycling Conditions for Amplification of Target Region

STEP	TEMPERATURE	TIME
Initial denaturation	98°C	30 seconds
Denaturation, annealing, extension for 35 cycles	98°C	10 seconds
	67°C (annealing)	15 seconds
	72°C	45 seconds
Final extension	72°C	5 minutes
Hold	4°C	Up to 24 hours

4. Extract PCR product using a PCR purification kit, then measure the concentration using a microvolume spectrophotometer.
5. Proceed with the T7 endonuclease I assay using ArciTect™ T7 Endonuclease I Kit (Catalog #76021), as described in the corresponding Product Information Sheet, available at www.stemcell.com or contact us to request a copy.

Related Products

For related products, including other genome editing tools, specialized cell culture and storage media, supplements, antibodies, cytokines, and small molecules, visit www.stemcell.com or contact us at techsupport@stemcell.com.

STEMCELL TECHNOLOGIES INC.'S QUALITY MANAGEMENT SYSTEM IS CERTIFIED TO ISO 13485. PRODUCTS ARE FOR RESEARCH USE ONLY AND NOT INTENDED FOR HUMAN OR ANIMAL DIAGNOSTIC OR THERAPEUTIC USES UNLESS OTHERWISE STATED.

Copyright © 2019 by STEMCELL Technologies Inc. All rights reserved including graphics and images. STEMCELL Technologies & Design, STEMCELL Shield Design, Scientists Helping Scientists, and ArciTect are trademarks of STEMCELL Technologies Canada Inc. All other trademarks are the property of their respective holders. The purchase of the ArciTect™ products conveys to the purchaser the limited, non-transferable right to use, in accordance with all applicable laws and regulations, the ArciTect™ products and any related material solely for research purposes only, not for any clinical, human, agricultural, veterinary, livestock, or commercial purpose, with no warranty (express or implied) from the owners or assignees of the Patents or STEMCELL Technologies each of whom expressly disclaim any warranty regarding results obtained through the use of the ArciTect™ products, and without any exception the owners or assignees of the Patents or STEMCELL Technologies, or as applicable a director, trustee, officer, employee, agent, faculty, official investigator or student thereof, will not suffer or be exposed to any liability, damages, loss, or expense of any kind whatsoever arising from or related to the use of the ArciTect™ products by a purchaser of same. Distribution of ArciTect™ products by STEMCELL Technologies is covered under at least US 8,697,359, US 8,771,945, US 8,795,965, US 8,865,406, US 8,871,445, US 8,889,356, US 8,889,418, US 8,895,308, US 8,906,616 and foreign equivalents ("Patents"). While STEMCELL has made all reasonable efforts to ensure that the information provided by STEMCELL and its suppliers is correct, it makes no warranties or representations as to the accuracy or completeness of such information.