

# StemSpan™-ACF Without Phenol Red

**cGMP, animal component-free medium for culture and expansion of human hematopoietic cells**

Catalog #100-0130

500 mL



Scientists Helping Scientists™ | [WWW.STEMCELL.COM](http://WWW.STEMCELL.COM)

TOLL FREE PHONE 1 800 667 0322 • PHONE +1 604 877 0713

[INFO@STEMCELL.COM](mailto:INFO@STEMCELL.COM) • [TECHSUPPORT@STEMCELL.COM](mailto:TECHSUPPORT@STEMCELL.COM)

FOR GLOBAL CONTACT DETAILS VISIT OUR WEBSITE

## Product Description

StemSpan™-ACF Without Phenol Red is an animal component-free (ACF) medium that has been developed for the in vitro culture and expansion of human hematopoietic cells, when appropriate growth factors and supplements are added. This allows users the flexibility to prepare medium that meets their requirements. StemSpan™-ACF Without Phenol Red contains only recombinant proteins and synthetic components, and does not contain serum or other human- or animal-derived components.

Using appropriate cytokines, StemSpan™-ACF Without Phenol Red may be used to expand CD34+ cells isolated from human cord blood, mobilized peripheral blood, or bone marrow samples, or to expand and differentiate lineage-committed progenitor cells to generate populations of myeloid or megakaryocyte progenitor cells.

StemSpan™-ACF Without Phenol Red is manufactured and tested following relevant cGMPs under a certified quality management system. For additional quality information, refer to [www.stemcell.com/compliance](http://www.stemcell.com/compliance).

## Properties

**Storage:** Store at 2 - 8°C. Do not freeze.

**Shelf Life:** Stable for 12 months from expiry date (EXP) on label.

**Contains:** This product contains only recombinant human proteins and synthetic components.

## Directions for Use

### PREPARATION OF COMPLETE MEDIUM

1. Add desired cytokines, growth factors and, optionally, other supplements to StemSpan™-ACF Without Phenol Red. Mix thoroughly.  
NOTE: Added supplements and cells in sterile cell culture medium should not exceed ~10% of total volume.
2. Add cells, mix thoroughly, and set up cultures as desired.

For a complete list of expansion supplements for use with StemSpan™-ACF Without Phenol Red, see Notes and Tips.

### ASSESSMENT OF HEMATOPOIETIC CELLS

Assessment of CD34+ cells before and after culture may be performed by flow cytometry using the following fluorochrome-conjugated antibody clones:

- Anti-Human CD34 Antibody, Clone 581 (Catalog #60013) or Clone 563 (Catalog #60119) or Clone 8G12 (Catalog #60121)
- Anti-Human CD45 Antibody, Clone HI30 (Catalog #60018) or Clone 2D1 (Catalog #60123)
- Anti-Human CD38 Antibody, Clone AT-1 (Catalog #60131) or Clone HIT2 (Catalog #60014)
- Anti-Human CD90 Antibody, Clone 5E10 (Catalog #60045)

## Notes and Tips

Selection of an optimal growth factor combination is dependent upon the source and type of cells and the experimental objectives of the researcher. StemSpan™ expansion supplements, described below, are recommended for use with StemSpan™-ACF Without Phenol Red. StemSpan™ expansion supplements are manufactured and tested under a quality management system certified to ISO 9001. Refer to the Product Information Sheet (PIS) for each supplement for recommended cell expansion protocols with StemSpan™ media.

- StemSpan™ CD34+ Expansion Supplement (10X) (Catalog #02691)
  - Culture and expansion of large numbers of human CD34+ progenitor cells
  - Contains: rh SCF, rh TPO, rh IL-3, rh IL-6, rh Flt3 ligand, other additives

- StemSpan™ Megakaryocyte Expansion Supplement (100X) (Catalog #02696)
  - Expansion and lineage-specific differentiation of human CD34+ cells into megakaryocyte progenitor cells
  - Contains: rh SCF, rh TPO, rh IL-6, rh IL-9
- StemSpan™ Myeloid Expansion Supplement (100X) (Catalog #02693)
  - Expansion and lineage-specific differentiation of human CD34+ cells into granulocytes
  - Contains: rh SCF, rh TPO, rh G-CSF, rh GM-CSF
- StemSpan™ Myeloid Expansion Supplement II (100X) (Catalog #02694)
  - Expansion and lineage-specific differentiation of human CD34+ cells into monocytes
  - Contains: rh Flt3 ligand, rh SCF, rh TPO, rh M-CSF, rh GM-CSF, other additives

SCF = stem cell factor; EPO = erythropoietin; TPO = thrombopoietin; rh = recombinant human; IL = interleukin; Flt = fms-like tyrosine kinase

#### RELATED PRODUCTS

For related products, including specialized culture and storage media, supplements, antibodies, cytokines, and small molecules, visit [www.stemcell.com/HSPCworkflow](http://www.stemcell.com/HSPCworkflow) or contact us at [techsupport@stemcell.com](mailto:techsupport@stemcell.com). For available fresh and cryopreserved peripheral blood, cord blood, and bone marrow products in your region, visit [www.stemcell.com/primarycells](http://www.stemcell.com/primarycells).

THIS PRODUCT IS MANUFACTURED AND TESTED FOLLOWING RELEVANT CGMPs UNDER A CERTIFIED QUALITY MANAGEMENT SYSTEM. PRODUCTS ARE FOR RESEARCH USE ONLY AND NOT INTENDED FOR HUMAN OR ANIMAL DIAGNOSTIC OR THERAPEUTIC USES UNLESS OTHERWISE STATED. FOR ADDITIONAL QUALITY INFORMATION, REFER TO [WWW.STEMCELL.COM/COMPLIANCE](http://WWW.STEMCELL.COM/COMPLIANCE).

Copyright © 2020 by STEMCELL Technologies Inc. All rights reserved including graphics and images. STEMCELL Technologies & Design, STEMCELL Shield Design, Scientists Helping Scientists, and StemSpan are trademarks of STEMCELL Technologies Inc. All other trademarks are the property of their respective holders. While STEMCELL has made all reasonable efforts to ensure that the information provided by STEMCELL and its suppliers is correct, it makes no warranties or representations as to the accuracy or completeness of such information.