

# Human Peripheral Blood Leukopak, Frozen

Primary human cells, frozen

<b>Catalog #200-0470</b>	Tenth Size
<b>Catalog #200-0132</b>	Quarter Size
<b>Catalog #200-0131</b>	Half Size
<b>Catalog #200-0130</b>	Full Size

## Product Description

Leukapheresis is performed on normal donors using Institutional Review Board (IRB)- or Research Ethics Committee (REC)-approved consent forms and protocols. Approximately two to three blood volumes are processed using the Spectra Optia<sup>®</sup> Apheresis System to produce a full-sized Leukopak. The collected product is then cryopreserved in a controlled-rate freezer.

<b>Donor Status:</b>	Normal
<b>Characterization Criteria:</b>	Age, Cell Count, Donor Viral Testing, Ethnicity, Height, Other Information, Sex, Smoking Status, Viability, Weight
<b>Product Format:</b>	Product is drawn into a sample collection bag containing anticoagulant and is frozen in CryoStor <sup>®</sup> CS10.
<b>Anticoagulant:</b>	Acid-citrate-dextrose solution A (ACDA)

For donor details, refer to the lot-specific Certificate of Analysis.

## Stability and Storage

Product is shipped on dry ice with liquid nitrogen shipping available upon request. Frozen leukopak products should be stored at -135°C or colder upon receipt. Thawed samples must be used immediately. As these are primary cells, they have a finite lifespan in culture.

## Precautions

**Donor Screening:** Donors are screened for HIV-1, HIV-2, hepatitis B, and hepatitis C.

Cryopreserved products are shipped with negative test results from donor screening that is performed within 90 days of collection.

Donors have been tested and found to be negative for HIV-1 and 2, hepatitis B, and hepatitis C prior to donation. As testing cannot completely guarantee that the donor was virus-free, **THIS PRODUCT SHOULD BE TREATED AS POTENTIALLY INFECTIOUS** and only used following appropriate handling precautions such as those described in biological safety level 2.

Storage of frozen cell products in the vapor phase of a liquid nitrogen storage tank is recommended. Storage in the liquid phase can result in cross-contamination if the seal is compromised or potential explosion when thawed if the liquid nitrogen manages to penetrate the product. Use of a face shield is required as a safety precaution when transferring cells from one container to another. When handling this product, do not use sharps such as needles and syringes.

STEMCELL cannot guarantee the specifications post thaw, the biological function or any other properties associated with performance of cells in a researcher's individual assay or culture systems. STEMCELL assures the cells will meet the specifications only when assessed immediately after thawing (before washing) by our test methods.

**FOR IN VITRO RESEARCH USE ONLY. NOT APPROVED FOR DIAGNOSTIC, THERAPEUTIC, OR CLINICAL APPLICATIONS. NOT APPROVED FOR HUMAN OR VETERINARY USE IN VIVO.**

## Directions for Use

**IMPORTANT:** The viability and number of cells in this product were assessed prior to cryopreservation, and some cell loss is expected upon thaw. To confirm the number of cells provided, a viable cell count must be done immediately after thawing. Work quickly once the Leukopak has been thawed to ensure high viability and recovery. Use sterile technique when processing thawed Leukopak.

**WARNING:** Handle frozen Leukopaks with care and inspect thoroughly upon receipt. Ensure the cryobag, tubing, and ports are intact and properly sealed prior to thawing. Freezing causes the plastic to become rigid, and while rare, can lead to breakage if improperly handled. Use extra care and avoid impacts when handling, transporting, and thawing the frozen product.

The following instructions are for thawing Leukopaks. Frozen leukopaks may come with external packaging; remove all external packaging prior to thawing for optimal performance. Instead of using a water bath (steps 3 - 4), frozen Leukopaks (Catalog #200-0130, #200-0131, #200-0132) can be thawed using ThawSTAR® CB Automated Thawing System (Catalog #100-1151).

**NOTE:** The ThawSTAR® CB Automated Thawing System is compatible with full, half, or quarter size frozen Leukopaks. It is not recommended for the tenth size frozen Leukopaks (Catalog #200-0470), and the water bath is preferred. For complete instructions, refer to the Product Information Sheet (Document #10000018736), available at [www.stemcell.com](http://www.stemcell.com), or contact us to request a copy.

1. Warm sufficient volume of HBSS Modified (Without Ca++ and Mg++) + 10% fetal bovine serum (FBS) in a 37°C water bath. Wipe with 70% ethanol and transfer to the biosafety cabinet.
2. Add DNase I to prepare thawing medium (HBSS Modified [Without Ca++ and Mg++] + 10% FBS + 0.1 mg/mL DNase I [final concentration]).  
**NOTE:** Do not add DNase I if the cells will be used for DNA or RNA extraction.  
**NOTE:** Thawing medium must be prepared fresh before each use.
3. Remove frozen Leukopak from liquid nitrogen storage and immediately place in a 37°C water bath. Submerge the Leukopak, but keep the ports elevated above the waterline, and do not agitate the bag while it is thawing.
4. When the Leukopak is mostly thawed (with a small amount of ice remaining), remove from the water bath and wipe the outside of the bag with 70% ethanol.
5. In a biosafety cabinet, slowly transfer the cell suspension into a sterile bottle that can contain 5 - 10 times the volume of the Leukopak, using the port on the bottom of the cryobag.
6. While gently swirling the cell suspension, add an equal volume of thawing medium dropwise to the bottle.
7. Optional: To improve recovery, rinse bag by adding one-half volume of thawing medium (relative to original Leukopak volume) to the Leukopak bag. Mix thoroughly and transfer to the cell suspension in the sterile bottle.
8. Slowly add an additional two volumes (relative to original Leukopak volume) of thawing medium to the cell suspension.
9. Thoroughly mix the sample (do not vortex).
10. Measure the volume precisely. Remove a small aliquot of cell suspension for counting and viability assessment. See Notes and Tips section

for performing cell counts with a hemocytometer.

11. Transfer the cell suspension to 50 mL conical tubes and centrifuge at 300 x g for 10 minutes at room temperature (15 - 25°C).
12. Carefully remove the supernatant (from step 11) with a pipette, leaving a small amount of medium to ensure the cell pellet is not disturbed. Resuspend the cell pellet by gently flicking the tube.
13. Optional: If cells are starting to clump, add 0.1 mg DNase I per mL of cell suspension and incubate for 15 minutes at room temperature.
14. Gently add at least 15 - 20 mL of thawing medium to each tube.
15. Centrifuge the cell suspension at 300 x g for 10 minutes at room temperature.
16. Carefully remove the supernatant with a pipette, leaving a small amount of medium to ensure the cell pellet is not disturbed. Resuspend the cell pellet by gently flicking the tube and adding desired medium.  
NOTE: Cell losses are expected during the wash steps.
17. Cells are now ready for use in downstream applications.

NOTE: Platelet and red blood cell (RBC) levels will vary between samples, and further processing may be required prior to use in downstream applications (e.g. slow spins to remove platelets, RBC lysis, or density gradient separation). See Notes and Tips section for an optional RBC Lysis and Platelet Removal protocol (after diluting and centrifuging the sample). For downstream EasySep™ isolations, it is recommended to perform an additional 1 - 2 washing steps before use..

## Notes and Tips

For a protocol on performing total nucleated cell counts using a hemocytometer, refer to <https://www.stemcell.com/technical-resources/how-to-count-cells-with-a-hemocytometer.html>

For a protocol on RBC Lysis and Platelet Removal, refer to "Processing a Leukopak for Downstream Cell Isolation" (Part II: Prepare Leukopak Contents, Option 2), available at <https://www.stemcell.com/leukopak-processing-protocol.html>.

## Accessory Products

- DNase I Solution (1 mg/mL; Catalog #07900)
- Falcon® Conical Tubes, 50 mL (Catalog #38010)
- Falcon® Serological Pipettes, 2 mL (Catalog #38002)
- Hausser Scientific™ Bright-Line Hemocytometer (Catalog #100-1181)
- HBSS, Modified (Without Ca++ and Mg++; Catalog #37250)
- ThawSTAR® CFT2 Automated Thawing System (Catalog #100-0650)
- Trypan Blue (Catalog #07050)

PRODUCTS ARE FOR RESEARCH USE ONLY AND NOT INTENDED FOR HUMAN OR ANIMAL DIAGNOSTIC OR THERAPEUTIC USES UNLESS OTHERWISE STATED.

Copyright © 2024 by STEMCELL Technologies Inc. All rights reserved including graphics and images. STEMCELL Technologies & Design, STEMCELL Shield Design, Scientists Helping Scientists, and EasySep are trademarks of STEMCELL Technologies Canada Inc. CryoStor and ThawStar are registered trademarks of BioLife Solutions. Falcon is a registered trademark of Corning Incorporated. All other trademarks are the property of their respective holders. While STEMCELL has made all reasonable efforts to ensure that the information provided by STEMCELL and its suppliers is correct, it makes no warranties or representations as to the accuracy or completeness of such information.