

Primary Cells

**Human Cancer - Whole Peripheral
Blood Collected with Sodium
Heparin, Fresh**



Scientists Helping Scientists™ | WWW.STEMCELL.COM

TOLL FREE PHONE 1 800 667 0322 • PHONE +1 604 877 0713

INFO@STEMCELL.COM • TECHSUPPORT@STEMCELL.COM

FOR GLOBAL CONTACT DETAILS VISIT OUR WEBSITE

Product Description

Whole peripheral blood (PB) is obtained from donors diagnosed with cancer using Institutional Review Board (IRB)-approved consent forms and protocols.

Donor Status: Cancer

Characterization Criteria: Primary cancer type, stage, donor ID, cell processing date, age, sex, ethnicity, weight, height, smoking status, blood type, anticoagulant, current treatments, other information

Format: Product is drawn directly into a Vacutainer® tube containing anticoagulant.

Anticoagulant: Sodium heparin (Na-Hep)

For donor details, refer to the lot-specific Certificate of Analysis.

Ordering Information

CANCER TYPE	CATALOG #
Bladder Cancer	200-0260
Breast Cancer	200-0261
Cervical Cancer	200-0262
Colorectal Cancer	200-0263
Endometrial Cancer	200-0264
Esophageal Cancer	200-0265
Gastric Cancer	200-0266
Head and Neck Cancer	200-0267
Kidney Cancer	200-0268
Liver Cancer	200-0269
Lung Cancer	200-0270
Melanoma	200-0271
Ovarian Cancer	200-0272
Pancreatic Cancer	200-0273
Prostate Cancer	200-0274
Solid Tumor Cancer	200-0400

Stability and Storage

Product is shipped at ambient temperature and should be used or processed immediately upon receipt.

Precautions

Donor Screening: Donors are screened for HIV-1, HIV-2, hepatitis B, and hepatitis C.

If the donor has been screened within 90 days of donation the product will be shipped with negative test results from donor screening.

If the donor has not been screened within 90 days of collection, a test sample will be taken at the time of donation and the product will be shipped before the screening results are available. In the unlikely event that a test result is positive, the customer will be contacted as soon as possible (usually within 24 - 72 hours from the time of shipment).

Donors have been tested and found to be negative for HIV-1, HIV-2, hepatitis B, and hepatitis C prior to donation. As testing cannot completely guarantee that the donor was virus-free, THIS PRODUCT SHOULD BE TREATED AS POTENTIALLY INFECTIOUS and only used following appropriate handling precautions such as those described in biological safety level 2. When handling this product, do not use sharps such as needles and syringes.

STEMCELL cannot guarantee the biological function or any other properties associated with performance of cells in a researcher's individual assay or culture systems. STEMCELL assures the cells will meet the specifications only when assessed, before washing, by our test methods.

FOR IN VITRO RESEARCH USE ONLY. NOT APPROVED FOR DIAGNOSTIC, THERAPEUTIC, OR CLINICAL APPLICATIONS. NOT APPROVED FOR HUMAN OR VETERINARY USE IN VIVO.

Handling / Directions for Use

IMPORTANT: To determine the number of cells provided, a cell count must be done upon receipt and before any processing steps (e.g. washing). Cell loss is expected during wash steps and may be up to 30%. Use sterile technique when processing cells.

Remove a 20 µL aliquot of cells for counting. Appropriately dilute in Trypan Blue (to assess viability) or 3% Acetic Acid with Methylene Blue (to assess nucleated cells). Adjust the dilution if there are more than 100 cells per square of the hemocytometer. See Tips section for more details on performing cell counts with a hemocytometer.

Tips

For a protocol on performing total nucleated cell counts using a hemocytometer, refer to <https://www.stemcell.com/how-to-count-cells-with-a-hemocytometer>.

Accessory Products

PRODUCT NAME	CATALOG #
Trypan Blue	07050
3% Acetic Acid with Methylene Blue	07060
Ammonium Chloride Solution	07800
Lymphoprep™	07801
SepMate™-50 (IVD)	85450
SepMate™-50 (RUO)	86450

PRODUCTS ARE FOR RESEARCH USE ONLY AND NOT INTENDED FOR HUMAN OR ANIMAL DIAGNOSTIC OR THERAPEUTIC USES UNLESS OTHERWISE STATED.

Copyright © 2021 by STEMCELL Technologies Inc. All rights reserved including graphics and images. STEMCELL Technologies & Design, STEMCELL Shield Design, Scientists Helping Scientists, and SepMate are trademarks of STEMCELL Technologies Canada Inc. Vacutainer is a registered trademark of Becton, Dickinson and Company. Lymphoprep is a trademark of Alere Technologies. All other trademarks are the property of their respective holders. While STEMCELL has made all reasonable efforts to ensure that the information provided by STEMCELL and its suppliers is correct, it makes no warranties or representations as to the accuracy or completeness of such information.