

# Primary Cells

## Human Cancer - Peripheral Blood Leukopak, Fresh



Scientists Helping Scientists™ | [WWW.STEMCELL.COM](http://WWW.STEMCELL.COM)

TOLL FREE PHONE 1 800 667 0322 • PHONE +1 604 877 0713

[INFO@STEMCELL.COM](mailto:INFO@STEMCELL.COM) • [TECHSUPPORT@STEMCELL.COM](mailto:TECHSUPPORT@STEMCELL.COM)

FOR GLOBAL CONTACT DETAILS VISIT OUR WEBSITE

## Product Description

Cancer-sourced Leukopaks are an enriched leukapheresis product consisting of a large number of single-donor mononuclear cells (MNCs) collected from peripheral blood of a donor who has been diagnosed with cancer. Leukapheresis is performed on donors diagnosed with cancer using Institutional Review Board (IRB)-approved consent forms and protocols.

**Donor Status:** Cancer

**Characterization Criteria:** Primary cancer type, stage, donor ID, cell processing date, age, sex, ethnicity, weight, height, smoking status, blood type, anticoagulant, current treatments, other information

**Format:** Product is drawn directly into a sample collection bag containing anticoagulant and transferred into a bottle, bag, or conical tube.

**Anticoagulant:** Acid-citrate-dextrose solution A (ACDA)

For donor details, refer to the lot-specific Certificate of Analysis.

## Ordering Information

CANCER TYPE	CATALOG #
Bladder Cancer	200-0290
Breast Cancer	200-0291
Cervical Cancer	200-0292
Colorectal Cancer	200-0293
Endometrial Cancer	200-0294
Esophageal Cancer	200-0295
Gastric Cancer	200-0296
Head and Neck Cancer	200-0297
Kidney Cancer	200-0298
Liver Cancer	200-0299
Lung Cancer	200-0300
Melanoma	200-0301
Ovarian Cancer	200-0302
Pancreatic Cancer	200-0303
Prostate Cancer	200-0304
Solid Tumor Cancer	200-0402

## Stability and Storage

Product is shipped at 2 - 8°C and should be used or processed immediately upon receipt.

## Precautions

Donor Screening: Donors are screened for HIV-1, HIV-2, hepatitis B, and hepatitis C.

If the donor has been screened within 90 days of donation the product will be shipped with negative test results from donor screening.

If the donor has not been screened within 90 days of collection, a test sample will be taken at the time of donation and the product will be shipped before the screening results are available. In the unlikely event that a test result is positive, the customer will be contacted as soon as possible (usually within 24 - 72 hours from the time of shipment).

Donors have been tested and found to be negative for HIV-1, HIV-2, hepatitis B, and hepatitis C prior to donation. As testing cannot completely guarantee that the donor was virus-free, THIS PRODUCT SHOULD BE TREATED AS POTENTIALLY INFECTIOUS and only used following appropriate handling precautions such as those described in biological safety level 2. When handling this product, do not use sharps such as needles and syringes.

STEMCELL cannot guarantee the biological function or any other properties associated with performance of cells in a researcher's individual assay or culture systems. STEMCELL assures the cells will meet the specifications only when assessed, before washing, by our test methods.

FOR IN VITRO RESEARCH USE ONLY. NOT APPROVED FOR DIAGNOSTIC, THERAPEUTIC, OR CLINICAL APPLICATIONS. NOT APPROVED FOR HUMAN OR VETERINARY USE IN VIVO.

## Directions for Use

**IMPORTANT:** To determine the number of cells provided, a cell count must be done upon receipt and before any processing steps (e.g. washing). Cell loss is expected during wash steps and may be up to 30%. Use sterile technique when processing cells.

Remove a 20 µL aliquot of cells for counting. Appropriately dilute in Trypan Blue (to assess viability) or 3% Acetic Acid with Methylene Blue (to assess nucleated cells). For most Leukopak samples, a dilution of 1 in 100 is sufficient. Adjust the dilution if there are more than 100 cells per square of the hemocytometer. See Tips section for more details on performing cell counts with a hemocytometer.

NOTE: SepMate™ tubes are not intended for use with leukapheresis samples.

## Tips

For a protocol on performing total nucleated cell counts using a hemocytometer, refer to <https://www.stemcell.com/how-to-count-cells-with-a-hemocytometer>.

## Accessory Products

PRODUCT NAME	CATALOG #
Trypan Blue	07050
3% Acetic Acid with Methylene Blue	07060
Ammonium Chloride Solution	07800
Lymphoprep™	07801

PRODUCTS ARE FOR RESEARCH USE ONLY AND NOT INTENDED FOR HUMAN OR ANIMAL DIAGNOSTIC OR THERAPEUTIC USES UNLESS OTHERWISE STATED.

Copyright © 2022 by STEMCELL Technologies Inc. All rights reserved including graphics and images. STEMCELL Technologies & Design, STEMCELL Shield Design, Scientists Helping Scientists, and SepMate are trademarks of STEMCELL Technologies Canada Inc. Lymphoprep is a trademark of Alere Technologies. All other trademarks are the property of their respective holders. While STEMCELL has made all reasonable efforts to ensure that the information provided by STEMCELL and its suppliers is correct, it makes no warranties or representations as to the accuracy or completeness of such information.