

# TeSR™-AOF 3D

**Animal origin-free suspension culture medium for human ES and iPS cells**

Catalog #100-0720 1 Kit



Scientists Helping Scientists™ | [WWW.STEMCELL.COM](http://WWW.STEMCELL.COM)

TOLL FREE PHONE 1 800 667 0322 • PHONE +1 604 877 0713

[INFO@STEMCELL.COM](mailto:INFO@STEMCELL.COM) • [TECHSUPPORT@STEMCELL.COM](mailto:TECHSUPPORT@STEMCELL.COM)

FOR GLOBAL CONTACT DETAILS VISIT OUR WEBSITE

## Product Description

TeSR™-AOF 3D is an animal origin-free culture medium developed for expansion and scale-up of undifferentiated human embryonic stem (ES) cells and human induced pluripotent stem (iPS) cells as aggregates in 3D suspension culture. Based on TeSR™-AOF, TeSR™-AOF 3D is a robust medium with an optimized fed-batch protocol for scaling up human ES and iPS cell production in suspension culture. The novel fed-batch workflow saves time and media; daily feeds replenish nutrients and eliminate the need for medium exchanges on non-passaging days.

The TeSR™-AOF 3D system is compatible with cells previously cultured in mTeSR™1 (Catalog #85850), mTeSR™ Plus (Catalog #100-0276), TeSR™-AOF (Catalog #100-0401), or TeSR™-E8™ (Catalog #05990) on Corning® Matrigel® hESC-Qualified Matrix (Corning Catalog #354277), Vitronectin XF™ (Catalog #07180), or CellAdhere™ Laminin-521 (Catalog #77003). TeSR™-AOF 3D is compatible with a range of suspension culture vessels.

No materials of animal or human origin are used in the manufacture of this medium or its components, to at least the secondary level of manufacturing. TeSR™-AOF 3D is animal component-free to the secondary level. Each lot of TeSR™-AOF 3D Seed or Feed Supplement that is used to prepare complete TeSR™-AOF medium is performance-tested in a culture assay using human pluripotent stem cells.

## Product Information

The following components are sold as part of the TeSR™-AOF 3D kit (Catalog #100-0720) and are not available for individual sale.

COMPONENT NAME	COMPONENT #	SIZE	STORAGE	SHELF LIFE
TeSR™-AOF 3D Seed Basal Medium	100-0721	480 mL	Store at 2 - 8°C.	Stable for 12 months from date of manufacture (MFG) on label.
TeSR™-AOF 3D Seed Supplement	100-0722	20 mL	Store at -20°C.	Stable for 12 months from date of manufacture (MFG) on label.
TeSR™-AOF 3D Feed Supplement	100-0723	20 mL	Store at -20°C.	Stable for 12 months from date of manufacture (MFG) on label.

## Preparation of Media

### A. TeSR™-AOF 3D Seed Medium

Use sterile technique to prepare TeSR™-AOF 3D Seed Medium (Basal Medium + Seed Supplement + Y-27632 [Catalog #72302]).

The following example is for preparing 500 mL of medium. If preparing other volumes, adjust accordingly.

1. Thaw TeSR™-AOF 3D Seed Supplement at room temperature (15 - 25°C) or overnight at 2 - 8°C. Do not thaw in a 37°C water bath. Mix thoroughly.

NOTE: Precipitate may be visible; warm to room temperature until dissolved. Performance will not be affected.

NOTE: If not used immediately, store TeSR™-AOF 3D Seed Supplement in one of the following containers:

- TeSR™-AOF 3D Seed Supplement bottle
- 50 mL polypropylene tubes (e.g. Catalog #38010)
- 2 mL or 0.5 mL polypropylene micro tubes (e.g. Sarstedt Catalog #72.694.006 or 72.730.005)

Do not use other storage containers. Store TeSR™-AOF 3D Seed Supplement at -20°C for up to 3 months. Do not exceed the shelf life of the supplement. After thawing the aliquots, use immediately. Do not re-freeze.

2. Add 20 mL of TeSR™-AOF 3D Seed Supplement to 480 mL of TeSR™-AOF 3D Seed Basal Medium. Mix thoroughly. If desired, the medium can be filtered using a 0.2 µm low-protein binding filter.

NOTE: If not used immediately, store TeSR™-AOF 3D Seed Medium (without Y-27632) in one of the following containers:

- TeSR™-AOF 3D Seed Basal Medium bottle
- 50 mL polypropylene tubes (e.g. Catalog #38010)

Do not use other storage containers. Store TeSR™-AOF 3D Seed Medium (without Y-27632) at 2 - 8°C for up to 2 weeks. Alternatively, aliquot and store at -20°C for up to 6 months. Do not exceed the shelf life of the individual components. After thawing the aliquots, proceed to step 3 and use immediately, or store at 2 - 8°C for up to 2 weeks. Do not re-freeze.

3. Immediately before use, add 10 µM Y-27632. Mix thoroughly.

## B. TeSR™-AOF 3D Feed Supplement

Thaw TeSR™-AOF 3D Feed Supplement at room temperature (15 - 25°C) or overnight at 2 - 8°C. Do not thaw in a 37°C water bath. Mix thoroughly.

NOTE: Precipitate may be visible, and/or supplement may appear cloudy. Warm the supplement to room temperature until all precipitate has dissolved and the solution no longer appears cloudy. Performance will not be affected.

NOTE: If not used immediately, store TeSR™-AOF 3D Feed Supplement in one of the following containers:

- TeSR™-AOF 3D Feed Supplement bottle
- 50 mL polypropylene tubes (e.g. Catalog #38010)
- 2 mL or 0.5 mL polypropylene micro tubes (e.g. Sarstedt Catalog #72.694.006 or 72.730.005)

Do not use other storage containers. Store TeSR™-AOF 3D Feed Supplement at 2 - 8°C for up to 2 weeks. Alternatively, aliquot and store at -20°C for up to 3 months. Do not exceed the shelf life of the supplement. After thawing the aliquots, use immediately or store at 2 - 8°C for up to 2 weeks. Do not re-freeze.

## Directions for Use

For complete instructions on how to culture human ES and iPS cells in TeSR™-AOF 3D, refer to the Technical Manual: Expansion of Human Pluripotent Stem Cells as Aggregates in Suspension Culture Using TeSR™-AOF 3D (Document #10000010775), available at [www.stemcell.com](http://www.stemcell.com), or contact us to request a copy.



TeSR™-AOF 3D is manufactured and sold under global exclusive license from Accellta for culture medium for hPSCs in suspension under feeder-free, non-adherent conditions.



This product was developed under license to intellectual property owned by WiCell™ Research Institute.

This product is sold for research use only (whether the buyer is an academic or for-profit entity) under a non-transferable, limited-use license. Purchase of this product does not include the right to sell, use or otherwise transfer this product for commercial purposes (i.e., any activity undertaken for consideration, such as use of this product for manufacturing, or resale of this product or any materials made using this product, or use of this product or any materials made using this product to provide services) or clinical use (i.e., administration of this product or any material using this product to humans) or the right to implant any material made using this product into an animal by, or in collaboration with, a for-profit entity, for purposes other than basic pre-clinical research applications (including without limitation teratoma assays) to validate the function of the cells. Purchasers who do not agree to the terms and conditions set forth above should return the product in acceptable conditions to the seller for a refund.

PRODUCTS ARE FOR RESEARCH USE ONLY AND NOT INTENDED FOR HUMAN OR ANIMAL DIAGNOSTIC OR THERAPEUTIC USES UNLESS OTHERWISE STATED. FOR ADDITIONAL INFORMATION ON QUALITY AT STEMCELL, REFER TO [WWW.STEMCELL.COM/COMPLIANCE](http://WWW.STEMCELL.COM/COMPLIANCE).

Copyright © 2022 by STEMCELL Technologies Inc. All rights reserved including graphics and images. STEMCELL Technologies & Design, STEMCELL Shield Design, CellAdhere, and Scientists Helping Scientists are trademarks of STEMCELL Technologies Canada Inc. Corning and Matrigel are registered trademarks of Corning Incorporated. mTeSR, TeSR, and E8 are trademarks of WARF. Vitronectin XF is developed and manufactured by Nucleus Biologics, and Vitronectin XF is a trademark of Nucleus Biologics. All other trademarks are the property of their respective holders. While STEMCELL has made all reasonable efforts to ensure that the information provided by STEMCELL and its suppliers is correct, it makes no warranties or representations as to the accuracy or completeness of such information.