

# Dyes and Stains

## Phalloidin iFluor™ Conjugates

Fluorescent dye conjugates for detection of F-actin in cell culture experiments



Scientists Helping Scientists™ | [WWW.STEMCELL.COM](http://WWW.STEMCELL.COM)

TOLL FREE PHONE 1 800 667 0322 • PHONE +1 604 877 0713

[INFO@STEMCELL.COM](mailto:INFO@STEMCELL.COM) • [TECHSUPPORT@STEMCELL.COM](mailto:TECHSUPPORT@STEMCELL.COM)

FOR GLOBAL CONTACT DETAILS VISIT OUR WEBSITE

|                   |           |
|-------------------|-----------|
| Catalog #100-0852 | 300 Tests |
| Catalog #100-0853 | 300 Tests |
| Catalog #100-0854 | 300 Tests |

## Product Description

Phalloidin is a heptapeptide derived from the poisonous mushroom *Amanita phalloides*. Phalloidin binds to F-actin filaments in cells with high specificity, preventing filament depolymerization. Fluorescently labeled phalloidin can be used in labeling, identification, and quantitation of F-actin in formaldehyde-fixed and permeabilized cells. Fluorescently labeled phalloidin has lower non-specific binding compared to actin antibodies across most species of plants and animals, which provides low background and high contrast during imaging.

## Product Information

| PRODUCT NAME            | CATALOG # | SIZE      | STORAGE   | SHELF LIFE                                   |
|-------------------------|-----------|-----------|---|--|
| Phalloidin, iFluor™ 488 | 100-0852  | 300 Tests | Store at -20°C. Protect product from prolonged exposure to light. | Stable until expiry date (EXP) on box label. |
| Phalloidin, iFluor™ 555 | 100-0853  | 300 Tests | Store at -20°C. Protect product from prolonged exposure to light. | Stable until expiry date (EXP) on box label. |
| Phalloidin, iFluor™ 647 | 100-0854  | 300 Tests | Store at -20°C. Protect product from prolonged exposure to light. | Stable until expiry date (EXP) on box label. |

## Physical & Spectral Properties

| PRODUCT NAME             | FORMULATION                                 | MOLECULAR WEIGHT (g/mol) | EXCITATION WAVELENGTH (nm) | EMISSION WAVELENGTH (nm) | EXTINCTION COEFFICIENT (cm <sup>-1</sup> M <sup>-1</sup> )** | QUANTUM YIELD** | CORRECTION FACTOR (280 nm) |
|--------------------------|---|--------------------------|----------------------------|--------------------------|--|-----------------|----------------------------|
| Phalloidin, iFluor™ 488* | 1000X solution in dimethyl sulfoxide (DMSO) | 1399.49                  | 491                        | 516                      | 75,000   | 0.9             | 0.11                       |
| Phalloidin, iFluor™ 555* |   | 1468.53                  | 557                        | 570                      | 100,000  | 0.64            | 0.14                       |
| Phalloidin, iFluor™ 647  | Powder                                      | 1632.87                  | 656                        | 670                      | 250,000  | 0.25            | 0.03                       |

\*Please refer to the Safety Data Sheet (SDS) for hazard information. This product contains components dissolved in DMSO. DMSO is a strong solvent and skin penetrant, and can transport many substances through the skin. DMSO can also penetrate some protective glove materials including latex and silicone. Extra caution should be utilized when handling this product.

\*\*Measured with aqueous buffer (pH 7.2).

## Directions for Use

Please read the entire protocol before proceeding. The following protocol is for staining cells in a black-wall/clear-bottom 96-well plate. If using other cultureware, adjust volumes accordingly. Staining may vary with cell type; the optimal concentration of working solution should be determined.

### Preparation of Phalloidin iFluor™ Working Solution

1. Warm the Phalloidin iFluor™ vial to room temperature (15 - 25°C) and centrifuge briefly before opening.
2. Prepare a Phalloidin iFluor™ working solution as follows:
  - iFluor™ 488 and 555 conjugates: Add 1 µL of solution from the vial to 1 mL of phosphate-buffered saline (PBS) with 1% bovine serum albumin (BSA). Mix thoroughly. Use the working solution immediately; do not store.

NOTE: If not used immediately, aliquot and store remaining solution from the vial at -20°C. Seal tubes tightly and protect from light. After thawing aliquots, use immediately; do not re-freeze.

- iFluor™ 647 conjugate:

a. To prepare a Phalloidin, iFluor™ 647 stock solution, add 30 µL of DMSO to the vial. Mix thoroughly.

NOTE: If not used immediately, aliquot and store remaining solution from the vial at -20°C. Seal tubes tightly and protect from light. After thawing aliquots, use immediately; do not re-freeze.

b. To prepare a Phalloidin, iFluor™ 647 working solution, add 1 µL of stock solution from the vial to 1 mL of PBS with 1% BSA. Mix thoroughly. Use the working solution immediately; do not store.

### Staining and Imaging Cells

1. Stain cells as follows:

a. Plate cells in a black-wall/clear-bottom 96-well plate and treat cells as desired.

b. Add 4% methanol-free formaldehyde in PBS to the cells, then incubate at room temperature for 10 - 30 minutes.

NOTE: Avoid using a methanol-containing fixative, since methanol can disrupt actin during the fixation process.

c. Remove fixative and rinse the fixed cells 2 - 3X with PBS.

d. If desired, add 0.1% Triton™ X-100 in PBS to fixed cells and incubate at room temperature for 3 - 5 minutes to increase permeability. Rinse cells 2 - 3X with PBS.

e. Add 100 µL/well of Phalloidin iFluor™ working solution to the fixed cells. Incubate at room temperature for 20 - 90 minutes; protect from light.

NOTE: The optimal incubation time should be determined for different cell types.

f. Rinse cells gently 2 - 3X with PBS to remove excess phalloidin conjugate.

2. Observe stained cells using a fluorescence microscope with the appropriate filter set, as follows:

- iFluor™ 488 conjugate: FITC

- iFluor™ 555 conjugate: Cy3/TRITC

- iFluor™ 647 conjugate: Cy5

PRODUCTS ARE FOR RESEARCH USE ONLY AND NOT INTENDED FOR HUMAN OR ANIMAL DIAGNOSTIC OR THERAPEUTIC USES UNLESS OTHERWISE STATED.

Copyright © 2021 by STEMCELL Technologies Inc. All rights reserved including graphics and images. STEMCELL Technologies & Design, STEMCELL Shield Design, and Scientists Helping Scientists are trademarks of STEMCELL Technologies Canada Inc. iFluor is a trademark of AAT Bioquest, Inc. Triton is a trademark of Dow Chemical Co. All other trademarks are the property of their respective holders. While STEMCELL has made all reasonable efforts to ensure that the information provided by STEMCELL and its suppliers is correct, it makes no warranties or representations as to the accuracy or completeness of such information.