

StemSpan™-AOF

cGMP, animal origin-free medium for culture and expansion of human hematopoietic cells

Catalog #100-0513

1 L Bag



Scientists Helping Scientists™ | WWW.STEMCELL.COM

TOLL FREE PHONE 1 800 667 0322 • PHONE +1 604 877 0713

INFO@STEMCELL.COM • TECHSUPPORT@STEMCELL.COM

FOR GLOBAL CONTACT DETAILS VISIT OUR WEBSITE

Product Description

StemSpan™-AOF is an animal origin-free (AOF) medium that has been developed for the in vitro culture and expansion of human hematopoietic cells, when appropriate growth factors and supplements are added. This allows users the flexibility to prepare medium that meets their requirements. StemSpan™-AOF contains only recombinant proteins and synthetic components, and does not contain serum or other human- or animal-derived components.

Using appropriate cytokines, StemSpan™-AOF may be used to expand CD34+ cells isolated from human cord blood, mobilized peripheral blood, or bone marrow samples, or to expand and differentiate lineage-committed progenitor cells to generate populations of myeloid or megakaryocyte progenitor cells.

StemSpan™-AOF is manufactured and tested following relevant cGMPs under a certified quality management system. For additional quality information, refer to www.stemcell.com/compliance.

Properties

Storage: Store at 2 - 8°C. Do not freeze. Do not stack. Store flat with tubing facing up. Outermost layers of packaging and bubble wrap may be removed if desired. Do not remove clear, vacuum-sealed polybag and inner protective bubble wrap until use. Innermost packaging is designed for optimal product storage, and to protect the flexible tubing and connectors.

Shelf Life: Stable until expiry date (EXP) on label.

Contains: This product contains only recombinant human proteins and synthetic components.

Directions for Use

PREPARATION OF COMPLETE MEDIUM

Instructions below serve as a general recommendation for the use of StemSpan™-AOF in a flexible 1 L bag. Protocol may vary depending on application and needs. Use sterile technique to prepare complete StemSpan™-AOF.

1. Remove the bag of StemSpan™-AOF from the vacuum-sealed polybag by grasping the flexible bag of medium by the handle at top of bag opposite tubing. Slide the tubing out of the protective bubble wrap, taking care to avoid leakage or damage.
2. Add desired cytokines, growth factors and, optionally, other supplements to StemSpan™-AOF using a method of choice, as follows:
 - a. Syringe transfer via needleless injection site
 - b. Appropriate connection to multipoint-connector (MPC)
 - c. Tube weld to allow aseptic connection

Mix thoroughly.

NOTE: Added supplements and cells in sterile cell culture medium should not exceed ~10% of total volume.

3. Connect the bag to the closed system or remove medium, as desired, by method of choice described in step 2.

For a complete list of expansion supplements for use with StemSpan™-AOF, see Notes and Tips.

ASSESSMENT OF HEMATOPOIETIC CELLS

Assessment of CD34+ cells before and after culture may be performed by flow cytometry using the following fluorochrome-conjugated antibody clones:

- Anti-Human CD34 Antibody, Clone 581 (Catalog #60013) or Clone 563 (Catalog #60119) or Clone 8G12 (Catalog #60121)
- Anti-Human CD45 Antibody, Clone HI30 (Catalog #60018) or Clone 2D1 (Catalog #60123)
- Anti-Human CD38 Antibody, Clone AT-1 (Catalog #60131) or Clone HIT2 (Catalog #60014)
- Anti-Human CD90 Antibody, Clone 5E10 (Catalog #60045)

Diagram




FIGURE 1. Overview of Connections and Tubing for StemSpan™-AOF (1 L bag; Catalog #100-0513)

A schematic showing the design of the StemSpan™-AOF, 1 L bag with the following lines:

- (1) Sealed line (1/4" ID x 3/8" OD) used for manufacturing
- (2) 18" of weldable/sealable tubing (1/4" ID x 3/8" OD) with MPC and ratchet pinch clamp on line
- (3) Consists of 2 segments: 3" of tubing (1/4" ID x 3/8" OD), and 18" of weldable/sealable tubing (1/8" ID x 1/4" OD) with swabbable valve/needleless injection site (luer) and ratchet pinch clamp on line.

OD = outer diameter, ID = inner diameter

Notes and Tips

Selection of an optimal growth factor combination is dependent upon the source and type of cells and the experimental objectives of the researcher. StemSpan™ expansion supplements*, described below, are recommended for use with StemSpan™-AOF. StemSpan™ expansion supplements are manufactured and tested under a quality management system certified to ISO 9001. Refer to the Product Information Sheet (PIS) for each supplement for recommended cell expansion protocols with StemSpan™ media.

- StemSpan™ CD34+ Expansion Supplement (10X) (Catalog #02691)
 - Culture and expansion of large numbers of human CD34+ progenitor cells
 - Contains: rh SCF, rh TPO, rh IL-3, rh IL-6, rh Flt3 ligand, other additives
- StemSpan™ Megakaryocyte Expansion Supplement (100X) (Catalog #02696)
 - Expansion and lineage-specific differentiation of human CD34+ cells into megakaryocyte progenitor cells
 - Contains: rh SCF, rh TPO, rh IL-6, rh IL-9
- StemSpan™ Myeloid Expansion Supplement (100X) (Catalog #02693)
 - Expansion and lineage-specific differentiation of human CD34+ cells into granulocytes
 - Contains: rh SCF, rh TPO, rh G-CSF, rh GM-CSF
- StemSpan™ Myeloid Expansion Supplement II (100X) (Catalog #02694)
 - Expansion and lineage-specific differentiation of human CD34+ cells into monocytes
 - Contains: rh Flt3 ligand, rh SCF, rh TPO, rh M-CSF, rh GM-CSF, other additives

SCF = stem cell factor; EPO = erythropoietin; TPO = thrombopoietin; rh = recombinant human; IL = interleukin; Flt = fms-like tyrosine kinase

* These supplements are serum-free and are not manufactured to relevant cGMPs. Their suitability depends on your experimental requirements.

Related Products

For related products, including specialized culture and storage media, supplements, antibodies, cytokines, and small molecules, visit www.stemcell.com/HSPCworkflow, or contact us at techsupport@stemcell.com. For available fresh and cryopreserved peripheral blood, cord blood, and bone marrow products in your region, visit www.stemcell.com/primarycells.

THIS PRODUCT IS MANUFACTURED AND TESTED FOLLOWING RELEVANT cGMPs UNDER A CERTIFIED QUALITY MANAGEMENT SYSTEM. PRODUCT IS FOR FURTHER MANUFACTURING OR RESEARCH USE. NOT INTENDED FOR HUMAN OR ANIMAL DIAGNOSTIC OR THERAPEUTIC USES UNLESS OTHERWISE STATED. FOR ADDITIONAL QUALITY INFORMATION, REFER TO WWW.STEMCELL.COM/COMPLIANCE.

Copyright © 2023 by STEMCELL Technologies Inc. All rights reserved including graphics and images. STEMCELL Technologies & Design, STEMCELL Shield Design, Scientists Helping Scientists, and StemSpan are trademarks of STEMCELL Technologies Inc. All other trademarks are the property of their respective holders. While STEMCELL has made all reasonable efforts to ensure that the information provided by STEMCELL and its suppliers is correct, it makes no warranties or representations as to the accuracy or completeness of such information.