



EasyEights™ EasySep™ Magnet

Magnet for column-free immunomagnetic separation

Simultaneously processes up to 8 samples per side, total of 16 samples

Catalog # 18103



Scientists Helping Scientists™ | WWW.STEMCELL.COM

TOLL FREE PHONE 1 800 667 0322 • PHONE +1 604 877 0713

INFO@STEMCELL.COM • TECHSUPPORT@STEMCELL.COM

FOR GLOBAL CONTACT DETAILS VISIT OUR WEBSITE

Product Description

The EasyEights™ EasySep™ Magnet is an anodized aluminum body housing rare earth magnetic elements.

The EasyEights™ EasySep™ Magnet is designed for cell separation procedures using selected EasySep™ reagents when processing up to 16 samples at a time. It is designed to hold up to:

- 8 standard 5 mL (12 x 75 mm) polystyrene round-bottom tubes (e.g. Catalog #38007), and
 - 8 standard 14 mL (17 x 95 mm) polystyrene round-bottom tubes (e.g. Catalog #38008)
- OR
- 8 standard 15 mL (17 x 118 mm) conical tubes (e.g. Catalog #38009)

Properties

Storage: Store at 15 - 25°C.

Shelf Life: Not applicable

STRONG MAGNETIC FIELD: Keep away from pacemakers, magnets, computer disks, watches, and other objects that respond to magnetic fields. Do not remove the screws from the bottom of the magnet.

Handling / Directions For Use

The EasyEights™ EasySep™ Magnet is only for use with certain EasySep™ kits. For more information, refer to the appropriate EasySep™ Product Information Sheet or contact us at techsupport@stemcell.com.

One side of the EasyEights™ EasySep™ Magnet (STEMCELL Technologies logo) is designed to hold up to 8 x 5 mL (12 x 75 mm) polystyrene round-bottom tubes (e.g. Catalog #38007). The other side (EasySep™ logo) is designed to hold up to 8 x 14 mL (17 x 95 mm) polystyrene round-bottom tubes (e.g. Catalog #38008) or 8 x 15 mL (17 x 118 mm) conical tubes (e.g. Catalog #38009).

After completing a cell separation, use a serological pipette to remove the supernatant; do not pour off. Use either a 2 mL (e.g. Catalog #38002) or a 10 mL (e.g. Catalog #38004) serological pipette, depending on the separation volume. Collect the entire supernatant, all at once, into a single pipette.

CARE OF MAGNET: After completing a cell separation, wipe the magnet clean with distilled water and then with alcohol to dry the magnet. To disinfect the magnet, wipe with 10% diluted household bleach and then wipe off with a water-dampened cloth to avoid corrosion of metal parts. Do not submerge the magnet in liquid.



STEMCELL TECHNOLOGIES INC.'S QUALITY MANAGEMENT SYSTEM IS CERTIFIED TO ISO 13485. PRODUCTS ARE FOR RESEARCH USE ONLY AND NOT INTENDED FOR HUMAN OR ANIMAL DIAGNOSTIC OR THERAPEUTIC USES UNLESS OTHERWISE STATED.

Copyright © 2017 by STEMCELL Technologies Inc. All rights reserved including graphics and images. STEMCELL Technologies & Design, STEMCELL Shield Design, Scientists Helping Scientists, and EasyEights, EasySep are trademarks of STEMCELL Technologies Canada Inc. All other trademarks are the property of their respective holders. While STEMCELL has made all reasonable efforts to ensure that the information provided by STEMCELL and its suppliers is correct, it makes no warranties or representations as to the accuracy or completeness of such information.