



**RosetteSep™ Human  
Cord Blood Progenitor  
Enrichment Kit**

Catalog #15276

For processing 500 mL cord blood



Scientists Helping Scientists™ | [WWW.STEMCELL.COM](http://WWW.STEMCELL.COM)

TOLL FREE PHONE 1 800 667 0322 • PHONE +1 604 877 0713  
 INFO@STEMCELL.COM • TECHSUPPORT@STEMCELL.COM  
 FOR GLOBAL CONTACT DETAILS VISIT OUR WEBSITE

Document #28365 | Version 2\_4\_1

## Description

Enrich untouched human cord blood progenitor cells directly from human cord blood by negative selection.

- Fast and easy-to-use
- Requires no special equipment or training
- Untouched, viable cells

This kit targets non-progenitor cells for removal with antibodies recognizing specific cell surface markers. The RosetteSep™ antibody cocktail crosslinks unwanted cells in human whole blood to multiple red blood cells (RBCs), forming immunorosettes. This increases the density of the unwanted (rosetted) cells, such that they pellet along with the free RBCs when centrifuged over a density gradient medium. Desired cells are never labeled with antibody and are easily collected as a highly enriched population at the interface between the plasma and the density gradient medium. Isolated cells are immediately available for downstream applications such as flow cytometry, cell culture, or DNA/RNA extraction.

NOTE: Disregard Document #28581 provided in the cocktail box for Catalog #15026. These instructions will replace those provided in the cocktail box when using kit #15276.

## Component Descriptions

COMPONENT NAME	COMPONENT #	QUANTITY	STORAGE	SHELF LIFE	FORMAT
RosetteSep™ Human Hematopoietic Progenitor Cell Enrichment Cocktail	15186C	2 x 2 mL	Store at 2 - 8°C. Do not freeze.	Stable until expiry date (EXP) on label.	A combination of monoclonal antibodies in PBS.
HetaSep™	07906	100 mL	Store at 15 - 25°C.	Stable until expiry date (EXP) on label.	6% Hetastarch Solution.

PBS - phosphate-buffered saline

Components may be shipped at room temperature (15 - 25°C) but should be stored as indicated above.

Precipitate may be observed in the cocktail vial but will not affect performance.

## Sample Preparation

For available cord blood products see [www.stemcell.com/primarycells](http://www.stemcell.com/primarycells).

### CORD BLOOD

For optimal performance use whole cord blood collected within the last 24 hours and stored at room temperature (15 - 25°C).

## Recommended Medium

Dulbecco's Phosphate Buffered Saline with 2% Fetal Bovine Serum (Catalog #07905).

## Density Gradient Medium

Lymphoprep™ (Catalog #07801) or other density gradient medium with a density of 1.077 g/mL.

## Directions for Use – RosetteSep™ Protocol

See page 1 for Sample Preparation and Recommended Medium.

Ensure that cord blood sample, HetaSep™, recommended medium, density gradient medium, and centrifuge are all at room temperature (15 - 25°C).

NOTE: This protocol was developed for separations in 50 mL centrifuge tubes but can be adapted for use with blood bags.

**Table 1. Concentrate Leukocytes**

STEP	INSTRUCTIONS	Standard Tube
1	Prepare sample.	Add 2 mL of HetaSep™ per 10 mL cord blood
		Mix well
2	Centrifuge.	50 x g for 5 minutes, <b>brake off</b>
		Retain supernatant plus top 10% of pellet volume and place into a new tube. Discard remainder of pellet.
3	Wash cells.	Top up with recommended medium
4	Centrifuge. Remove and discard supernatant using a pipette.	300 x g for 10 minutes brake low <b>Do not pour off</b> , as pellet is not firmly packed.
5	Resuspend cells in recommended medium.	Add 0.5 mL per 10 mL of <b>original</b> cord blood volume. Final volume should be approximately 10 - 15% of original blood volume.

**Table 2. RosetteSep™ Human Cord Blood Progenitor Enrichment Kit Protocol**

		ROSETTESEP™
STEP	INSTRUCTIONS	Standard Tube
1	Add RosetteSep™ Cocktail to sample.	75 µL/10 mL of <b>original</b> cord blood volume
	Mix and incubate.	RT for 10 minutes
2	Dilute sample with recommended medium and mix gently.	Use twice as much volume of recommended medium
3	Add density gradient medium to required tube.	See Table 3 for volumes and tubes
4	Add diluted sample to the tube containing the density gradient medium. Use at least half as much density gradient medium as the total volume of the diluted sample	Layer diluted sample on density gradient medium, being careful to minimize their mixing
5	Centrifuge.	1200 x g for 20 minutes, <b>brake off</b>
6	Collect enriched cells. * For platelet removal see footnote below.	Harvest enriched cell layer with a pipette and transfer to new tube**
7	Wash enriched cells.	Top up with recommended medium
8	Centrifuge.	300 x g for 10 minutes brake low***
		Discard supernatant
9	Resuspend cells in recommended medium.	The enriched cells are ready for use

RT - room temperature (15 - 25°C)

\* To minimize platelet contamination, remove and discard the top third of the plasma layer before collecting the cells at the density gradient medium:plasma interface. Platelets may also be removed by including an extra wash after step 8, by centrifugation at 120 x g for 10 minutes at room temperature with no brake.

\*\* Sometimes it is difficult to see the cells at the interface. It is recommended to remove some of the density gradient medium along with the pre-enriched cells in order to ensure complete recovery.

\*\*\* The wash step can be done with Ammonium Chloride Solution (Catalog #07800) prior to flow cytometry analysis or if residual RBCs will interfere with subsequent assays.

**Table 3. Recommended Volumes and Tube Sizes**

WHOLE BLOOD VOLUME	RECOMMENDED MEDIUM VOLUME	STANDARD TUBE	
		TUBE SIZE	DENSITY GRADIENT MEDIUM VOLUME
0.5 mL	1 mL	5 mL	0.75 mL
1 mL	2 mL	5 mL	1.5 mL
2 mL	4 mL	14 mL	3 mL
3 mL	6 mL	14 mL	4.5 mL
4 mL	8 mL	50 mL	6 mL
5 mL	10 mL	50 mL	7.5 mL
10 mL	20 mL	50 mL	15 mL

\* Small bubbles may be present in the density gradient medium after pipetting. This will not affect performance.

## Notes and Tips

### CONVERSION OF g TO RPM

To convert g to RPM, use the following formula:

$$RPM = \sqrt{\frac{RCF}{(1.118 \times 10^{-5}) \times (\text{Radius})}}$$

Where:  
 RCF = relative centrifugal force (g)  
 RPM = centrifuge speed in revolutions per minute  
 Radius = radius of rotor in cm

### ASSESSING PURITY

For purity assessment of human cord blood progenitor cells by flow cytometry, use one of the following fluorochrome-conjugated antibodies:

- Anti-Human CD34 Antibody, Clone 581 (Catalog #60013), or
- Anti-Human CD34 Antibody, Clone 8G12 (Catalog #60121)

STEMCELL TECHNOLOGIES INC.'S QUALITY MANAGEMENT SYSTEM IS CERTIFIED TO ISO 13485. PRODUCTS ARE FOR RESEARCH USE ONLY AND NOT INTENDED FOR HUMAN OR ANIMAL DIAGNOSTIC OR THERAPEUTIC USES UNLESS OTHERWISE STATED.

Copyright © 2019 by STEMCELL Technologies Inc. All rights reserved including graphics and images. STEMCELL Technologies & Design, STEMCELL Shield Design, Scientists Helping Scientists, HetaSep, and RosetteSep are trademarks of STEMCELL Technologies Canada Inc. Lymphoprep is a trademark of Alere Technologies. All other trademarks are the property of their respective holders. While STEMCELL has made all reasonable efforts to ensure that the information provided by STEMCELL and its suppliers is correct, it makes no warranties or representations as to the accuracy or completeness of such information.