



EasySep™ Human Pan-DC Pre-Enrichment Kit

Negative Selection
Catalog #19251

For processing 1 x 10⁹ cells



Scientists Helping Scientists™ | WWW.STEMCELL.COM

TOLL FREE PHONE 1 800 667 0322 • PHONE +1 604 877 0713
INFO@STEMCELL.COM • TECHSUPPORT@STEMCELL.COM
FOR GLOBAL CONTACT DETAILS VISIT OUR WEBSITE

Document #29038 | Version 3_2_1

Description

Isolate untouched pre-enriched dendritic cells (including myeloid and plasmacytoid dendritic cells; DCs) from fresh human peripheral blood mononuclear cells (PBMCs) by immunomagnetic negative selection.

- Fast, easy-to-use and column-free
- Up to 80% purity
- Untouched, viable cells

This kit targets non-DCs for removal with antibodies recognizing specific cell surface markers. Unwanted cells are labeled with antibodies and magnetic particles, and separated without columns using an EasySep™ magnet. Desired cells are simply poured off into a new tube. Isolated cells are immediately available for downstream applications such as flow cytometry, culture or DNA/RNA extraction.

Component Descriptions

COMPONENT NAME	COMPONENT #	QUANTITY	STORAGE	SHELF LIFE	FORMAT
EasySep™ Human Pan-DC Pre-Enrichment Cocktail	19251C.1	1 x 1 mL	Store at 2 - 8°C. Do not freeze.	Stable until expiry date (EXP) on label.	A combination of monoclonal antibodies in PBS.
EasySep™ D Magnetic Particles	19250	2 x 1 mL	Store at 2 - 8°C. Do not freeze.	Stable until expiry date (EXP) on label.	A suspension of magnetic particles in TBS.
Anti-Human CD32 (Fc gamma RII) Blocker	14551C	1 x 0.8 mL	Store at 2 - 8°C. Do not freeze.	Stable until expiry date (EXP) on label.	A combination of monoclonal antibodies in PBS.

PBS - phosphate-buffered saline; TBS - TRIS-buffered saline

Components may be shipped at room temperature (15 - 25°C) but should be stored as indicated above.

Sample Preparation

For available fresh and frozen samples, see www.stemcell.com/primarycells.

PERIPHERAL BLOOD

Prepare a PBMC suspension from whole peripheral blood by centrifugation over a density gradient medium (e.g. Lymphoprep™, Catalog #07801). For more rapid PBMC preparation without the need for careful sample layering, use the SepMate™-15 (Catalog #15415) or SepMate™-50 (Catalog #15450) cell isolation tube.

If using previously frozen PBMCs, thaw up to 5 x 10⁸ cells in 50 mL of PBS containing 20% fetal bovine serum (FBS) and 1 mM EDTA. Centrifuge at 200 x g for 10 minutes. Carefully remove and discard the supernatant and resuspend the cells in recommended medium. Centrifuge at 200 x g for 10 minutes and carefully remove and discard the supernatant. DNase I treatment is not recommended when preparing the cells for separation.

After preparation, resuspend cells at 5 x 10⁷ cells/mL in recommended medium.



Recommended Medium

EasySep™ Buffer (Catalog #20144), RoboSep™ Buffer (Catalog #20104), or PBS containing 2% FBS and 1 mM EDTA. Medium should be free of Ca⁺⁺ and Mg⁺⁺.

Directions for Use – Manual EasySep™ Protocols

See page 1 for Sample Preparation and Recommended Medium. Refer to Tables 1 and 2 for detailed instructions regarding the EasySep™ procedure for each magnet.

Table 1. EasySep™ Human Pan- DC Pre-Enrichment Kit Protocol


		EASYSEP™ MAGNETS	
STEP	INSTRUCTIONS	 EasySep™ (Catalog #18000)	 “The Big Easy” (Catalog #18001)
1	Prepare sample at the indicated cell concentration within the volume range.	5 x 10 ⁷ cells/mL 0.25 - 2 mL	5 x 10 ⁷ cells/mL 0.5 - 8.5 mL
	Add sample to required tube.	5 mL (12 x 75 mm) polystyrene round-bottom tube (e.g. Corning Catalog #352058)	14 mL (17 x 100 mm) polystyrene round-bottom tube (e.g. Corning Catalog #352057)
2	Add Blocker to sample.‡	30 µL/mL of sample	30 µL/mL of sample
3	Add Enrichment Cocktail to sample.	50 µL/mL of sample	50 µL/mL of sample
	Mix and incubate.	RT for 30 minutes	RT for 30 minutes
4	Vortex Magnetic Particles.	30 seconds	30 seconds
5	Add Magnetic Particles to sample.	100 µL/mL of sample	100 µL/mL of sample
	Mix and incubate.	RT for 10 minutes	RT for 10 minutes
6	Add recommended medium to top up the sample to the indicated volume. Mix by gently pipetting up and down 2 - 3 times.	Top up to 2.5 mL	<ul style="list-style-type: none"> • Top up to 5 mL for samples < 2 mL • Top up to 10 mL for samples ≥ 2 mL
	Place the tube (without lid) into the magnet and incubate.	RT for 5 minutes	RT for 5 minutes
7	Pick up the magnet, and in one continuous motion invert the magnet and tube,* pouring off the enriched cell suspension into a new tube.	Use a new 5 mL tube	Use a new 14 mL tube
OPTIONAL ADDITIONAL SEPARATION Note: This will increase purity but may decrease recovery. Use of this protocol will require additional Magnetic Particles. Please contact Technical Support to request more.		---	---
8	Add Magnetic Particles to sample.	100 µL/mL of sample	100 µL/mL of sample
	Mix and incubate.	RT for 10 minutes	RT for 10 minutes
9	Remove the tube from the magnet and place the new tube (without lid) into the magnet and incubate for a second separation.	RT for 5 minutes	RT for 5 minutes
10	Pick up the magnet, and in one continuous motion invert the magnet and tube,* pouring the enriched cell suspension into a new tube.	Isolated cells are ready for use	Isolated cells are ready for use

RT - room temperature (15 - 25°C)

‡ Addition of Anti-Human CD32 (Fc gamma RII) Blocker is optional but may improve product performance by preventing non-specific depletion of DCs. Use of the Anti-Human CD32 (Fc gamma RII) Blocker may prevent subsequent attempts at cross-linking CD32 molecules on the surface of enriched cells to trigger signaling through these receptors. It may therefore be desirable to not add the Anti-Human CD32 (Fc gamma RII) Blocker to the cell suspension for such studies.

* Leave the magnet and tube inverted for 2 - 3 seconds, then return upright. Do not shake or blot off any drops that may remain hanging from the mouth of the tube.

Table 2. EasySep™ Human Pan- DC Pre-Enrichment Kit Protocol

		EASYSEP™ MAGNETS
STEP	INSTRUCTIONS	Easy 50 (Catalog #18002) 
1	Prepare sample at the indicated cell concentration within the volume range.	5 x 10 ⁷ cells/mL 1 - 35 mL
	Add sample to required tube.	50 mL conical tube (e.g. Corning Catalog #352070)
2	Add Blocker to sample.‡	30 µL/mL of sample
3	Add Enrichment Cocktail to sample.	50 µL/mL of sample
	Mix and incubate.	RT for 30 minutes
4	Vortex Magnetic Particles.	30 seconds
5	Add Magnetic Particles to sample.	250 µL/mL of sample
	Mix and incubate.	RT for 10 minutes
6	Add recommended medium to top up sample to the indicated volume. Mix by gently pipetting up and down 2 - 3 times.	<ul style="list-style-type: none"> • Top up to 10 mL for samples < 5 mL • Top up to 20 mL for samples ≥ 5 - 10 mL • Top up to 30 mL for samples > 10 - 15 mL • Top up to 40 mL for samples > 15 - 20 mL • Top up to 50 mL for samples > 20 mL
	Place the tube (without lid) into the magnet and incubate.	RT for 5 minutes
7	Carefully pipette** (do not pour) the enriched cell suspension into a new tube.	Use a new 50 mL tube
8	Remove the tube from the magnet and place the new tube (without lid) into the magnet and incubate for a second round of separation.	RT for 5 minutes
9	Carefully pipette** (do not pour) the enriched cell suspension into a new tube	Isolated cells are ready for use

RT - room temperature (15 - 25°C)

‡ Addition of Anti-Human CD32 (Fc gamma RII) Blocker is optional but may improve product performance by preventing non-specific depletion of DCs. Use of the Anti-Human CD32 (Fc gamma RII) Blocker may prevent subsequent attempts at cross-linking CD32 molecules on the surface of enriched cells to trigger signaling through these receptors. It may therefore be desirable to not add the Anti-Human CD32 (Fc gamma RII) Blocker to the cell suspension for such studies.

** Collect the entire supernatant, all at once, into a single pipette.

Directions for Use – Fully Automated RoboSep™ Protocol

See page 1 for Sample Preparation and Recommended Medium. Refer to Table 3 for detailed instructions regarding the RoboSep™ procedure.

Table 3. RoboSep™ Human Pan-DC Pre-Enrichment Kit Protocol

STEP	INSTRUCTIONS	RoboSep™ (Catalog #20000 and #21000)
1	Prepare sample at the indicated cell concentration within the volume range.	5 x 10 ⁷ cells/mL 0.5 - 6.5 mL
	Add sample to required tube.	14 mL (17 x 100 mm) polystyrene round-bottom tube (e.g. Corning Catalog #352057)
2	Add Blocker to sample.	30 µL/mL of sample
3	Select protocol.	Human panDC pre-Enrichment Negative Selection 19251-high recovery
4	Vortex Magnetic Particles .	30 seconds
5	Load the carousel.	Follow on-screen prompts
	Start the protocol.	Press the green "Run" button
6	Unload the carousel when the run is complete.	Isolated cells are ready for use

Notes and Tips


ASSESSING PURITY

For purity assessment of pan-DCs by flow cytometry use an appropriate combination of fluorochrome-conjugated:

- Anti-Human CD14 Antibody, Clone M5E2 (Catalog #60004), and
- Anti-Human CD16 Antibody, Clone 3G8 (Catalog #60041), and
- Anti-Human CD20 Antibody, Clone 2H7 (Catalog #60008), and
- Anti-Human CD56 (NCAM) Antibody, Clone HCD56 (Catalog #60021), and
- Anti-human HLA-DR antibody, and
- Anti-human TCR alpha/beta antibody, and
- Anti-Human CD123 (IL-3Ra) Antibody, Clone 6H6 (Catalog #60110), or anti-human CD11c antibody

DCs are described as either Lineage- (CD14, CD16, CD20, CD56, TCRab), HLA-DR+CD123+ or Lineage- (CD14, CD16, CD20, CD56, TCRab), HLA-DR+CD11c+.

Data



Starting with peripheral blood mononuclear cells, the dendritic cell content (Lin-/HLA-DR+/CD123+ or CD11c+) of the enriched fraction typically ranges from 40 - 80%. In the above example, the final purities of the start and enriched fractions are 1.25% and 72.7%, respectively.

STEMCELL TECHNOLOGIES INC.'S QUALITY MANAGEMENT SYSTEM IS CERTIFIED TO ISO 13485. PRODUCTS ARE FOR RESEARCH USE ONLY AND NOT INTENDED FOR HUMAN OR ANIMAL DIAGNOSTIC OR THERAPEUTIC USES UNLESS OTHERWISE STATED.

Copyright © 2015 by STEMCELL Technologies Inc. All rights reserved including graphics and images. STEMCELL Technologies & Design, STEMCELL Shield Design and Scientists Helping Scientists, EasySep, RoboSep, and SepMate are trademarks of STEMCELL Technologies Inc. Lymphoprep is a trademark of AXIS-SHIELD. All other trademarks are the property of their respective holders. While STEMCELL has made all reasonable efforts to ensure that the information provided by STEMCELL and its suppliers is correct, it makes no warranties or representations as to the accuracy or completeness of such information.