

EasySep™ Human Eosinophil Enrichment Kit

Negative Selection
Catalog #19256

For processing 1 x 10⁹ cells



Scientists Helping Scientists™ | WWW.STEMCELL.COM

TOLL FREE PHONE 1 800 667 0322 • PHONE +1 604 877 0713
INFO@STEMCELL.COM • TECHSUPPORT@STEMCELL.COM
FOR GLOBAL CONTACT DETAILS VISIT OUR WEBSITE

Document #29039 | Version 1_0_3

Description

Isolate untouched and highly purified eosinophils (CD16-CD66b+CD45+) from fresh human peripheral blood polymorphonuclear cells (PMNCs) immunomagnetic negative selection.

- Fast, easy-to-use and column-free
- Up to 99% purity
- Untouched, viable cells

This kit targets non-eosinophils for removal with antibodies recognizing specific cell surface markers. Unwanted cells are labeled with antibodies and magnetic particles, and separated without columns using an EasySep™ magnet. Desired cells are simply poured off into a new tube. Isolated cells are immediately available for downstream applications such as flow cytometry, culture, or DNA/RNA extraction.

Component Descriptions

COMPONENT NAME	COMPONENT #	QUANTITY	STORAGE	SHELF LIFE	FORMAT
EasySep™ Human Eosinophil Enrichment Cocktail	19256C	1 x 1 mL	Store at 2 - 8°C. Do not freeze.	Stable until expiry date (EXP) on label.	A combination of monoclonal antibodies in PBS.
EasySep™ Magnetic Nanoparticles	19150.1	3 x 1 mL	Store at 2 - 8°C. Do not freeze.	Stable until expiry date (EXP) on label.	A suspension of magnetic particles in water.

PBS - phosphate-buffered saline

Components may be shipped at room temperature (15 - 25°C) but should be stored as indicated above.

Sample Preparation

Important: Do not use dextran sedimentation to prepare cells.

WHOLE BLOOD USING RED BLOOD CELL (RBC) LYSIS (preferred for slightly higher purity)

1. Collect whole blood in a blood collection tube containing anticoagulant.
2. Carefully perform a standard density gradient separation (e.g. using Lymphoprep™; Catalog #07801). Do not use SepMate™.
3. Remove and discard the plasma layer, the band of mononuclear cells and the density gradient medium leaving the RBC pellet intact.
4. Add Ammonium Chloride Solution (Catalog #07800) to the RBC pellet and mix well.
5. Incubate on ice for 10 minutes then centrifuge at 300 x g for 8 minutes.
6. Discard supernatant and wash pellet with cold recommended medium, centrifuging at 250 x g for 10 minutes.
7. Discard supernatant and resuspend cells at 5 x 10⁷ cells/mL in recommended medium.

WHOLE BLOOD USING HETASEP™ RBC SEDIMENTATION (preferred for faster, lysis-free sample processing)

1. Collect whole blood in a blood collection tube containing anticoagulant.
2. Add 1 part HetaSep™ (Catalog #07906) to 5 parts whole blood and mix well. Use the minimum-sized tube for the total volume of HetaSep™ : blood sample. A 14 mL tube is the maximum size recommended for optimal leukocyte recovery.
3. Centrifuge sample at 50 x g for 5 minutes at room temperature (15 - 25°C) with the brake off.
4. Remove tube from centrifuge and let sit undisturbed (maximum 10 minutes) until the RBC : plasma interface is approximately 40% of the total volume.
5. Harvest the leukocyte-rich plasma (everything above the RBC fraction) into a 50 mL tube and add 4 parts recommended medium to 1 part harvested cells/plasma.
6. Centrifuge at 500 x g for 10 minutes at room temperature with the brake on low.
7. Discard supernatant and wash pellet to remove excess platelets, centrifuging at 120 x g for 10 minutes at room temperature with the brake off.
8. Discard supernatant and resuspend cells at 5 x 10⁷ cells/mL in recommended medium.

Recommended Medium

EasySep™ Buffer (Catalog #20144), RoboSep™ Buffer (Catalog #20104), or PBS containing 2% fetal bovine serum (FBS) and 1 mM EDTA. Medium should be free of Ca⁺⁺ and Mg⁺⁺.

Directions for Use – Manual EasySep™ Protocols

See page 1 for Sample Preparation and Recommended Medium. Refer to Table 1 for detailed instructions regarding the EasySep™ procedure for each magnet.

Table 1. EasySep™ Human Eosinophil Enrichment Kit Protocol

		EASYSEP™ MAGNETS	
STEP	INSTRUCTIONS	 EasySep™ (Catalog #18000)	“The Big Easy” (Catalog #18001) 
1	Prepare sample at the indicated cell concentration within the volume range.	5 x 10 ⁷ cells/mL 1 - 2 mL	5 x 10 ⁷ cells/mL 1 - 6.5 mL
	Add sample to required tube.	5 mL (12 x 75 mm) polystyrene round-bottom tube (e.g. Corning Catalog #352058)	14 mL (17 x 100 mm) polystyrene round-bottom tube (e.g. Corning Catalog #352057)
2	Add Selection Cocktail to sample.	50 µL/mL of sample	50 µL/mL of sample
	Mix and incubate.	2 - 8°C for 10 minutes	2 - 8°C for 10 minutes
3	Mix Magnetic Particles. NOTE: Particles should appear evenly dispersed.	Pipette up and down more than 5 times	Pipette up and down more than 5 times
4	Add Magnetic Particles to sample.	100 µL/mL of sample	100 µL/mL of sample
	Mix and incubate.	2 - 8°C for 10 minutes	2 - 8°C for 10 minutes
5	Add recommended medium to top up the sample to the indicated volume. Mix by gently pipetting up and down 2 - 3 times.	Top up to 2.5 mL	<ul style="list-style-type: none"> • Top up to 5 mL for samples < 2 mL • Top up to 10 mL for samples ≥ 2 mL
	Place the tube (without lid) into the magnet and incubate.	RT for 5 minutes	RT for 10 minutes
6	Pick up the magnet, and in one continuous motion invert the magnet and tube,* pouring the enriched cell suspension into a new tube.	Use a new 5 mL tube Isolated cells are ready for use	Use a new 14 mL tube Isolated cells are ready for use
OPTIONAL ADDITIONAL SEPARATION NOTE: This will improve purity but may reduce recovery.		---	---
7	Remove the tube from the magnet and place the new tube (without lid) into the magnet and incubate for a second separation.	RT for 5 minutes	RT for 10 minutes
8	Pick up the magnet, and in one continuous motion invert the magnet and tube,* pouring the enriched cell suspension into a new tube.	Isolated cells are ready for use	Isolated cells are ready for use

RT - room temperature (15 - 25°C)

* Leave the magnet and tube inverted for 2 - 3 seconds, then return upright. Do not shake or blot off any drops that may remain hanging from the mouth of the tube.

Directions for Use – Fully Automated RoboSep™ Protocol

See page 1 for Sample Preparation and Recommended Medium. Refer to Table 2 for detailed instructions regarding the RoboSep™ procedure.

Table 2. RoboSep™ Human Eosinophil Enrichment Kit Protocol

STEP	INSTRUCTIONS	RoboSep™ (Catalog #20000 and #21000)
1	Prepare sample at the indicated cell concentration within the volume range.	5 x 10 ⁷ cells/mL 1 - 6.5 mL
	Add sample to required tube.	14 mL (17 x 100 mm) polystyrene round-bottom tube (e.g. Corning Catalog #352057)
2	Select protocol.	Human Eosinophil Negative Selection 19256-high purity
3	Mix Magnetic Particles. NOTE: Particles should appear evenly dispersed.	Pipette up and down more than 5 times
4	Load the carousel.	Follow on-screen prompts
	Start the protocol.	Press the green "Run" button
5	Unload the carousel when the run is complete.	Isolated cells are ready for use

Notes and Tips

ASSESSING PURITY


For purity assessment of eosinophils by flow cytometry use the following fluorochrome-conjugated antibody clones:

- Anti-Human CD16 Antibody, Clone 3G8 (Catalog #60041), and
- Anti-Human CD66b Antibody, Clone G10F5 (Catalog #60086)
- Anti-Human CD45 Antibody, Clone HI30 (Catalog #60018; optional)

Eosinophils are CD16-CD66b+ and are low in forward scatter but high in side scatter.

Alternatively, purity may be assessed by performing a cytospin on the enriched cells followed by Wright's or May-Grunwald staining (e.g. Sigma-Aldrich Catalog #W0625 or #205435, respectively).

Data



Starting with whole blood prepared using RBC lysis or HetaSep™, the eosinophil content (CD16-CD66b+CD45+) of the enriched fraction typically ranges from 86 - 99%. In the above examples, the purities of the start and final enriched fractions prepared using RBC lysis are 1.2% and 98% (gated on CD45+), respectively and prepared using HetaSep™ are 1.2% and 96% (gated on CD45+), respectively.

STEMCELL TECHNOLOGIES INC.'S QUALITY MANAGEMENT SYSTEM IS CERTIFIED TO ISO 13485. PRODUCTS ARE FOR RESEARCH USE ONLY AND NOT INTENDED FOR HUMAN OR ANIMAL DIAGNOSTIC OR THERAPEUTIC USES UNLESS OTHERWISE STATED.

Copyright © 2016 by STEMCELL Technologies Inc. All rights reserved including graphics and images. STEMCELL Technologies & Design, STEMCELL Shield Design, Scientists Helping Scientists, HetaSep, EasySep, RoboSep, and SepMate are trademarks of STEMCELL Technologies Inc. Lymphoprep is a trademark of AXIS-SHIELD. All other trademarks are the property of their respective holders. While STEMCELL has made all reasonable efforts to ensure that the information provided by STEMCELL and its suppliers is correct, it makes no warranties or representations as to the accuracy or completeness of such information.