

# ClonaCell™-HY Medium E



**Hybridoma growth medium with hypoxanthine and thymidine (serum-containing)**

Catalog # 03805      500 mL

Scientists Helping Scientists™ | [WWW.STEMCELL.COM](http://WWW.STEMCELL.COM)

TOLL FREE PHONE 1 800 667 0322 • PHONE +1 604 877 0713

[INFO@STEMCELL.COM](mailto:INFO@STEMCELL.COM) • [TECHSUPPORT@STEMCELL.COM](mailto:TECHSUPPORT@STEMCELL.COM)

FOR GLOBAL CONTACT DETAILS VISIT OUR WEBSITE

## Product Description

ClonaCell™-HY Medium E is a serum-containing, nutritionally rich liquid medium optimized for hybridoma expansion after hypoxanthine, aminopterin, thymidine (HAT) selection. The medium contains hypoxanthine and thymidine (HT) and is used to wean hybridomas off aminopterin used during the selection process. This medium has been verified for use in mouse and rat hybridoma development and monoclonal antibody production, and reportedly is compatible for production, cloning, and expansion of hybridomas using lymphocytes from a variety of host animals including human, mouse, rat, and hamster. ClonaCell™-HY Medium E has proven effective for subcloning and stabilization of unstable hybridoma cell lines.

## Properties

- Storage:** Store at -20°C.
- Shelf Life:** Stable until expiry date (EXP) on label.
- Contains:**
- DMEM
  - Pre-selected serum
  - Hypoxanthine and thymidine (HT)
  - Gentamicin
  - 2-Mercaptoethanol
  - Phenol red
  - L-Glutamine and other supplements
  - Other ingredients

## Handling / Directions For Use

1. Thaw ClonaCell™-HY Medium E at room temperature (15 - 25°C) or overnight at 2 - 8°C. Mix well.  
NOTE: Do not thaw ClonaCell™-HY Medium E in a 37°C water bath.
2. If ClonaCell™-HY Medium E is not used immediately, store at 2 - 8°C for up to 2 weeks. Alternatively, aliquot and store at -20°C until expiry date as indicated on the label.

For further information, refer to the Technical Manual: ClonaCell™-HY: A Complete Workflow for Hybridoma Generation (Document #28411), available at [www.stemcell.com](http://www.stemcell.com) or contact us to request a copy.

## References

- Chen ZC et al. (2000) Genes coding evolutionary novel anti-carbohydrate antibodies: Studies on anti-Gal production in alpha 1,3galactosyltransferase knock out mice. *Mol Immunol* 37(8): 455–66.
- Flyak AI et al. (2015) Mechanism of human antibody-mediated neutralization of Marburg virus. *Cell* 160(5): 893–903.
- García-Barreno B et al. (2014) Characterization of an enhanced antigenic change in the pandemic 2009 H1N1 influenza virus haemagglutinin. *J Gen Virol* 95(5): 1033–42.
- Kaabinejadian S et al. (2016) Immunodominant West Nile virus T cell epitopes are fewer in number and fashionably late. *J Immunol* 196(10): 4263–73.
- Loveless BC et al. (2011) Structural characterization and epitope mapping of the glutamic acid/alanine-rich protein from *Trypanosoma congolense*: Defining assembly on the parasite cell surface. *J Biol Chem* 286(23): 20658–65.
- Okai S et al. (2016) High-affinity monoclonal IgA regulates gut microbiota and prevents colitis in mice. *Nat Microbiol* 1(9): 16103.
- Sapparapu G et al. (2016) Neutralizing human antibodies prevent Zika virus replication and fetal disease in mice. *Nature* 540(7633): 443–7.

- Smith SA et al. (2012) Persistence of circulating memory B cell clones with potential for Dengue virus disease enhancement for decades following infection. *J Virol* 86(5): 2665–75.
- Spanier JA et al. (2016) Efficient generation of monoclonal antibodies against peptide in the context of MHCII using magnetic enrichment. *Nat Commun* 7: 11804.
- Tan GS et al. (2014) Characterization of a broadly neutralizing monoclonal antibody that targets the fusion domain of Group 2 influenza A virus hemagglutinin. *J Virol* 88(23): 13580–92.
- Wittman VP et al. (2006) Antibody targeting to a class I MHC-peptide epitope promotes tumor cell death. *J Immunol* 177(6): 4187–95.
- Yew CW & Tan YJ. (2016) Generation of mouse monoclonal antibodies specific to chikungunya virus using ClonaCell-HY Hybridoma Cloning Kit. *Methods Mol Biol* 1426: 225–33.

STEMCELL TECHNOLOGIES INC.'S QUALITY MANAGEMENT SYSTEM IS CERTIFIED TO ISO 13485. PRODUCTS ARE FOR RESEARCH USE ONLY AND NOT INTENDED FOR HUMAN OR ANIMAL DIAGNOSTIC OR THERAPEUTIC USES UNLESS OTHERWISE STATED.

Copyright © 2019 by STEMCELL Technologies Inc. All rights reserved including graphics and images. STEMCELL Technologies & Design, STEMCELL Shield Design, Scientists Helping Scientists, and ClonaCell are trademarks of STEMCELL Technologies Canada Inc. All other trademarks are the property of their respective holders. While STEMCELL has made all reasonable efforts to ensure that the information provided by STEMCELL and its suppliers is correct, it makes no warranties or representations as to the accuracy or completeness of such information.