

MethoCult™ SF H4236

Serum-Free Methylcellulose-Based Medium Without Cytokines for Human Cells

Catalog #04236

80 mL



Scientists Helping Scientists™ | WWW.STEMCELL.COM

TOLL FREE PHONE 1 800 667 0322 • PHONE +1 604 877 0713

INFO@STEMCELL.COM • TECHSUPPORT@STEMCELL.COM

FOR GLOBAL CONTACT DETAILS VISIT OUR WEBSITE

FOR RESEARCH USE ONLY. NOT INTENDED FOR HUMAN OR ANIMAL DIAGNOSTIC OR THERAPEUTIC USES.

Product Description

Incomplete Methylcellulose-Based Medium for Colony-Forming Unit (CFU) Assays for Human Cells

MethoCult™ SF H4236 is recommended as a base medium for the culture of human cells in defined serum-free conditions. It is optimized for the detection and quantification of hematopoietic progenitor cells in human bone marrow (BM), mobilized peripheral blood (MPB), peripheral blood (PB) and cord blood (CB) samples using CFU assays. This formulation allows for the addition of an exogenous source of erythropoietin and other cytokines, and is ideal for testing cytokine effects where the presence of fetal bovine serum (FBS) is not desired.

Properties

Storage: Store at -20°C.

Shelf Life: Stable until expiry date (EXP) on label.

Contains:

- Methylcellulose in Iscove's MDM
- Bovine serum albumin
- Recombinant human insulin
- Human transferrin (iron-saturated)
- 2-Mercaptoethanol
- Supplements

Human transferrin used in this product has been derived from human plasma. Donors have been tested and found negative for hepatitis B surface antigen (HBsAg) and HIV-1 antibodies and/or HIV-1 antigen. However, this product should be considered potentially infectious and treated in accordance with universal handling precautions.

Handling/Directions for Use

NOTE: If product is received partially thawed, place immediately at -20°C or thaw and aliquot as described below.

PREPARATION OF COMPLETE METHOCULT™ SF H4236 MEDIUM

MethoCult™ SF H4236 base medium does not contain EPO or other cytokines. These can be added directly to the bottle or to each tube after aliquoting. Refer to Table 1 for volumes required to prepare complete MethoCult™ SF H4236 medium per bottle or per tube. The 4:1 (v:v) ratio of MethoCult™ to other components in the liquid medium (e.g. cytokines) gives the correct viscosity to ensure optimal CFU growth and morphology.

Use sterile techniques to prepare complete MethoCult™ SF H4236 medium (MethoCult™ SF H4236 base medium + desired components).

NOTE: Do not use pipettes to dispense methylcellulose as the volume dispensed will not be accurate. Syringes and large-bore blunt-end needles should be used for accurate dispensing of viscous methylcellulose medium and to prevent needle-stick injuries.

A. TO PREPARE 100 mL BOTTLE

1. Thaw 80 mL bottle of MethoCult™ at room temperature (15 - 25°C) or overnight at 2 - 8°C. Do not thaw MethoCult™ at 37°C.
2. Prepare desired components in Iscove's MDM (IMDM) with 25 mM HEPES (Catalog #36150) in a volume of 20 mL and add to MethoCult™ (total volume of 100 mL).
3. Shake vigorously for 1 - 2 minutes and then let stand for at least 5 minutes to allow bubbles to rise to the top before aliquoting.
4. Using a 3 or 6 mL luer lock syringe attached to a 16 gauge Blunt-End Needle (Catalog #28110), aliquot 3 mL per tube for 1.1 mL duplicate cultures or 4 mL per tube for 1.1 mL triplicate cultures. Complete MethoCult™ medium is now ready for use.

B. TO PREPARE INDIVIDUAL TUBES

1. Thaw 80 mL bottle of MethoCult™ at room temperature (15 - 25°C) or overnight at 2 - 8°C. Do not thaw MethoCult™ at 37°C.
2. Shake vigorously for 1 - 2 minutes and then let stand for at least 5 minutes to allow bubbles to rise to the top before aliquoting.
3. Using a 3 or 6 mL luer lock syringe attached to a 16 gauge Blunt-End Needle (Catalog #28110), aliquot MethoCult™ SF H4236 base medium into tubes (refer to Table 1 for required volumes).
NOTE: Before adding components, tubes of incomplete MethoCult™ medium may be stored at -20°C until expiry date as indicated on label. After thawing aliquoted tubes, add desired components and mix well.
4. Add desired growth factors, supplements, and IMDM with 25 mM Hepes (Catalog #36150) to tubes of MethoCult™ SF H4236 (see Table 1 for required volumes).
5. Vortex tubes to mix well. Complete MethoCult™ medium is now ready for use.
6. Aliquot any remaining MethoCult™ SF H4236 base medium for duplicate or triplicate cultures (see Table 1 for required volumes), store at -20°C, then add desired components after thawing. Mix well before use.

Table 1: Volumes Required for Preparation of Complete MethoCult™ SF H4236 Medium

COMPONENT	PER BOTTLE	PER TUBE (duplicate 1.1 ml cultures)	PER TUBE (triplicate 1.1 ml cultures)
MethoCult™ SF H4236	80 mL	2.4 mL	3.2 mL
IMDM with cytokines*	20 mL	0.6 mL	0.8 mL
TOTAL VOLUME	100 mL	3.0 mL	4.0 mL

*For a complete list of available cytokines, refer to our website at www.stemcell.com.

For recommended cell plating concentrations, set-up of human CFU assays, and counting and classification of human colonies, refer to the Technical Manual: Human Colony-Forming Unit Assays Using MethoCult™ (Document #28404), available on our website at www.stemcell.com or contact us to request a copy.

Notes and Tips

STEMCELL Technologies recommends the use of human Human LDL (Catalog #02698) as a culture supplement. It has been prescreened for the culture, expansion, and colony assay of human hematopoietic and non-hematopoietic cells in serum-free culture media. It promotes the proliferation and survival of human hematopoietic and other progenitor cells in culture, resulting in increased cell output in expansion cultures and increased colony numbers and/or colony size in colony assays.

References

- Eaves CJ, Eaves AC: Anatomy and physiology of hematopoiesis. In: Childhood Leukemias (2nd Edition) (CH Pui, ed.), Cambridge University Press, Cambridge, pp 69-105, 2006
- Eaves C, Lambie K: Atlas of Human Hematopoietic Colonies, STEMCELL Technologies Inc., Vancouver, 1995 (Catalog #28700)
- Wognum B, Yuan N, Lai B, & Miller CL: Colony Forming Cell Assays for Human Hematopoietic Progenitor Cells. In: Basic Cell Culture Protocols (CD Helgason, CL Miller, eds.), Humana Press Inc., Totowa, New Jersey, pp 267-283, 2013
- Nissen-Druey C, Tichelli A, Meyer-Monard S: Human Hematopoietic Colonies in Health and Disease, S. Karger Medical and Scientific Publishers, Basel, 2005. Reprint of Acta Haematol 113: 5-96, 2005 (Catalog #28760)
- Atlas of Hematopoietic Colonies from Cord Blood, STEMCELL Technologies Inc., Vancouver, 2010

STEMCELL TECHNOLOGIES INC.'S QUALITY MANAGEMENT SYSTEM IS CERTIFIED TO ISO 13485 MEDICAL DEVICE STANDARDS.

Copyright ©2014 by STEMCELL Technologies Inc. All rights reserved including graphics and images. STEMCELL Technologies & Design, STEMCELL Shield Design, Scientists Helping Scientists, and MethoCult are trademarks of STEMCELL Technologies Inc. While STEMCELL has made all reasonable efforts to ensure that the information provided by STEMCELL and its suppliers is correct, it makes no warranties or representations as to the accuracy or completeness of such information.