

MethoCult™ H4330

Methylcellulose-Based Medium With EPO for Human Cells

Catalog #04330

90 mL



Scientists Helping Scientists™ | WWW.STEMCELL.COM

TOLL FREE PHONE 1 800 667 0322 • PHONE +1 604 877 0713

INFO@STEMCELL.COM • TECHSUPPORT@STEMCELL.COM

FOR GLOBAL CONTACT DETAILS VISIT OUR WEBSITE

FOR RESEARCH USE ONLY. NOT INTENDED FOR HUMAN OR ANIMAL DIAGNOSTIC OR THERAPEUTIC USES.

Product Description

Incomplete Methylcellulose-Based Medium for Colony-Forming Unit (CFU) Assays for Human Cells

MethoCult™ H4330 is recommended as a base medium to detect and quantify hematopoietic progenitors in human bone marrow (BM), mobilized peripheral blood (MPB), peripheral blood (PB) and cord blood (CB) samples using CFU assays. This formulation allows for the addition of an exogenous source of cytokines to meet the specific requirements of investigators. It can also be used without the addition of exogenous cytokines for the assessment of colony stimulating factor activity/burst promoting activity in unknown samples.

Properties

Storage: Store at -20°C.

Shelf Life: Stable until expiry date (EXP) on label.

Contains:

- Methylcellulose in Iscove's MDM
- Fetal bovine serum
- Bovine serum albumin
- 2-Mercaptoethanol
- Recombinant human erythropoietin (EPO)
- Supplements

Directions for Use

NOTE: If product is received partially thawed, place immediately at -20°C or thaw and aliquot as described below.

PREPARATION OF COMPLETE METHOCULT™ H4330 MEDIUM

MethoCult™ H4330 does not contain cytokines other than EPO. Cytokines can be added directly to the bottle or to each tube after aliquoting. Refer to Table 1 for volumes required to prepare complete MethoCult™ H4330 medium per bottle or per tube. The ratio (v:v) of MethoCult™ to other components in liquid medium (e.g. cytokines) is important for viscosity, which ensures optimal CFU growth and morphology.

Use sterile techniques to prepare complete MethoCult™ H4330 medium (MethoCult™ H4330 base medium + desired components).

NOTE: Do not use pipettes to dispense methylcellulose as the volume dispensed will not be accurate. Syringes and large bore blunt-end needles should be used for accurate dispensing of viscous methylcellulose medium and to prevent needle-stick injuries.

A. TO PREPARE 100 ML BOTTLE

1. Thaw 90 mL bottle at room temperature (15 - 25°C) or overnight at 2 - 8°C. Do not thaw MethoCult™ at 37°C.
2. Prepare desired growth factors, supplements, and Iscove's MDM (IMDM) with 25mM HEPES (Catalog #36150) in a volume of 10 mL and add to MethoCult™ (for a total volume of 100 mL). Shake vigorously for 1 - 2 minutes and then let stand for at least 5 minutes to allow bubbles to rise to the top before aliquoting.
3. Using a 3 or 6 mL luer lock syringe attached to a 16 gauge Blunt-End Needle (Catalog #28110), aliquot 3 mL per tube for 1.1 mL duplicate cultures or 4 mL per tube for 1.1 mL triplicate cultures. Complete MethoCult™ medium is now ready for use.

B. TO PREPARE INDIVIDUAL TUBES

1. Thaw 90 mL bottle at room temperature (15 - 25°C) or overnight at 2 - 8°C. Do not thaw MethoCult™ at 37°C.
2. Shake vigorously for 1 - 2 minutes and then let stand for at least 5 minutes to allow bubbles to rise to the top before aliquoting.

3. Using a 3 or 6 mL luer lock syringe attached to a 16 gauge Blunt-End Needle (Catalog #28110), aliquot 2.7 mL per tube for 1.1 mL duplicate cultures or 3.6 mL per tube for 1.1 mL triplicate cultures.
NOTE: Before adding components, tubes of incomplete MethoCult™ medium may be stored at -20°C until expiry date as indicated on label. After thawing aliquoted tubes, add desired components (see step 4) and mix well.
4. Add desired growth factors, supplements, and Iscove's MDM (IMDM) with 25mM HEPES (Catalog #36150) to tubes of MethoCult™ H4330 (see Table 1 for required volumes).
5. Vortex tubes to mix well. Complete MethoCult™ medium is now ready for use.
6. Aliquot any remaining MethoCult™ H4330 base medium for duplicate or triplicate cultures (see Table 1 for required volumes), store at -20°C, then add desired components after thawing. Mix well before use.

Table 1: Volumes Required for Preparation of Complete MethoCult™ H4330 Medium

COMPONENT	PER BOTTLE	PER TUBE (duplicate 1.1 ml cultures)	PER TUBE (triplicate 1.1 ml cultures)
MethoCult™ H4330	90 mL	2.7 mL	3.6 mL
IMDM with cytokines*	10 mL	0.3 mL	0.4 mL
TOTAL VOLUME	100 mL	3.0 mL	4.0 mL

*For a complete list of available cytokines, refer to our website at www.stemcell.com.

For recommended cell plating concentrations, set-up of human CFU assays, and counting and classification of colonies, refer to the Technical Manual: Human Colony-Forming Cell Assays Using MethoCult™ (Document #28404), available on our website at www.stemcell.com, or contact us to request a copy.

References

- Eaves CJ, Eaves AC: Anatomy and physiology of hematopoiesis. In: Childhood Leukemias (2nd Edition) (CH Pui, ed.), Cambridge University Press, Cambridge, pp 69-105, 2006
- Eaves C, Lambie K: Atlas of Human Hematopoietic Colonies, STEMCELL Technologies Inc., Vancouver, 1995 (Document #28700)
- Wognum B, Yuan N, Lai B, & Miller CL: Colony Forming Cell Assays for Human Hematopoietic Progenitor Cells. In: Basic Cell Culture Protocols (CD Helgason, CL Miller, eds.), Humana Press Inc., Clifton, New Jersey, p267-83, 2013.
- Nissen-Druey C, Tichelli A, Meyer-Monard S: Human Hematopoietic Colonies in Health and Disease, S. Karger Medical and Scientific Publishers, Basel, 2005. Reprint of Acta Haematol 113: 5-96, 2005 (Document #28760)
- Atlas of Hematopoietic Colonies from Cord Blood, STEMCELL Technologies Inc., Vancouver, 2010

STEMCELL TECHNOLOGIES INC.'S QUALITY MANAGEMENT SYSTEM IS CERTIFIED TO ISO 13485 MEDICAL DEVICE STANDARDS.

Copyright ©2014 by STEMCELL Technologies Inc. All rights reserved including graphics and images. STEMCELL Technologies & Design, STEMCELL Shield Design, Scientists Helping Scientists, and MethoCult are trademarks of STEMCELL Technologies Inc. While STEMCELL has made all reasonable efforts to ensure that the information provided by STEMCELL and its suppliers is correct, it makes no warranties or representations as to the accuracy or completeness of such information.