# STEMdiff<sup>™</sup> Definitive Endoderm Kit

Defined animal component-free medium for the differentiation of human ES and iPS cells to definitive endoderm

| Catalog #05110 | 1 Kit |
|----------------|-------|
| Catalog #05115 | 1 Kit |



Scientists Helping Scientists<sup>™</sup> | WWW.STEMCELL.COM

TOLL FREE PHONE 1 800 667 0322 • PHONE +1 604 877 0713 INFO@STEMCELL.COM • TECHSUPPORT@STEMCELL.COM FOR GLOBAL CONTACT DETAILS VISIT OUR WEBSITE

| 0.0  |   |    |
|------|---|----|
| 0.0  |   |    |
| 12.3 |   |    |
| 1.5  |   |    |
| 10   | $ ( \circ $ | ij |
| 1    |   | 1  |

### **Product Description**

STEMdiff<sup>TM</sup> Definitive Endoderm Kit is a serum-free and animal component-free combination of a basal medium and supplements for the differentiation of human embryonic stem (ES) cells and induced pluripotent stem (iPS) cells to definitive endoderm. Cells differentiated to definitive endoderm using this kit can be used to generate more specified cells of endodermal lineage, including hepatocytes and pancreatic progenitors.

The purity of definitive endoderm cells (CXCR4+SOX17+ or CXCR4+c-Kit+) obtained with the STEMdiff™ Definitive Endoderm Kit is typically in the range of 75 - 99%.

STEMdiff™ Definitive Endoderm Kit (Catalog #05110) and STEMdiff™ Definitive Endoderm Kit (TeSR™-E8™ Optimized; Catalog #05115) have been optimized for the differentiation of human ES and iPS cells cultured in mTeSR™1 and TeSR™, respectively.

# Product Information

| PRODUCT NAME   | CATALOG # | SIZE  | COMPONENTS   |
|--|-----------|-------|--|
| STEMdiff™ Definitive Endoderm Kit                          | 05110     | 1 Kit | <ul> <li>STEMdiff™ Definitive Endoderm Basal Medium (100 mL)</li> <li>STEMdiff™ Definitive Endoderm Supplement A (100X; 0.35 mL)</li> <li>STEMdiff™ Definitive Endoderm Supplement B (100X; 1.1 mL)</li> </ul>   |
| STEMdiff™ Definitive Endoderm Kit<br>(TeSR™-E8™ Optimized) | 05115     | 1 Kit | <ul> <li>STEMdiff<sup>™</sup> Definitive Endoderm Basal Medium (100 mL)</li> <li>STEMdiff<sup>™</sup> Definitive Endoderm Supplement A (100X; 0.35 mL)</li> <li>STEMdiff<sup>™</sup> Definitive Endoderm Supplement B (100X; 1.1 mL)</li> <li>STEMdiff<sup>™</sup> Definitive Endoderm TeSR<sup>™</sup>-E8<sup>™</sup> Supplement (20X; 7 mL)</li> </ul> |

# **Component Storage and Stability**

The following components are sold as part of the STEMdiff™ Definitive Endoderm Kits (see Product Information) and are not available for individual sale.

| COMPONENT NAME  | COMPONENT # | STORAGE           | SHELF LIFE  |
|---|-------------|-------------------|---|
| STEMdiff™ Definitive Endoderm<br>Basal Medium               | 05111       | Store at -20°C.   | Stable until expiry date (EXP) on label.                      |
| STEMdiff™ Definitive Endoderm<br>Supplement A (100X)        | 05112       | Store at -20°C.   | Stable for 12 months from date of manufacture (MFG) on label. |
| STEMdiff™ Definitive Endoderm<br>Supplement B (100X)        | 05113       | Store at -20°C.   | Stable for 12 months from date of manufacture (MFG) on label. |
| STEMdiff™ Definitive Endoderm<br>TeSR™-E8™ Supplement (20X) | 05116       | Store at 2 - 8°C. | Stable for 12 months from date of manufacture (MFG) on label. |



## Handling Frozen Components

### 05111 STEMdiff™ Definitive Endoderm Basal Medium

Thaw entire bottle at room temperature (15 - 25°C) or overnight at 2 - 8°C, and mix thoroughly. Once thawed, use immediately or store at 2 - 8°C for up to 2 months. Alternatively, aliquot and store at -20°C until the expiry date as indicated on the label. After thawing the aliquots, use immediately or store at 2 - 8°C for up to 2 weeks. Do not re-freeze.

### 05112 STEMdiff™ Definitive Endoderm Supplement A (100X) OR 05113 STEMdiff™ Definitive Endoderm Supplement B (100X)

• Thaw on ice and mix thoroughly. Once thawed, use immediately or aliquot and store at -20°C for up to 12 months from the date of manufacture as indicated on the label. After thawing the aliquots, use immediately. Do not re-freeze.

| Materials I | Required | But Not | Included |
|-------------|----------|---------|----------|
|-------------|----------|---------|----------|

| PRODUCT NAME  | CATALOG #                     |
|---|-------------------------------|
| mTeSR™1<br>OR<br>TeSR™-E8™  | 05850<br>OR<br>05940          |
| Corning® Matrigel® hESC-qualified matrix<br>OR<br>Vitronectin XF™ | Corning 354277<br>OR<br>07180 |
| DMEM/F-12 with 15 mM HEPES  | 36254                         |
| Gentle Cell Dissociation Reagent                                  | 07174                         |
| D-PBS (Without Ca++ and Mg++)                                     | 37350                         |
| Y-27632   | 72302                         |

# Schematic of STEMdiff<sup>™</sup> Definitive Endoderm Kit Procedure



# **Directions for Use**

Please read the entire protocol before proceeding.

NOTE: For complete instructions on coating plates with Corning® Matrigel® or Vitronectin XF<sup>™</sup>, and maintaining high quality human ES and iPS cells for use in differentiation, please refer to the Technical Manual: Maintenance of Human Pluripotent Stem Cells in mTeSR<sup>™</sup>1 (Document #29106) or TeSR<sup>™</sup>-E8<sup>™</sup> (Document #29267) available on our website at www.stemcell.com, or contact us to request a copy. Matrix-coated plates should be prepared in advance and be brought to room temperature (15 - 25°C) for at least 30 minutes prior to use.

Use sterile techniques when performing the following protocols. The following are instructions for use with 6-well plates. Indicated volumes are for a single well. If using alternative cultureware, adjust volumes accordingly.



### 1. Passaging Cells for Definitive Endoderm Induction

For optimal product performance, passage human ES or iPS cells using the specific passaging protocols for cells cultured in mTeSR™1 or TeSR™-E8™ as outlined in this section, before proceeding with differentiation to definitive endoderm (section 2).

NOTE: Human ES and iPS cells are ready for passage when cultures are approximately 70% confluent.

#### <u>mTeSR™1 Cultures</u>

This protocol is specific to human ES and iPS cells cultured in mTeSR™1 medium.

- On Day 0, warm (15 25°C) sufficient volumes of mTeSR™1, DMEM/F-12, and Gentle Cell Dissociation Reagent for passaging. Prepare Single-Cell Passaging Medium by adding Y-27632 to mTeSR™1 to reach a final concentration of 10 µM.
- 2. Wash the well to be passaged with 1 mL of D-PBS (Without Ca++ and Mg++).
- 3. Aspirate wash medium and add 1 mL of Gentle Cell Dissociation Reagent.
- 4. Incubate at 37°C for 8 10 minutes.
- 5. Dislodge cells by pipetting up and down 1 3 times using a pipette with a P1000 tip. Ensure all remaining cell aggregates are broken up into single cells.
- 6. Immediately transfer cells to a tube containing an equal volume of DMEM/F-12. Wash the well once with 1 mL of DMEM/F-12 to collect any remaining cells and transfer to the tube. Centrifuge the tube at 300 x g for 5 minutes.
- 7. Resuspend cells in 1 mL of Single-Cell Passaging Medium and count the number of live cells using a hemocytometer.
- 8. Plate cells at a density of 2.1 x 10^5 cells/cm<sup>2</sup> (i.e. 2 x 10^6 cells/well) onto pre-coated plates. Adjust density if necessary, so that the cells are approximately 90 100% confluent on Day 1.
- 9. Incubate at 37°C for 24 hours.
- 10. Continue to section 2 (Differentiating Monolayer Cultures to Definitive Endoderm).

#### TeSR<sup>™</sup>-E8<sup>™</sup> Cultures

This protocol is specific to human ES and iPS cells cultured in TeSR™-E8™ medium.

1. Follow a standard passaging protocol to passage TeSR<sup>™</sup>-E8<sup>™</sup> cultures into one well of a 6-well plate, and perform daily medium changes for 4 days.

NOTE: Refer to the Technical Manual: Maintenance of Human Pluripotent Stem Cells in TeSR<sup>™</sup>-E8<sup>™</sup> (Document #29267) for recommended passaging protocols using TeSR<sup>™</sup>-E8<sup>™</sup>.

 Four days after passaging TeSR<sup>™</sup>-E8<sup>™</sup> cultures, prepare complete TeSR<sup>™</sup>-E8<sup>™</sup> Pre-Differentiation Medium by diluting cold (2 - 8°C) STEMdiff<sup>™</sup> Definitive Endoderm TeSR<sup>™</sup>-E8<sup>™</sup> Supplement 1 in 20 in cold (2 - 8°C) TeSR<sup>™</sup>-E8<sup>™</sup> medium (e.g. add 1 mL of Supplement to 19 mL of TeSR<sup>™</sup>-E8<sup>™</sup>). Prepare sufficient complete TeSR<sup>™</sup>-E8<sup>™</sup> Pre-Differentiation Medium to be used until step 6 (i.e. at least 4 mL per well).

NOTE: Complete Pre-Differentiation Medium can be stored at 2 - 8°C for up to 2 weeks.

- 3. Warm (15 25°C) only the volume of complete TeSR<sup>™</sup>-E8<sup>™</sup> Pre-Differentiation Medium required on this day (i.e. 2 mL per well). Store remaining medium at 2 8°C.
- 4. Aspirate medium from the culture well and add 2 mL of complete TeSR<sup>™</sup>-E8<sup>™</sup> Pre-Differentiation Medium.
- 5. Incubate at 37°C and perform daily medium changes (steps 3 and 4) until cultures are approximately 70% confluent, and are ready to be passaged.

NOTE: For optimal differentiation performance, cells must be exposed to complete TeSR<sup>™</sup>-E8<sup>™</sup> Pre-Differentiation Medium for at least 24 hours before the next passaging step.

- 6. Passage cells (Day 0):
  - i. Warm (15 25°C) sufficient volumes of complete TeSR<sup>™</sup>-E8<sup>™</sup> Pre-Differentiation Medium, DMEM/F-12, and Gentle Cell Dissociation Reagent for passaging. Prepare Single-Cell Passaging Medium by adding Y-27632 to TeSR<sup>™</sup>-E8<sup>™</sup> Pre-Differentiation Medium to reach a final concentration of 10 μM.
  - ii. Wash the well to be passaged with 1 mL of D-PBS (Without Ca++ and Mg++).
  - iii. Aspirate wash medium and add 1 mL of Gentle Cell Dissociation Reagent.
  - iv. Incubate at 37°C for 8 10 minutes.
  - v. Dislodge cells by pipetting up and down 1 3 times using a pipette with a P1000 tip. Ensure all remaining cell aggregates are broken up into single cells.
  - vi. Immediately transfer cells to a tube containing an equal volume of DMEM/F-12. Wash the well once with 1 mL of DMEM/F-12 to collect any remaining cells and transfer to the tube. Centrifuge the tube at 300 x g for 5 minutes.
  - vii. Resuspend cells in 1 mL of Single-Cell Passaging Medium and count the number of live cells using a hemocytometer.

### STEMdiff<sup>™</sup> Definitive Endoderm Kit



- viii. Plate cells at a density of 2.1 x 10^5 cells/cm<sup>2</sup> (i.e. 2 x 10^6 cells/well) onto pre-coated plates. Adjust density if necessary, so that the cells are approximately 90 - 100% confluent on Day 1.
- ix. Incubate at 37°C for 24 hours.
- x. Continue to section 2 (Differentiating Monolayer Cultures to Definitive Endoderm).

### 2. Differentiating Monolayer Cultures to Definitive Endoderm

- 1. On Day 1, warm (37°C) sufficient volumes of DMEM/F-12 and STEMdiff™ Definitive Endoderm Basal Medium for Day 1 use.
- Prepare Medium 1 by diluting both STEMdiff<sup>™</sup> Definitive Endoderm Supplement A and STEMdiff<sup>™</sup> Definitive Endoderm Supplement B 1 in 100 in STEMdiff<sup>™</sup> Definitive Endoderm Basal Medium (e.g. add 10 µL of Supplement A and 10 µL of Supplement B to 980 µL of Basal Medium).

NOTE: Supplements should be thawed on ice and kept cold until added to STEMdiff™ Definitive Endoderm Basal Medium.

- 3. Aspirate medium and wash with 1 mL DMEM/F-12.
- 4. Aspirate wash medium and replace with 2 mL of Medium 1.
- 5. Incubate at 37°C for 24 hours.
- 6. On Day 2, prepare Medium 2 by diluting STEMdiff<sup>™</sup> Definitive Endoderm Supplement B 1 in 100 in STEMdiff<sup>™</sup> Definitive Endoderm Basal Medium (e.g. add 10 µL of Supplement B to 990 µL of Basal Medium). Prepare sufficient Medium 2 to be used on Days 2, 3 and 4 (i.e. 6 mL per well).

NOTE: STEMdiff<sup>™</sup> Definitive Endoderm Supplement B should be thawed on ice and added to cold (2 - 8°C) STEMdiff<sup>™</sup> Definitive Endoderm Basal Medium.

- 7. Warm (37°C) only the volume of Medium 2 required for Day 2 use (i.e. 2 mL per well). Store remaining Medium 2 at 2 8°C.
- Aspirate medium from the well and add 2 mL of Medium 2.
   NOTE: A wash step with DMEM/F-12 is not required at this step or during subsequent medium changes.
- 9. Incubate at 37°C for 24 hours.
- 10. On Day 3, warm (37°C) only the volume of Medium 2 required for Day 3 use (i.e. 2 mL per well). Store remaining Medium 2 at 2 8°C.
- 11. Aspirate medium from the well and add 2 mL of Medium 2.
- 12. Incubate at 37°C for 24 hours.
- 13. On Day 4, warm (37°C) only the volume of Medium 2 required for the Day 4 medium change (i.e. 2 mL per well).
- 14. Aspirate medium from the well and add 2 mL of Medium 2.
- 15. Incubate at 37°C for 24 hours.
- 16. On Day 5, cells are ready to be assayed for the formation of definitive endoderm or carried forward into more specialized lineage differentiation protocols.

NOTE: Expression of definitive endoderm markers may peak by Day 4 in some cell lines.

### Assessment of Definitive Endoderm Cells

Purity of definitive endoderm cells can be measured by flow cytometry after labeling with fluorochrome-conjugated anti-CXCR4 (e.g. Anti-Human CD184 [CXCR4] Antibody, Clone 12G5, Catalog #60089) and anti-c-Kit (e.g. Anti-Human CD117 [c-Kit] Antibody, Clone 104D2, Catalog #60087) or anti-SOX17 antibodies. Results may vary depending on cell line used.

# **Related Products**

For related products, including specialized media, matrices, antibodies, cytokines, and small molecules, visit ww.stemcell.com/DEworkflow or contact us at techsupport@stemcell.com.

STEMCELL TECHNOLOGIES INC.'S QUALITY MANAGEMENT SYSTEM IS CERTIFIED TO ISO 13485. PRODUCTS ARE FOR RESEARCH USE ONLY AND NOT INTENDED FOR HUMAN OR ANIMAL DIAGNOSTIC OR THERAPEUTIC USES UNLESS OTHERWISE STATED.

Copyright © 2016 by STEMCELL Technologies Inc. All rights reserved including graphics and images. STEMCELL Technologies & Design, STEMCELL Shield Design, Scientists Helping Scientists, and STEMdiff are trademarks of STEMCELL Technologies Canada Inc. mTeSR, TeSR and E8 are trademarks of WARF. Corning and Matrigel are registered trademarks of Corning Incorporated. Vitronectin XF is developed and manufactured by Primorigen Biosciences Inc., and Vitronectin XF is a trademark of Primorigen Biosciences, Inc. All other trademarks are the property of their respective holders. While STEMCELL and its suppliers is correct, it makes no warranties or representations as to the accuracy or completeness of such information.