

3% Acetic Acid with Methylene Blue

Reagent for manual counting of nucleated mammalian cells

Catalog # 07060

100 mL



Scientists Helping Scientists™ | WWW.STEMCELL.COM

TOLL FREE PHONE 1 800 667 0322 • PHONE +1 604 877 0713

INFO@STEMCELL.COM • TECHSUPPORT@STEMCELL.COM

FOR GLOBAL CONTACT DETAILS VISIT OUR WEBSITE

Product Description

3% Acetic Acid with Methylene Blue is recommended for manual counting of nucleated mammalian cells. It will lyse red and white blood cell membranes. The remaining white blood cell nuclei will stain lightly with methylene blue.

Properties

Storage: Store at 15 - 25°C.

Shelf Life: Not applicable.

Contains:

- Acetic acid (3%)
- Methylene blue
- Distilled water

Please refer to the Safety Data Sheet (SDS) for hazard information.

Handling / Directions For Use

1. Prepare an appropriate dilution of cells in 3% Acetic Acid with Methylene Blue.
2. Prepare a hemocytometer by first cleaning the chamber surface with alcohol. Wipe dry.
3. Position the coverslip over the chambers. Carefully transfer sufficient volume of the 3% Acetic Acid with Methylene Blue/cell solution to each chamber using a capillary tube or pipetman. Do not over- or underfill.
4. Begin by counting the cells in one chamber. Count all cells in each 1 mm square of each chamber. If cells are on the border outlining each square, count only the cells on the top and left border of the square.

NOTE: Each square has a total volume of 0.1 mm³ (or 10⁻⁴ cm³, which is approximately equivalent to 10⁻⁴ mL).

5. Determine the cell count (cells per mL) as follows:

$$\text{average cell count per square} \times \text{dilution factor} \times 10^4 = \text{cell count per mL}$$

STEMCELL TECHNOLOGIES INC.'S QUALITY MANAGEMENT SYSTEM IS CERTIFIED TO ISO 13485. PRODUCTS ARE FOR RESEARCH USE ONLY AND NOT INTENDED FOR HUMAN OR ANIMAL DIAGNOSTIC OR THERAPEUTIC USES UNLESS OTHERWISE STATED.

Copyright © 2018 by STEMCELL Technologies Inc. All rights reserved including graphics and images. STEMCELL Technologies & Design, STEMCELL Shield Design, and Scientists Helping Scientists are trademarks of STEMCELL Technologies Canada Inc. All other trademarks are the property of their respective holders. While STEMCELL has made all reasonable efforts to ensure that the information provided by STEMCELL and its suppliers is correct, it makes no warranties or representations as to the accuracy or completeness of such information.