

DNase I Solution (1 mg/mL)

Cell dissociation reagent

Catalog # 07900 1 mL



Scientists Helping Scientists™ | WWW.STEMCELL.COM

TOLL FREE PHONE 1 800 667 0322 • PHONE +1 604 877 0713

INFO@STEMCELL.COM • TECHSUPPORT@STEMCELL.COM

FOR GLOBAL CONTACT DETAILS VISIT OUR WEBSITE

Product Description

Deoxyribonuclease (DNase) I Solution (1 mg/mL) is useful for reducing or preventing aggregation of concentrated and/or cryopreserved cell suspensions following thawing. The activity range is ≥ 2000 Kunitz units per mg of protein. One Kunitz unit refers to the amount of DNase I required to act on 1 mg/mL of DNA (pH 5.0, 25°C) to produce an increase in absorbance of 0.001 per minute at a wavelength of 260 nm.

Properties

- Storage:** Store at -20°C.
- Shelf Life:** Stable until expiry date (EXP) on label.
- Contains:**
- Bovine pancreatic deoxyribonuclease I (1 mg/mL)
 - Phosphate-buffered saline

Handling / Directions For Use

Storage at 2 - 8°C is not recommended. Use the protocol below when cell samples appear to be aggregating; this can occur when cells have been frozen and thawed or after enzymatic tissue digestion. For optimal performance, use cells diluted in a medium that does not contain a chelator of calcium and magnesium ions (e.g. EDTA).

1. Thaw DNase I Solution at room temperature (15 - 25°C) or overnight at 2 - 8°C.
 2. Centrifuge cells and carefully remove the supernatant.
 3. Resuspend cells in 0.1 mg/mL of DNase I Solution.
 4. Incubate at room temperature (15 - 25°C) for 15 minutes.
- NOTE: For optimal cell separation results, filter aggregated suspensions through a 37 μ m Reversible Strainer (Catalog #27215/27250) and then resuspend at the appropriate cell concentration in desired medium.
- NOTE: For downstream applications that are sensitive to DNase, wash cells with appropriate DNase-free medium.

STEMCELL TECHNOLOGIES INC.'S QUALITY MANAGEMENT SYSTEM IS CERTIFIED TO ISO 13485. PRODUCTS ARE FOR RESEARCH USE ONLY AND NOT INTENDED FOR HUMAN OR ANIMAL DIAGNOSTIC OR THERAPEUTIC USES UNLESS OTHERWISE STATED.

Copyright © 2018 by STEMCELL Technologies Inc. All rights reserved including graphics and images. STEMCELL Technologies & Design, STEMCELL Shield Design, and Scientists Helping Scientists are trademarks of STEMCELL Technologies Canada Inc. All other trademarks are the property of their respective holders. While STEMCELL has made all reasonable efforts to ensure that the information provided by STEMCELL and its suppliers is correct, it makes no warranties or representations as to the accuracy or completeness of such information.