

0.1% Gelatin in Water

Coating for cultureware

Catalog # 07903 500 mL



Scientists Helping Scientists™ | WWW.STEMCELL.COM

TOLL FREE PHONE 1 800 667 0322 • PHONE +1 604 877 0713

INFO@STEMCELL.COM • TECHSUPPORT@STEMCELL.COM

FOR GLOBAL CONTACT DETAILS VISIT OUR WEBSITE

Product Description

This product is recommended for the coating of culture dishes or flasks to grow mouse embryonic stem (ES) cells, mouse induced pluripotent stem (iPS) cells, or primary embryonic fibroblasts.

Properties

Storage: Store at 15 - 25°C.

Shelf Life: Stable until expiry date (EXP) on label.

Handling / Directions For Use

NOTE: The gelatin solution may be cloudy upon receipt due to exposure to low temperatures during shipping. If this occurs, place bottle in a 37°C water bath for 1 hour (maximum of 2 hours). Ensure the temperature of the water bath is $\geq 35^{\circ}\text{C}$ and $< 40^{\circ}\text{C}$. The gelatin will dissolve, resulting in a clear solution. Product performance will not be affected.

TO COAT CULTURE VESSELS:

1. Aseptically dispense enough gelatin into the dish or flask to completely cover the bottom of the vessel. Recommended volumes are:
 - 3 mL per 25 cm² flask or 60 mm tissue culture dish
 - 7 - 8 mL per 100 mm tissue culture dish
2. Leave gelatin solution in contact with the cultureware for at least 20 minutes at room temperature (15 - 25°C).
3. Aspirate off the excess gelatin solution and leave cultureware sitting open in the biosafety cabinet until dry.

NOTE: Gelatin-coated dishes or flasks can be stored for up to 2 weeks at room temperature (15 - 25°C) or at 2 - 8°C until required.

For further information, refer to the Technical Manual: In Vitro Hematopoietic Differentiation of Mouse ES & iPS Cells Using ES-Cult™ (Document #28415), available on our website at www.stemcell.com or contact us to request a copy.

STEMCELL TECHNOLOGIES INC.'S QUALITY MANAGEMENT SYSTEM IS CERTIFIED TO ISO 13485. PRODUCTS ARE FOR RESEARCH USE ONLY AND NOT INTENDED FOR HUMAN OR ANIMAL DIAGNOSTIC OR THERAPEUTIC USES UNLESS OTHERWISE STATED.

Copyright © 2016 by STEMCELL Technologies Inc. All rights reserved including graphics and images. STEMCELL Technologies & Design, STEMCELL Shield Design, Scientists Helping Scientists, and ES-Cult are trademarks of STEMCELL Technologies Canada Inc. All other trademarks are the property of their respective holders. While STEMCELL has made all reasonable efforts to ensure that the information provided by STEMCELL and its suppliers is correct, it makes no warranties or representations as to the accuracy or completeness of such information.