

Hydrocortisone

Cell culture supplement

Catalog # 07904 100 mg



Scientists Helping Scientists™ | WWW.STEMCELL.COM

TOLL FREE PHONE 1 800 667 0322 • PHONE +1 604 877 0713

INFO@STEMCELL.COM • TECHSUPPORT@STEMCELL.COM

FOR GLOBAL CONTACT DETAILS VISIT OUR WEBSITE

Product Description

Hydrocortisone is required for supplementation of the following media to give a final concentration of 10^{-6} M:

- MyeloCult™ H5100 (Human, Catalog #05100) or MyeloCult™ M5300 (Mouse, Catalog #05300), for long-term cultures and the long-term culture-initiating cell assay
- EpiCult™-B Human Medium Kit (Catalog #05601), for culture of human mammary luminal and myoepithelial cells
- MammoCult™ Human Medium Kit (Catalog #05620), for culture of mammospheres and tumorspheres


Properties

Storage: Store at -20°C .

Shelf Life: Store powder with a desiccant at -20°C .

Handling / Directions For Use

1. To prepare a 10^{-3} M solution of hydrocortisone, accurately weigh out 2.42 mg and dilute to 5 mL in an appropriate medium (e.g. Alpha MEM with Nucleosides [Catalog #36450]). Alternatively, weigh out a small amount of powder and dissolve at a ratio of 1 mg of powder to 2 mL of medium.
2. Filter sterilize using a using a 0.2 - 0.22 μm low protein binding polyethersulfone (PES) filter unit (e.g. Fisher 09-741-04 [0.2 μm , 250 mL]; Fisher SCGP00525 [0.22 μm , 50 mL]).
3. Perform a further 1 in 10 dilution of this stock in medium to obtain a working stock solution of 10^{-4} M.
NOTE: If not used immediately, store the working stock solution at 2 - 8°C for up to 1 week.
4. Add 1 mL of 10^{-4} M hydrocortisone to 100 mL of medium to achieve a final concentration of 10^{-6} M hydrocortisone.



Hydrocortisone

STEMCELL TECHNOLOGIES INC.'S QUALITY MANAGEMENT SYSTEM IS CERTIFIED TO ISO 13485. PRODUCTS ARE FOR RESEARCH USE ONLY AND NOT INTENDED FOR HUMAN OR ANIMAL DIAGNOSTIC OR THERAPEUTIC USES UNLESS OTHERWISE STATED.

Copyright © 2018 by STEMCELL Technologies Inc. All rights reserved including graphics and images. STEMCELL Technologies & Design, STEMCELL Shield Design, Scientists Helping Scientists, MyeloCult, EpiCult, and MammoCult are trademarks of STEMCELL Technologies Canada Inc. All other trademarks are the property of their respective holders. While STEMCELL has made all reasonable efforts to ensure that the information provided by STEMCELL and its suppliers is correct, it makes no warranties or representations as to the accuracy or completeness of such information.