

# ACCUTASE™

## Cell detachment solution

Catalog # 07920                      100 mL



Scientists Helping Scientists™ | [WWW.STEMCELL.COM](http://WWW.STEMCELL.COM)

TOLL FREE PHONE 1 800 667 0322 • PHONE +1 604 877 0713

[INFO@STEMCELL.COM](mailto:INFO@STEMCELL.COM) • [TECHSUPPORT@STEMCELL.COM](mailto:TECHSUPPORT@STEMCELL.COM)

FOR GLOBAL CONTACT DETAILS VISIT OUR WEBSITE

## Product Description

ACCUTASE™ is a cell detachment solution of proteolytic and collagenolytic enzymes. It is useful for the routine detachment of cells from standard tissue culture plasticware and adhesion-coated plasticware. ACCUTASE™ does not contain mammalian- or bacterial-derived products. Each lot of ACCUTASE™ is tested for sterility (by USP membrane filtration method), enzymatic activity (tested with synthetic chromagenic tetrapeptides), and cell detachment from tissue culture plastic.

ACCUTASE™ has been shown to efficiently detach or dissociate fibroblasts, neurospheres, keratinocytes, vascular endothelial cells, hepatocytes, vascular smooth muscle cells, hepatocyte progenitor cells, bone marrow stem cells, macrophages, 293 cells, L929 cells, immortalized mouse testicular germ cells, 3T3, Vero, COS, HeLa, NT2, MG63, M24 and A375 metastatic melanomas, gliomas U251 and D54, and T1080 fibrosarcoma cells.

## Properties

**Storage:** Store at -20°C.

**Shelf Life:** Stable until expiry date (EXP) on label.

**Contains:**

- 1X ACCUTASE™ enzymes in Dulbecco's phosphate-buffered saline (PBS)
- 0.5 mM EDTA•4Na
- 3 mg/L Phenol red

## Handling / Directions For Use

1. Thaw ACCUTASE™ overnight at 2 - 8°C, at room temperature (15 - 25°C), or in a container of cool water. Do not thaw at 37°C. Mix thoroughly.  
NOTE: Once thawed, use immediately or store at 2 - 8°C for up to 2 months. Alternatively, aliquot and store at -20°C. Do not exceed the expiry date on the label.
2. Wash plate, flask, or beads with sterile PBS.
3. Add ACCUTASE™ to culture dish or flask at 10 mL per 75 cm<sup>2</sup> surface area.
4. Incubate at 37°C for 5 - 10 minutes, allowing cells to detach.
5. Count cells and passage as usual; no additional washes or enzyme inhibitors are required.

STEMCELL TECHNOLOGIES INC.'S QUALITY MANAGEMENT SYSTEM IS CERTIFIED TO ISO 13485. PRODUCTS ARE FOR RESEARCH USE ONLY AND NOT INTENDED FOR HUMAN OR ANIMAL DIAGNOSTIC OR THERAPEUTIC USES UNLESS OTHERWISE STATED.

Copyright © 2016 by STEMCELL Technologies Inc. All rights reserved including graphics and images. STEMCELL Technologies & Design, STEMCELL Shield Design, and Scientists Helping Scientists are trademarks of STEMCELL Technologies Canada Inc. All other trademarks are the property of their respective holders. ACCUTASE is a trademark of Innovative Cell Technologies, Inc. While STEMCELL has made all reasonable efforts to ensure that the information provided by STEMCELL and its suppliers is correct, it makes no warranties or representations as to the accuracy or completeness of such information.