## RosetteSep™ DM-L Density Medium

REF 15705

100 mL

### ENGLISH

#### INTENDED USE

RosetteSep™ DM-L Density Medium is a density separation medium designed specifically for use with RosetteSep™ cocktails for the enrichment of specific human lymphocyte populations (CD3+ T cells, CD4+ T cells, CD8+ T cells, and B cells) from whole blood.

#### PRODUCT SPECIFICATIONS

RosetteSep™ DM-L Density Medium

REF #15705

Density: 1.081 g/mL

#### QUALITY CONTROL

RosetteSep™ DM-L Density Medium is manufactured using aseptic technique and tightly controlled processes.

Each lot of RosetteSep™ DM-L Density Medium is sterility tested according to USP methods.

#### STORAGE AND STABILITY

Store at 15 - 25°C. Storage at 2 - 8°C is acceptable, but ensure that the medium equilibrates to 15 - 25°C and invert bottle to mix contents before use. Keep protected from direct light. Product stable at 15 - 25°C until expiry date as indicated on label.

#### WARNINGS AND PRECAUTIONS

1. For professional users only.
2. This product is for in vitro diagnostic use.
3. Do not use bottle if contents have leaked. Unused density medium may be disposed of according to standard laboratory procedures for non-hazardous liquids.
4. This product should be handled by trained personnel observing good laboratory practices. Once this product is added to human cells, treat the suspension as potentially biohazardous. Handling of reagents and disposal of wastes should observe all local, state, or national regulations.
5. This product is a potential irritant to eyes, respiratory system, and skin. This product may also be harmful if ingested. Avoid exposure through skin, eye contact, inhalation, and ingestion.

#### HANDLING AND DIRECTIONS FOR USE

For a list of RosetteSep™ cell enrichment kits this product can be used with, see the Related Products section on this Product Information Sheet. For detailed instructions, refer to the appropriate RosetteSep™ Product Information Sheet.

**NOTE: Do not use RosetteSep™ DM-L Density Medium with RosetteSep™ Human Monocyte Enrichment Cocktail (Catalog #15028) or RosetteSep™ Human Progenitor Cocktails (e.g. Catalog #15276, #15026, #15126, #15027, or #15128). The density of RosetteSep™ DM-L Density Medium is not optimized for the isolation of these cell types.**

#### RELATED PRODUCTS

PRODUCT	CATALOG #
<i>Included in:</i>	
RosetteSep™ HLA Human Lymphoid Cell Enrichment Kit	15271HLA
<i>Recommended for use with:</i>	
RosetteSep™ Human T Cell Enrichment Cocktail	15021/15061
RosetteSep™ Human CD4+ T Cell Enrichment Cocktail	15022/15062
RosetteSep™ Human CD8+ T Cell Enrichment Cocktail	15023/15063
RosetteSep™ Human B Cell Enrichment Cocktail	15024/15064
RosetteSep™ Human Total Lymphocyte Enrichment Cocktail	15223/15263
RosetteSep™ HLA T Cell Enrichment Cocktail	15061HLA/ 15081HLA
RosetteSep™ HLA B Cell Enrichment Cocktail	15064HLA/ 15084HLA
RosetteSep™ HLA Total Lymphocyte Enrichment Cocktail	15263HLA/ 15283HLA


STEMCELL Technologies Inc | 1618 Station Street, Vancouver, BC | V6A 1B6 | Canada | [www.stemcell.com](http://www.stemcell.com)

#### For Technical Assistance

Tel: +1.604.877.0713

European toll-free number: 00800 7836 2355


E-mail: [techsupport@stemcell.com](mailto:techsupport@stemcell.com)

MDSS GmbH

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30175 Hannover • Germany

Document #29684

Version 3.4.0

2017

Page 1 of 10


## TECHNICAL ASSISTANCE

For technical support, contact us by email at [techsupport@stemcell.com](mailto:techsupport@stemcell.com) or call either **+1.604.877.0713** or the European toll-free number 00800 7836 2355. For more information, visit [www.stemcell.com](http://www.stemcell.com).

If you require a printed copy or a translated version of this document in a certain language, contact technical support.

 **MDSS GmbH**

Schiffgraben 41  
30175 Hannover, Germany

 Catalog or reference number	 Batch code	 Use by: YYYY-MM
 Caution, consult accompanying documents	 In Vitro Diagnostic Medical Device	 For storage within temperature limits
 CE Mark	 Manufacturer's identification (name & address)	 Authorized EC representative in the European Community

STEMCELL TECHNOLOGIES INC.'S QUALITY MANAGEMENT SYSTEM IS CERTIFIED TO ISO 13485. PRODUCTS ARE FOR PROFESSIONAL IN VITRO DIAGNOSTIC USE.

Copyright © 2017 by STEMCELL Technologies Inc. All rights reserved including graphics and images. STEMCELL Technologies & Design, STEMCELL Shield Design, Scientists Helping Scientists, and RosetteSep are trademarks of STEMCELL Technologies Canada Inc. All other trademarks are the property of their respective holders. While STEMCELL has made all reasonable efforts to ensure that the information provided by STEMCELL and its suppliers is correct, it makes no warranties or representations as to the accuracy or completeness of such information.

 **STEMCELL Technologies Inc** | 1618 Station Street, Vancouver, BC | V6A 1B6 | Canada | [www.stemcell.com](http://www.stemcell.com)

**For Technical Assistance**  
Tel: +1.604.877.0713  
European toll-free number: 00800 7836 2355  
E-mail: [techsupport@stemcell.com](mailto:techsupport@stemcell.com)

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
STEMCELL Technologies Inc | 1618 Station Street, Vancouver, BC | V6A 1B6 | Canada | www.stemcell.com

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
 **STEMCELL Technologies Inc | 1618 Station Street, Vancouver, BC | V6A 1B6 | Canada | [www.stemcell.com](http://www.stemcell.com)**

**For Technical Assistance**  
Tel: +1.604.877.0713  
European toll-free number: 00800 7836 2355  
E-mail: [techsupport@stemcell.com](mailto:techsupport@stemcell.com)

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