

# StemSpan™-ACF



Scientists Helping Scientists™ | [WWW.STEMCELL.COM](http://WWW.STEMCELL.COM)

TOLL FREE PHONE 1 800 667 0322 • PHONE +1 604 877 0713

[INFO@STEMCELL.COM](mailto:INFO@STEMCELL.COM) • [TECHSUPPORT@STEMCELL.COM](mailto:TECHSUPPORT@STEMCELL.COM)

FOR GLOBAL CONTACT DETAILS VISIT OUR WEBSITE

## Animal component-free medium for culture and expansion of human hematopoietic cells

Catalog # 09805	100 mL
09855	500 mL

## Product Description

StemSpan™-ACF is an animal component-free medium for the in vitro culture and expansion of human hematopoietic cells. This medium contains recombinant proteins and synthetic components. It does not contain serum or other human- or animal-derived components. Recombinant hematopoietic growth factors, required for the optimal growth and expansion of hematopoietic cells, have not been added to StemSpan™-ACF. This allows users the flexibility to add growth factors and other supplements according to their requirements.

## Properties

**Storage:** Store at 2 - 8°C.  
**Shelf Life:** Stable for 12 months from date of manufacture (MFG) on label.  
**Contains:** This product contains only recombinant proteins and synthetic components.

## Handling / Directions For Use

1. Add desired cytokines, growth factors and other components to StemSpan™-ACF and mix well.  
NOTE: Added components and cells in sterile cell culture medium (e.g. Iscove's MDM or DMEM) should not exceed ~10% of total volume.
2. Add cells, mix well, and set up cultures as desired.

### ASSESSMENT OF HEMATOPOIETIC CELLS

Assessment of CD34+ cells before and after culture may be performed by flow cytometry using the following fluorochrome-conjugated antibody clones:

- Anti-Human CD34 Antibody, Clone 581 (Catalog #60013) or Clone 563 (Catalog #60119) or Clone 8G12 (Catalog #60121)
- Anti-Human CD45 Antibody, Clone HI30 (Catalog #60018) or Clone 2D1 (Catalog #60123)

## Notes and Tips

Selection of an optimal growth factor combination is dependent upon the source and type of cells and the experimental objectives of the researcher. StemSpan™ expansion supplements, described below, are suitable for use with StemSpan™-ACF.

- StemSpan™ CD34+ Expansion Supplement (10X) (Catalog #02691)
  - Culture and expansion of large numbers of human CD34+ progenitor cells
  - Contains: rh SCF, rh TPO, rh IL-3, rh IL-6, rh Flt3 ligand, other additives
- StemSpan™ CC100 (Catalog #02690)
  - Culture and expansion of human hematopoietic cells
  - Contains: rh Flt3 ligand, rh SCF, rh IL-3, rh IL-6
- StemSpan™ CC110 (Catalog #02697)
  - Culture and expansion of human hematopoietic cells
  - Contains: rh Flt3 ligand, rh SCF, rh TPO
- StemSpan™ Erythroid Expansion Supplement (100X) (Catalog #02692)
  - Expansion and lineage-specific differentiation of human CD34+ cells into erythroid progenitor cells
  - Contains: rh SCF, rh IL-3, rh EPO

- StemSpan™ Megakaryocyte Expansion Supplement (100X) (Catalog #02696)
  - Expansion and lineage-specific differentiation of human CD34+ cells into megakaryocyte progenitor cells
  - Contains: rh SCF, rh TPO, rh IL-6, rh IL-9
- StemSpan™ Myeloid Expansion Supplement (100X) (Catalog #02693)
  - Expansion and lineage-specific differentiation of human CD34+ cells into myeloid progenitor cells
  - Contains: rh SCF, rh TPO, rh G-CSF, rh GM-CSF

SCF = Stem cell factor; EPO = Erythropoietin; TPO = Thrombopoietin; rh = recombinant human; IL = interleukin; Flt = fms-like tyrosine kinase

#### RELATED PRODUCTS

For related products, including specialized culture and storage media, supplements, antibodies, cytokines, and small molecules, visit [www.stemcell.com/HSPCworkflow](http://www.stemcell.com/HSPCworkflow) or contact us at [techsupport@stemcell.com](mailto:techsupport@stemcell.com). For available fresh and cryopreserved peripheral blood, cord blood, and bone marrow products in your region, visit [www.stemcell.com/primarycells](http://www.stemcell.com/primarycells).

STEMCELL TECHNOLOGIES INC.'S QUALITY MANAGEMENT SYSTEM IS CERTIFIED TO ISO 13485. PRODUCTS ARE FOR RESEARCH USE ONLY AND NOT INTENDED FOR HUMAN OR ANIMAL DIAGNOSTIC OR THERAPEUTIC USES UNLESS OTHERWISE STATED.

Copyright © 2017 by STEMCELL Technologies Inc. All rights reserved including graphics and images. STEMCELL Technologies & Design, STEMCELL Shield Design, Scientists Helping Scientists, and StemSpan are trademarks of STEMCELL Technologies Canada Inc. All other trademarks are the property of their respective holders. While STEMCELL has made all reasonable efforts to ensure that the information provided by STEMCELL and its suppliers is correct, it makes no warranties or representations as to the accuracy or completeness of such information. While STEMCELL has made all reasonable efforts to ensure that the information provided by STEMCELL and its suppliers is correct, it makes no warranties or representations as to the accuracy or completeness of such information.