

CellAdhere™ Type I Collagen, Human, Lyophilized

Type I collagen, human, for tissue engineering research and cell culture

Catalog # 07005

15 mg



Scientists Helping Scientists™ | WWW.STEMCELL.COM

TOLL FREE PHONE 1 800 667 0322 • PHONE +1 604 877 0713

INFO@STEMCELL.COM • TECHSUPPORT@STEMCELL.COM

FOR GLOBAL CONTACT DETAILS VISIT OUR WEBSITE

Product Description

Type I Collagen is the most abundant type of collagen in the human body and is a major structural protein prevalent in skin, bone, tendon, and other fibrous connective tissues (Bornstein et al.). All collagens are characterized by a secondary protein structure consisting of three polypeptide chains wrapped around one another in a triple helical arrangement (Tanzer). The different types of collagen vary in the primary amino acid sequence of their polypeptide chains, with type I collagen comprising two identical $\alpha 1(I)$ chains and one $\alpha 2(I)$ chain that spontaneously form the triple helix in solution (Heino).

CellAdhere™ Type I Collagen, Human, Lyophilized is prepared from the extracellular matrix secreted by normal human fibroblasts and is xeno-free. When reconstituted with 5 mL of 0.01 N HCl, the resulting concentration is 3 mg/mL. This product is purified using a highly controlled manufacturing process with validated steps to ensure inactivation of possible viral contaminants, and it is sterility tested. The endotoxin level is 5.0 EU/mL or less, as assayed by Limulus amoebocyte lysate (LAL) assay.

This product is ideal for coating surfaces with a thin layer of protein to support cell attachment, but is not suitable for 3D cell culture applications. The optimal protein concentration may vary depending on the cell type being used, and therefore must be titrated for best results.

Properties

Storage: Store at -20°C.

Shelf Life: Stable for 6 months from date of receipt.

This product contains material derived from human fibroblasts. The cell line has been tested and found negative for HIV-1 and -2, hepatitis B, and hepatitis C. However, this product should be considered potentially infectious and treated in accordance with universal handling precautions.

Handling / Directions For Use

PREPARING COLLAGEN-COATED TISSUE CULTUREWARE

1. Reconstitute CellAdhere™ Type I Collagen, Human, Lyophilized (15 mg) with 5 mL of sterile 0.01 N HCl solution to obtain a 3 mg/mL solution.
2. Mix gently until everything has been solubilized.
NOTE: Reconstituted solution may be slightly opaque. If not used immediately, store at 2 - 8°C for up to 3 months.
3. Dilute reconstituted Type I Collagen solution with sterile 0.01 N HCl solution to obtain desired concentration.
NOTE: Different dilutions will need to be tested to determine the optimal concentration required for each culture system. Typical coating concentrations range between 10 - 100 µg/mL.
4. Add desired amount of diluted Type I Collagen solution on the surface of the cultureware to be coated.
For example, use 1 mL to coat a 35 mm Culture Dish (Catalog #27100).
5. Cover coated cultureware to protect from contamination. Incubate at room temperature (15 - 25°C) for 1 - 2 hours.
7. Aspirate excess solution. Avoid scratching the coated surface.
8. Rinse coated cultureware with sterile D-PBS (Without Ca⁺⁺ and Mg⁺⁺) (Catalog #37350) or culture medium.
9. Use coated cultureware immediately. Alternatively, keep sterile and store at 2 - 8°C damp or air dried.

References

- Bornstein P & Sage H. (1980) Structurally distinct collagen types. *Annu Rev Biochem* 49: 957–1003.
Heino J. (2007) The collagen family members as cell adhesion proteins. *Bioessays* 29(10): 1001–10.
Tanzer ML. (1973) Cross-linking of collagen. *Science* 180(4086): 561–6.

STEMCELL TECHNOLOGIES INC.'S QUALITY MANAGEMENT SYSTEM IS CERTIFIED TO ISO 13485. PRODUCTS ARE FOR RESEARCH USE ONLY AND NOT INTENDED FOR HUMAN OR ANIMAL DIAGNOSTIC OR THERAPEUTIC USES UNLESS OTHERWISE STATED.

Copyright © 2018 by STEMCELL Technologies Inc. All rights reserved including graphics and images. STEMCELL Technologies & Design, STEMCELL Shield Design, Scientists Helping Scientists, and CellAdhere are trademarks of STEMCELL Technologies Canada Inc. All other trademarks are the property of their respective holders. While STEMCELL has made all reasonable efforts to ensure that the information provided by STEMCELL and its suppliers is correct, it makes no warranties or representations as to the accuracy or completeness of such information.