

3 cc Syringes



Scientists Helping Scientists™ | WWW.STEMCELL.COM

TOLL FREE PHONE 1 800 667 0322 • PHONE +1 604 877 0713
INFO@STEMCELL.COM • TECHSUPPORT@STEMCELL.COM
FOR GLOBAL CONTACT DETAILS VISIT OUR WEBSITE

Catalog # 28230	30 Syringes/Bag
28240	100 Syringes/Bag

Product Description

Sterile 3 cc Syringes are ideal for accurate dispensing of methylcellulose-based medium (e.g. MethoCult™), when used with Blunt-End Needles, 16 Gauge (Catalog #28110) purchased from STEMCELL Technologies to ensure compatibility and accuracy. Manufactured from non-pyrogenic materials.

Properties

Storage: Store at 15 - 25°C.
Shelf Life: Stable until expiry date (EXP) on label.

Handling / Directions For Use

Use a 3 cc Syringe and Blunt-End Needle, 16 Gauge (Catalog #28110) to dispense MethoCult™ methylcellulose-based medium. This is for accurate dispensing of viscous methylcellulose medium and to prevent needle-stick injuries.

Use 3 cc Syringes with Blunt-End Needles, 16 Gauge purchased from STEMCELL Technologies to ensure that they are compatible. STEMCELL Technologies cannot guarantee that 3 cc Syringes will be compatible with blunt-end needles from other suppliers.

3 cc Syringes and Blunt-End Needles, 16 Gauge are single-use only. Dispose according to standard procedures for biohazardous materials.

For complete instructions, refer to one of the following Technical Manuals, available at www.stemcell.com or contact us to request a copy:

- Human Colony-Forming Unit Assays Using MethoCult™ (Document #28404)
- Mouse Colony-Forming Unit Assays Using MethoCult™ (Document #28405)

STEMCELL TECHNOLOGIES INC.'S QUALITY MANAGEMENT SYSTEM IS CERTIFIED TO ISO 13485. PRODUCTS ARE FOR RESEARCH USE ONLY AND NOT INTENDED FOR HUMAN OR ANIMAL DIAGNOSTIC OR THERAPEUTIC USES UNLESS OTHERWISE STATED.

Copyright © 2019 by STEMCELL Technologies Inc. All rights reserved including graphics and images. STEMCELL Technologies & Design, STEMCELL Shield Design, Scientists Helping Scientists, and MethoCult are trademarks of STEMCELL Technologies Canada Inc. All other trademarks are the property of their respective holders. While STEMCELL has made all reasonable efforts to ensure that the information provided by STEMCELL and its suppliers is correct, it makes no warranties or representations as to the accuracy or completeness of such information.