

ImmunoCult™ Mouse Th2 Differentiation Supplement

Supplement for the differentiation of mouse naïve CD4+ T cells into Th2 cells

Catalog # 10955 1 mL



Scientists Helping Scientists™ | WWW.STEMCELL.COM

TOLL FREE PHONE 1 800 667 0322 • PHONE +1 604 877 0713

INFO@STEMCELL.COM • TECHSUPPORT@STEMCELL.COM

FOR GLOBAL CONTACT DETAILS VISIT OUR WEBSITE

Product Description

ImmunoCult™ Mouse Th2 Differentiation Supplement contains a combination of recombinant mouse cytokines and a blocking monoclonal antibody optimized to promote the differentiation of mouse naïve CD4+ T cells into Th2 cells.

This supplement is intended for use with RPMI Medium (Catalog #36750) containing fetal bovine serum and other additives, as well as anti-mouse CD3 and anti-mouse CD28 monoclonal antibodies as activating agents.

- Optimized for the induction of Th2 cells from naïve CD4+ T cells isolated from the spleen of a BALB/c mouse
- Supplied as a 100X concentrate; after thawing and mixing, the tube contents can be added directly to medium.

Properties

- Storage:** Store at -20°C.
- Shelf Life:** Stable for 12 months from date of manufacture (MFG) on label.
- Contains:**
- Recombinant mouse interleukin 2 (IL-2)
 - Recombinant mouse interleukin 4 (IL-4)
 - Rat anti-mouse interferon-gamma (anti-IFN- γ)

Handling / Directions For Use

Please read the entire protocol before proceeding. Use sterile techniques when performing the following protocols.

A. COATING CULTUREWARE WITH ANTI-MOUSE CD3 ANTIBODY

Coat a flat-bottom tissue culture-treated plate (e.g. Catalog #38015) with anti-mouse CD3 antibody (e.g. Anti-Mouse CD3e Antibody, Clone 145-2C11, Catalog #60015) at a density of 312.5 ng/cm². Cover plate with plastic wrap and store at 2 - 8°C overnight.

B. PREPARATION OF Th2 DIFFERENTIATION MEDIUM

1. Thaw ImmunoCult™ Mouse Th2 Differentiation Supplement on ice until just thawed. Mix thoroughly.
NOTE: If necessary, centrifuge vial for 30 seconds to recover liquid from cap.
NOTE: If not used immediately, store at 2 - 8°C for up to 1 month. Do not re-freeze.
2. Add the following components to RPMI 1640 Medium (Catalog #36750) and mix thoroughly:
 - 5 - 10% fetal bovine serum
 - 2 mM L-Glutamine (Catalog #07100)
 - 10 mM HEPES Buffer Solution (Catalog #07200)
 - 1 mM Sodium Pyruvate (Catalog #07000)
 - 100 μ M MEM Non-Essential Amino Acid Solution (Catalog #07600)
 - 50 μ M β -mercaptoethanol
3. Add ImmunoCult™ Mouse Th2 Differentiation Supplement at a 1 in 100 dilution. Mix thoroughly.
NOTE: If not used immediately, store at 2 - 8°C for up to 1 month.

C. DIFFERENTIATION TO Th2 CELLS

1. DAY 0: Remove and discard antibody from coated wells (prepared in section A). Wash coated wells twice with sterile phosphate-buffered saline.
2. Add anti-mouse CD28 antibody (recommended clone 37.51) to Th2 Differentiation Medium (prepared in section B) to a final concentration of 0.5 μ g/mL. Mix thoroughly.

3. Isolate naïve CD4+ T cells using EasySep™ Mouse Naïve CD4+ T Cell Isolation Kit (Catalog #19765). Dilute cells to 2.5×10^5 cells/mL in Th2 Differentiation Medium + anti-mouse CD28 antibody (prepared in step 2).
NOTE: ImmunoCult™ Mouse Th2 Differentiation Supplement is optimized for the polarization of naïve CD4+ T cells isolated from the spleen of BALB/c mice.
4. Add cell suspension to coated wells at a density of 1.56×10^5 cells/cm². Incubate at 37°C and 5% CO₂ for 2 days.
5. DAY 2: Add cells to fresh tissue culture-treated wells at a split ratio of 1:2 to 1:4. Add fresh Th2 Differentiation Medium (without anti-mouse CD28 antibody). Incubate at 37°C and 5% CO₂ for 2 days.
NOTE: The cells do not need to be re-stimulated with anti-mouse CD3 and anti-mouse CD28 antibodies.
6. DAY 4: Add cells to fresh tissue culture-treated wells at a split ratio of 1:2 to 1:4. Add fresh Th2 Differentiation Medium (without anti-mouse CD28 antibody). Incubate at 37°C and 5% CO₂ for 2 days.
7. DAY 6: Th2 cells are ready to be assayed in the desired application.

Notes and Tips

For assessment of Th2 cells by flow cytometry the following fluorochrome-conjugated antibody clones can be used:

- Anti-mouse CD4 antibody, clone GK1.5
- Anti-mouse GATA3 antibody, clone L50-823
- Anti-mouse IL-4 antibody, clone BVD4-1D11

STEMCELL TECHNOLOGIES INC.'S QUALITY MANAGEMENT SYSTEM IS CERTIFIED TO ISO 13485. PRODUCTS ARE FOR RESEARCH USE ONLY AND NOT INTENDED FOR HUMAN OR ANIMAL DIAGNOSTIC OR THERAPEUTIC USES UNLESS OTHERWISE STATED.

Copyright © 2018 by STEMCELL Technologies Inc. All rights reserved including graphics and images. STEMCELL Technologies & Design, STEMCELL Shield Design, Scientists Helping Scientists, EasySep, and ImmunoCult are trademarks of STEMCELL Technologies Canada Inc. All other trademarks are the property of their respective holders. While STEMCELL has made all reasonable efforts to ensure that the information provided by STEMCELL and its suppliers is correct, it makes no warranties or representations as to the accuracy or completeness of such information.