

Handling and Storage

Upon receipt, immediately transfer the cryovial to liquid nitrogen storage.

Cell Culture Surfaces

For most applications, use tissue culture vessels pre-coated with Poly-D-Lysine.

Other cell culture surfaces can influence adherence, morphology and function; refer to the iCell® Microglia User's Guide for more information and protocols.

Preparing the Maintenance Medium

1. Prepare maintenance medium (see **Table 1**).
2. Filter maintenance medium using a 0.2 µm PES filter unit.
3. Store maintenance medium at 4°C for up to 2 weeks.
4. Single-use aliquots can be stored at -20°C for up to 3 months.
5. Equilibrate maintenance medium to room temperature before use.

Thawing the Cells

1. Transfer 8 ml of maintenance medium to a 15 ml centrifuge tube.
2. Thaw iCell Microglia cryovial in a 37°C water bath for 3 minutes. Clean with 70% ethanol.
3. Transfer the cells to the 15 ml centrifuge tube containing 8 ml of maintenance medium.
4. Rinse the cryovial with 1 ml of maintenance medium and add it to the centrifuge tube.
5. Gently mix by inverting the centrifuge tube or slowly pipetting.
6. Centrifuge the cells at 1000 x g for 10 minutes.
7. Carefully remove the supernatant leaving approximately 200-300 µl above the pellet to avoid disturbing the cell pellet.

Plating the Cells

1. Dilute the cell suspension with maintenance medium to obtain the desired cell plating density using the total viable cells from the Certificate of Analysis. See table below for plating density examples.

Culture Vessel	Surface Area (cm ²)	Plating Volume (ml)	Cell Number
12-well cell culture plate	3.8	1.2	1.8 x 10 ⁵
96-well cell culture plate	0.32	0.1	1.5 x 10 ⁴

2. Dispense the cells into the cell culture vessel.
3. Culture the cells at 37°C, 5% CO₂.

Maintaining the Cells

1. Replace 50% of the medium every 2-3 days.
2. Culture the cells at 37°C, 5% CO₂.

Contacting Technical Support

Email: techsupport@stemcell.com




Figure 1: iCell Microglia cells, 20X

Table 1: Maintenance Medium Formulation for iCell Microglia

Component	Vendor Catalog #	Volume (ml)
DMEM/F-12, HEPES, no phenol red	ThermoFisher #110390	93.3
N-2 Supplement, 100X	ThermoFisher #175020	0.5
B-27 Supplement, 50X	ThermoFisher #175040	1
10% BSA in DPBS ¹	MilliporeSigma #A1470	0.5
1-Thioglycerol (MTG) 11.5M	MilliporeSigma #M6145	0.004
Ascorbic Acid, 20mg/ml ²	FUJIFILM Wako Chemicals #013-19641	0.25
Penicillin-Streptomycin	ThermoFisher #15140	1
GlutaMAX Supplement	ThermoFisher #350500	1
MEM Non-essential Amino Acids, 100X	ThermoFisher #111400	1
Insulin-Transferrin-Selenium, 100X	ThermoFisher #41400045	1
Human Insulin Solution	MilliporeSigma #I9278	0.05
rhM-CSF, 100µg/ml ²	PeproTech #300-25	0.025
rhTGF-β 1, 100 µg/ml ²	R&D Systems #240-B	0.05
rhIL-34, 100µg/ml ²	PeproTech #200-34	0.1
rhFractalkine, 100 µg/ml ²	PeproTech #300-31	0.1
rhCD200, 100 µg/ml ²	Acro Biosystems #OX2-H5228	0.1


¹ DPBS no calcium, no magnesium, ThermoFisher #14190

² Reconstitute according to manufacturer's recommendations

Conditions of Use

The cells are for RESEARCH USE ONLY. See www.fujifilmcdi.com/product-warranty/ for USE RESTRICTIONS applicable to the cells and other terms and conditions related to the cells and their use.

Trademarks

iCell and MyCell are registered trademarks, and Cellular Dynamics and the  logo are trademarks of FUJIFILM Cellular Dynamics, Inc. All other brands, product names, company names, trademarks, and service marks are the properties of their respective owners.

Copyright Notice

© 2018 FUJIFILM Cellular Dynamics, Inc. All rights reserved. This document may not be reproduced, distributed, modified or publicly displayed without the prior express written permission of FUJIFILM Cellular Dynamics, Inc.

Copyright © 2019 by STEMCELL Technologies Inc. All rights reserved including graphics and images. STEMCELL Technologies & Design, STEMCELL Shield Design, and Scientists Helping Scientists are trademarks of STEMCELL Technologies Canada Inc. All other trademarks are the property of their respective holders. While STEMCELL has made all reasonable efforts to ensure that the information provided by STEMCELL and its suppliers is correct, it makes no warranties or representations as to the accuracy or completeness of such information.

Revision History

Document ID: X1022

Version 1.0 December 2018

Scientists Helping Scientists™ | WWW.STEMCELL.COM