

EBV (LMP2) Peptide Pool

Epstein-Barr virus (LMP2) peptide pool for immune cell activation

Catalog #100-0671

1 Unit

~25 µg (15 nmol)/peptide

Product Description

EBV (LMP2) Peptide Pool is a lyophilized mixture of 122 peptides from latent membrane protein 2 (LMP2) of Epstein-Barr virus (EBV; strain B95-8). LMP2 is expressed in EBV latency, as well as in EBV-associated pathologies (Lee et al.; Longnecker) and it blocks reactivation of EBV from latency (Miller et al.). The pool consists of 15-mer peptides with 11-amino-acid overlaps that cover amino acids 1 - 497 on LMP2. One unit of this product (i.e. ~25 µg/peptide) is sufficient for stimulating 2.5×10^8 cells.

Product Information

Amino Acid Sequence:	MGSLEMVPMGAGPPSPGGDPDGYDGGNNSQYPSASGSSGNTPTPPNDEERESNEEPPPPYEDPYW GNGDRHSDYQPLGTQDQSLYLGLQHDGNDGLPPPYSPRDDSSQHIYEEAGRGSMNPVCLPVIVAPY LFWLAAIAASCFTASVSTVVTATGLALSLLLLAAVASSYAAAQRKLLTPVTVLTAVTFFAICLTWRIEDPP FNSLLFALLAAAGGLQGIYVLVMLVLLILAYRRRWRLTVC GGIMFLACVLVLIVDAVLQLSPLLGA VTVVS MTLLLLAFVLWLSSPGGLGT LG AALLT LAAALALLASLILGT LNLTTM FLLMLLWTLVLLICSSCSCPLSK ILLARLFLYALALLLASALIAGGSILQTNFKSLSSTEFIPNLFCMLLLIVAGILFILAILTEWGSNGNRTYGPVFM CLGGLTMVAGAVWLTVM SNTLLSAWILTAGFLIFLIGFALFGVIRCCRYCCYCLTLESEERPPTPYRNTV
Product Formulation:	Lyophilized as trifluoroacetate salts
Source:	Epstein-Barr virus (strain B95-8)
Number of Peptides:	122
Protein ID:	P13285
Protein Name:	Latent membrane protein 2 (LMP2)
Gene Name:	LMP2
Purity:	Average 70%

Preparation and Storage

Stability and Storage:	Store at -20°C. Stable as supplied until expiry date (EXP) on label.
Preparation:	Warm to room temperature (15 - 25°C) before reconstitution. Add pure dimethyl sulfoxide (DMSO; ~40 µL) and dilute with water to the desired concentration. Final concentration of DMSO must be below 1% (v/v) to avoid toxicity in the biological system. If not used immediately, aliquot and store at -20°C. Protect from light. Avoid repeated freeze-thaw cycles.

Related Products

For a complete list of cytokines or peptide pools, as well as related products available from STEMCELL Technologies, visit www.stemcell.com/cytokines or contact us at techsupport@stemcell.com.

References

- Lee M-A et al. (1999) Genetic evidence that EBNA-1 is needed for efficient, stable latent infection by Epstein-Barr virus. *J Virol* 73(4): 2974–82.
- Longnecker R. (2000) Epstein-Barr virus latency: LMP2, a regulator or means for Epstein-Barr virus persistence. *Adv Cancer Res* 79: 175–200.
- Miller CL et al. (1994) An integral membrane protein (LMP2) blocks reactivation of Epstein-Barr virus from latency following surface immunoglobulin crosslinking. *Proc Natl Acad Sci USA* 91(2): 772–6.

PRODUCTS ARE FOR RESEARCH USE ONLY AND NOT INTENDED FOR HUMAN OR ANIMAL DIAGNOSTIC OR THERAPEUTIC USES UNLESS OTHERWISE STATED.

Copyright © 2024 by STEMCELL Technologies Inc. All rights reserved including graphics and images. STEMCELL Technologies & Design, STEMCELL Shield Design, and Scientists Helping Scientists are trademarks of STEMCELL Technologies Canada Inc. All other trademarks are the property of their respective holders. While STEMCELL has made all reasonable efforts to ensure that the information provided by STEMCELL and its suppliers is correct, it makes no warranties or representations as to the accuracy or completeness of such information.