

Anti-Human IFN-gamma Antibody, Clone 1-D1K

Mouse monoclonal antibody against human IFN-gamma, unconjugated

Catalog #100-1460

250 µg

1 mg/mL

Product Description

This mouse monoclonal antibody (clone 1-D1K) reacts with human interferon-gamma (IFN- γ). IFN- γ is an ~17 kDa protein produced in response to mitogenic and antigenic stimuli by T and NK cells, and in smaller amounts by dendritic cells and macrophages. IFN- γ is controlled by cytokines such as IL-12 and IL-18 secreted in response to infection. IFN- γ functions as a homodimer composed of interconnecting alpha helices which bind to the alpha subunit of the IFN- γ receptor (IFNGR) on target cells, resulting in receptor crosslinking and activation of JAK-STAT signaling. IFN- γ mediates many cellular processes, including upregulation of class I and II antigen presentation pathways, antiviral and antimicrobial mechanisms, antiproliferative and apoptotic effects, immunomodulation, and leukocyte trafficking. The 1-D1K antibody is suitable for the detection of intracellular IFN- γ by flow cytometry.

Target Antigen:	IFN γ
Alternative Names:	IFN- γ , immune interferon, interferon gamma, type II interferon
Gene ID:	3458
Species Reactivity:	Human, Non Human Primate
Host Species:	Mouse
Clonality:	Monoclonal
Clone:	1-D1K
Isotype:	IgG1
Immunogen:	Recombinant human IFN- γ
Conjugate:	Unconjugated

Applications

Verified Applications:	ELISA, ELISpot
Reported Applications:	ELISA, ELISpot, FC, IHC

Abbreviations: CellSep: Cell separation; ChIP: Chromatin immunoprecipitation; FA: Functional assay; FACS: Fluorescence-activated cell sorting; FC: Flow cytometry; ICC: Immunocytochemistry; IF: Immunofluorescence microscopy; IHC: Immunohistochemistry; IHC-P: Immunohistochemistry (paraffin-embedded); IP: Immunoprecipitation; RIA: Radioimmunoassay; WB: Western blotting

Properties

Purification:	The antibody was purified by affinity chromatography.
Stability and Storage:	Product stable at 2 - 8°C when stored undiluted. Do not freeze. Stable until expiry date (EXP) on label.
Directions for Use:	The suggested use of this antibody is: ELISA, 1:500; ELISpot, 1:67. It is recommended that the antibody be titrated for optimal performance for each application.

Related Products

For a complete list of antibodies, including other conjugates, sizes, and clones, as well as related products available from STEMCELL Technologies, visit www.stemcell.com/antibodies, or contact us at techsupport@stemcell.com.

References

- Chen H et al. (2009) Anti-CTLA-4 therapy results in higher CD4+ ICOShi T cell frequency and IFN- γ levels in both nonmalignant and malignant prostate tissues. *Proc Natl Acad Sci U S A* 106(8): 2729–34.
- Kak G et al. (2018). Interferon-gamma (IFN- γ): exploring its implications in infectious diseases. *Biomol Concepts* 9(1): 64–79.
- Nilsen EM et al. (1998) Gluten induces an intestinal cytokine response strongly dominated by interferon gamma in patients with celiac disease. *Gastroenterol* 115(3): 551–63.
- Schroder K et al. (2004) Interferon- γ : an overview of signals, mechanisms and functions. *J Leukoc Biol* 75(2): 163–89.
- Schwinn N et al. (2009) Interferon- γ down-regulates NKG2D ligand expression and impairs the NKG2D-mediated cytolysis of MHC class I-deficient melanoma by natural killer cells. *Int J Cancer* 124(7): 1594–1604.
- Sun P et al. (2003) Protective immunity induced with malaria vaccine, RTS, S, is linked to Plasmodium falciparum circumsporozoite protein-specific CD4+ and CD8+ T cells producing IFN- γ . *J Immunol* 171(12): 6961–67.
- Thiel DJ et al. (2000) Observation of an unexpected third receptor molecule in the crystal structure of human interferon- γ receptor complex. *Structure* 8(9): 927–36.
- Wahlström J et al. (2001) Analysis of intracellular cytokines in CD4+ and CD8+ lung and blood T cells in sarcoidosis. *Am J Respir Crit Care Med* 163(1): 115–21.

PRODUCTS ARE FOR RESEARCH USE ONLY AND NOT INTENDED FOR HUMAN OR ANIMAL DIAGNOSTIC OR THERAPEUTIC USES UNLESS OTHERWISE STATED.

Copyright © 2024 by STEMCELL Technologies Inc. All rights reserved including graphics and images. STEMCELL Technologies & Design, STEMCELL Shield Design, and Scientists Helping Scientists are trademarks of STEMCELL Technologies Canada Inc. All other trademarks are the property of their respective holders. While STEMCELL has made all reasonable efforts to ensure that the information provided by STEMCELL and its suppliers is correct, it makes no warranties or representations as to the accuracy or completeness of such information.