



Release Notes

RoboSep™-S Software Version 1.5.2.

These release notes describe the new functionality and/or fixes to the previous RoboSep™-S Software Version 1.4.0.0.

Software	New Functionality/Fix
RoboSep™-S Software Version 1.5.2.	<ul style="list-style-type: none"> ● User Interface Changes, including: <ul style="list-style-type: none"> ○ User Profile Selection screen upon startup of instrument ○ Separated User Selection and User Management ○ End of Run Report detail view replacing PDF viewer ○ Added Lab Manager Menu ● Implemented failed protocol run recovery feature and downloadable report ● Addition of reagent expiry date verification during reagent vial barcode scanning ● Addition of reagent expiry date to End of Run Report ● Changes to Diagnostic Package including <ul style="list-style-type: none"> ○ Option to include video logs and reports ○ Include Windows events by default ○ Removal of FTP option ● Improved accuracy of required buffer volume displayed at Carousel Load screen ● Improved accuracy of estimated protocol run time displayed at Run Progress screen ● Change to carousel homing logic for improved positional accuracy ● Bug Fixes

Copyright © 2026 by STEMCELL Technologies Inc. All rights reserved including graphics and images. STEMCELL Technologies & Design, STEMCELL Shield Design, Scientists Helping Scientists, and RoboSep are trademarks of STEMCELL Technologies Canada Inc. All other trademarks are the property of their respective holders. While STEMCELL has made all reasonable efforts to ensure that the information provided by STEMCELL and its suppliers is correct, it makes no warranties or representations as to the accuracy or completeness of such information.

UNLESS OTHERWISE STATED, PRODUCTS ARE FOR RESEARCH USE ONLY AND NOT INTENDED FOR HUMAN OR ANIMAL DIAGNOSTIC OR THERAPEUTIC USES. FOR PRODUCT-SPECIFIC COMPLIANCE AND INTENDED USE INFORMATION, REFER TO THE PRODUCT INFORMATION SHEET. GENERAL INFORMATION ON QUALITY AT STEMCELL MAY BE FOUND AT WWW.STEMCELL.COM/COMPLIANCE.

Scientists Helping Scientists™ | WWW.STEMCELL.COM

TOLL FREE PHONE 1 800 667 0322 · PHONE +1 604 877 0713 · INFO@STEMCELL.COM · TECHSUPPORT@STEMCELL.COM

FOR GLOBAL CONTACT DETAILS VISIT WWW.STEMCELL.COM **DOCUMENT #27293, VERSION 1.0.0** MARCH 2026