



Release Notes

RoboSep™-S Software Version 1.3.0.0

This release note describes the new functionality and/or fixes to the previous RoboSep™-S Software Version.

RoboSep™-S Software Version	New Functionality/Fix
1.2.1.1	<ul style="list-style-type: none"> • Pre-set user profiles removed • Updated branding • Barcode scanning updated to include reagent lot number • Bug fixes
1.3.0.0	<ul style="list-style-type: none"> • Customization of reagent names displayed on graphical user interface (GUI): Reagent names can now be customized via the Protocol Editor. • Update to barcode scanning logic, as required for new format of catalog numbers • Update to End of Run report template • Bug fixes

Copyright © 2020 by STEMCELL Technologies Inc. All rights reserved including graphics and images. STEMCELL Technologies & Design, STEMCELL Shield Design, Scientists Helping Scientists, and RoboSep are trademarks of STEMCELL Technologies Canada Inc. All other trademarks are the property of their respective holders. While STEMCELL has made all reasonable efforts to ensure that the information provided by STEMCELL and its suppliers is correct, it makes no warranties or representations as to the accuracy or completeness of such information.

STEMCELL TECHNOLOGIES INC.'S QUALITY MANAGEMENT SYSTEM IS CERTIFIED TO ISO 13485. PRODUCTS ARE FOR RESEARCH USE ONLY AND NOT INTENDED FOR HUMAN OR ANIMAL DIAGNOSTIC OR THERAPEUTIC USES UNLESS OTHERWISE STATED.

Scientists Helping Scientists™ | WWW.STEMCELL.COM

TOLL FREE PHONE 1 800 667 0322 • PHONE +1 604 877 0713 • INFO@STEMCELL.COM • TECHSUPPORT@STEMCELL.COM
 FOR GLOBAL CONTACT DETAILS VISIT WWW.STEMCELL.COM DOCUMENT #27037 VERSION 3.0.0 MAY 2020