



## Release Notes

# STEMvision™ Software Version 2023.08.08.00

These release notes describe the new functionality and/or fixes to the previous STEMvision™ Software Version 2021.05.07.00.

Software	New Functionality/Fix
STEMvision™ Acquisition Version 2.6.3.0	<p>Compatibility with upcoming hardware revisions</p> <p>Bug fixes for improved stability:</p> <ul style="list-style-type: none"> <li>• "Automatic ID" shows the correct values.</li> </ul>
STEMvision™ Analyzer (Human) Version 2.6.7.0	<p>New features:</p> <ul style="list-style-type: none"> <li>• Report type can be selected from the menu</li> <li>• Support the new STEMvision™ Custom Report Generator</li> </ul> <p>Bug fixes for improved stability:</p> <ul style="list-style-type: none"> <li>• CFU-E analysis crash fix</li> <li>• Enable "abort" function with CFR mode</li> <li>• Various cosmetic changes (clarification of labels, etc.)</li> </ul>
STEMvision™ Analyser (Mouse) Version 2.5.0.0	No updates
STEMvision™ Colony Marker Version 2.5.0.0	No updates
STEMvision™ Report Generator Version 1.0.7.0	Initial release

Copyright © 2023 by STEMCELL Technologies Inc. All rights reserved including graphics and images. STEMCELL Technologies & Design, STEMCELL Shield Design, Scientists Helping Scientists, and STEMvision are trademarks of STEMCELL Technologies Canada Inc. All other trademarks are the property of their respective holders. While STEMCELL has made all reasonable efforts to ensure that the information provided by STEMCELL and its suppliers is correct, it makes no warranties or representations as to the accuracy or completeness of such information.

PRODUCTS ARE FOR RESEARCH USE ONLY AND NOT INTENDED FOR HUMAN OR ANIMAL DIAGNOSTIC OR THERAPEUTIC USES UNLESS OTHERWISE STATED. FOR ADDITIONAL INFORMATION ON QUALITY AT STEMCELL, REFER TO [WWW.STEMCELL.COM/COMPLIANCE](http://WWW.STEMCELL.COM/COMPLIANCE).

Scientists Helping Scientists™ | [WWW.STEMCELL.COM](http://WWW.STEMCELL.COM)

TOLL FREE PHONE 1 800 667 0322 · PHONE +1 604 877 0713 · [INFO@STEMCELL.COM](mailto:INFO@STEMCELL.COM) · [TECHSUPPORT@STEMCELL.COM](mailto:TECHSUPPORT@STEMCELL.COM)

FOR GLOBAL CONTACT DETAILS VISIT [WWW.STEMCELL.COM](http://WWW.STEMCELL.COM) DOCUMENT #27251 VERSION 1.0.0 NOV 2023