

Anti-Adherence Rinsing Solution

Rinsing solution for cultureware to prevent cell adhesion

Catalog #07010

100 mL

Product Description

Anti-Adherence Rinsing Solution is a surfactant solution for pre-treating cultureware to reduce surface tension and prevent cell adhesion. It can be used with plates, flasks, pipettes, strainers, and other cultureware.

Anti-Adherence Rinsing Solution is required for use with the following:

- AggreWell[™] plates (e.g. Catalog #34811)
- Organoid growth media, including HepatiCult[™] (Catalog #06030) and PancreaCult[™] (Catalog #06040)

Properties

Stability and Storage: Store at 15 - 25°C. Stable until expiry date (EXP) on label.

Directions for Use

For use with AggreWell[™] plates, refer to the Product Information Sheet (PIS) for AggreWell[™]400 (Document #10000001001) or AggreWell[™]800 (Document #10000001005), available at www.stemcell.com, or contact us to request a copy.

For use with organoid growth media, refer to the appropriate PIS.

For other applications, coat the desired cultureware surface with Anti-Adherence Rinsing Solution, then wash with phosphate-buffered saline (e. g. Catalog #37350) or an appropriate medium.

PRODUCTS ARE FOR RESEARCH USE ONLY AND NOT INTENDED FOR HUMAN OR ANIMAL DIAGNOSTIC OR THERAPEUTIC USES UNLESS OTHERWISE STATED. FOR ADDITIONAL INFORMATION ON QUALITY AT STEMCELL, REFER TO WWW.STEMCELL.COM/COMPLIANCE.

Copyright © 2024 by STEMCELL Technologies Inc. All rights reserved including graphics and images. STEMCELL Technologies & Design, STEMCELL Shield Design, Scientists Helping Scientists, AggreWell, HepatiCult, and PancreaCult are trademarks of STEMCELL Technologies Canada Inc. All other trademarks are the property of their respective holders. While STEMCELL has made all reasonable efforts to ensure that the information provided by STEMCELL and its suppliers is correct, it makes no warranties or representations as to the accuracy or completeness of such information.