

StemSpan™ Lymphoid Progenitor Expansion Supplement

StemSpan™ Lymphoid Differentiation Coating Material

Supplement and coating material for generation of lymphoid progenitor cells from human CD34⁺ hematopoietic progenitor cells



Scientists Helping Scientists™ | WWW.STEMCELL.COM

TOLL FREE PHONE 1 800 667 0322 • PHONE +1 604 877 0713

INFO@STEMCELL.COM • TECHSUPPORT@STEMCELL.COM

FOR GLOBAL CONTACT DETAILS VISIT OUR WEBSITE

Product Description

StemSpan™ Lymphoid Progenitor Expansion Supplement contains recombinant human cytokines and other additives formulated to selectively promote the expansion and differentiation of CD34⁺ cells to lymphoid progenitor cells when used in combination with StemSpan™ SFEM II (Catalog #09605), and on plates coated with StemSpan™ Lymphoid Differentiation Coating Material.

Product Information

PRODUCT NAME	CATALOG #	SIZE	STORAGE	SHELF LIFE
StemSpan™ Lymphoid Progenitor Expansion Supplement (10X)*	09915	5 mL	Store at -20°C.	Stable for 1.5 years from date of manufacture (MFG) on label.
StemSpan™ Lymphoid Differentiation Coating Material (100X)*	09925	250 µL	Store at -20°C.	Stable for 1.5 years from date of manufacture (MFG) on label.

*This product contains material derived from human plasma. Donors have been tested and found negative for HIV-1 and -2, hepatitis B, and hepatitis C prior to donation. However, this product should be considered potentially infectious and treated in accordance with universal handling precautions.

Directions for Use

For preparation of StemSpan™ Lymphoid Progenitor Expansion Medium and StemSpan™ Lymphoid Differentiation Coating Material, and for complete directions for use, refer to the Technical Manual for T cells (Document #10000007541) or NK cells (Document #10000007537), available at www.stemcell.com or contact us to request a copy.

PRODUCTS ARE FOR RESEARCH USE ONLY AND NOT INTENDED FOR HUMAN OR ANIMAL DIAGNOSTIC OR THERAPEUTIC USES UNLESS OTHERWISE STATED. FOR ADDITIONAL INFORMATION ON QUALITY AT STEMCELL, REFER TO WWW.STEMCELL.COM/COMPLIANCE.

Copyright © 2020 by STEMCELL Technologies Inc. All rights reserved including graphics and images. STEMCELL Technologies & Design, STEMCELL Shield Design, Scientists Helping Scientists, and StemSpan are trademarks of STEMCELL Technologies Canada Inc. All other trademarks are the property of their respective holders. While STEMCELL has made all reasonable efforts to ensure that the information provided by STEMCELL and its suppliers is correct, it makes no warranties or representations as to the accuracy or completeness of such information.