

StemSpan™ T Cell Progenitor Maturation Supplement (10X)

Supplement for the maturation of lymphoid progenitors to T cells

Catalog #09930

12.5 mL



Scientists Helping Scientists™ | WWW.STEMCELL.COM

TOLL FREE PHONE 1 800 667 0322 • PHONE +1 604 877 0713

INFO@STEMCELL.COM • TECHSUPPORT@STEMCELL.COM

FOR GLOBAL CONTACT DETAILS VISIT OUR WEBSITE

Product Description

StemSpan™ T Cell Progenitor Maturation Supplement (10X) enables further maturation of lymphoid progenitor cells to double-positive (DP) T cells. It also supports maturation of DP T cells to CD8 SP T cells with additional stimuli. StemSpan™ T Cell Progenitor Maturation Supplement (10X) is intended for use in combination with StemSpan™ SFEM II (Catalog #09605), and on plates coated with StemSpan™ Lymphoid Differentiation Coating Material (100X; Catalog #09925).

Properties

Storage: Store at -20°C.

Shelf Life: Stable for 12 months from date of manufacture (MFG) on label.

Directions For Use

For complete instructions, refer to the Technical Manual for STEMdiff™ or StemSpan™ T Cell Kits (Document #10000007541), available at www.stemcell.com or contact us to request a copy.

PRODUCTS ARE FOR RESEARCH USE ONLY AND NOT INTENDED FOR HUMAN OR ANIMAL DIAGNOSTIC OR THERAPEUTIC USES UNLESS OTHERWISE STATED. FOR ADDITIONAL INFORMATION ON QUALITY AT STEMCELL, REFER TO WWW.STEMCELL.COM/COMPLIANCE.

Copyright ©2020 by STEMCELL Technologies Inc. All rights reserved including graphics and images. STEMCELL Technologies & Design, STEMCELL Shield Design, Scientists Helping Scientists, STEMdiff, and StemSpan are trademarks of STEMCELL Technologies Canada Inc. All other trademarks are the property of their respective holders. While STEMCELL has made all reasonable efforts to ensure that the information provided by STEMCELL and its suppliers is correct, it makes no warranties or representations as to the accuracy or completeness of such information.