

Falcon® Serological Pipettes, 1 mL

Sterile polystyrene pipettes for cell culture

Catalog # 38001

1000 Pipettes



Scientists Helping Scientists™ | WWW.STEMCELL.COM

TOLL FREE PHONE 1 800 667 0322 • PHONE +1 604 877 0713

INFO@STEMCELL.COM • TECHSUPPORT@STEMCELL.COM

FOR GLOBAL CONTACT DETAILS VISIT OUR WEBSITE

Product Description

Falcon® Serological Pipettes, 1 mL can be used for tissue culture or with any non-toxic solution that requires aliquoting or aspirating/dispensing in a sterile environment. These pipettes have been tested for dispensing and aspirating mTeSR™1 (Catalog #85850) and TeSR™-AOF (Catalog #100-0401) media. They have also been used in cell separation with EasySep™ Magnet (Catalog #18000) and EasyRights™ EasySep™ Magnet (Catalog #18103).

Catalog #38001, 1000 Pipettes; 100 Pipettes/box, individually wrapped (Corning Catalog #357521)

APPLICATIONS

- Cell separation
- Tissue culture
- Aliquoting or aspirating/dispensing in a sterile environment

Properties

Storage: Store at 15 - 25°C.

Shelf Life: Stable until expiry date on manufacturer label.

Features

MATERIALS

- Pipette: clear and transparent polystyrene
- Plug: polyester fiber
- Color-coded yellow
- Non-pyrogenic, low-endotoxin

CAPACITY

- 1 mL
- Increments of 0.01 mL with reverse graduations
- Overfill is 0.04 mL
- Reliable and accurate to $\pm 2\%$ at full volume

STERILITY

- Irradiation

Directions for Use

Use as directed in protocol of choice.

Notes and Tips

For chemical resistance or thermoplastic property information, contact us at techsupport@stemcell.com.

Related Products

For a complete list of related products available from STEMCELL Technologies, visit www.stemcell.com, or contact us at techsupport@stemcell.com.

PRODUCTS ARE FOR RESEARCH USE ONLY AND NOT INTENDED FOR HUMAN OR ANIMAL DIAGNOSTIC OR THERAPEUTIC USES UNLESS OTHERWISE STATED.

Copyright © 2023 by STEMCELL Technologies Inc. All rights reserved including graphics and images. STEMCELL Technologies & Design, STEMCELL Shield Design, Scientists Helping Scientists, EasyEights, and EasySep are trademarks of STEMCELL Technologies Canada Inc. Falcon is a registered trademark of Corning Incorporated. mTeSR and TeSR are trademarks of WARF. All other trademarks are the property of their respective holders. While STEMCELL has made all reasonable efforts to ensure that the information provided by STEMCELL and its suppliers is correct, it makes no warranties or representations as to the accuracy or completeness of such information.