

Round-Bottom Polystyrene Tubes with Caps, 5 mL



Sterile polystyrene round-bottom tubes for cell separation and cell culture

Scientists Helping Scientists™ | WWW.STEMCELL.COM

TOLL FREE PHONE 1 800 667 0322 • PHONE +1 604 877 0713

INFO@STEMCELL.COM • TECHSUPPORT@STEMCELL.COM

FOR GLOBAL CONTACT DETAILS VISIT OUR WEBSITE

Catalog # 100-0088	25 Tubes
38007	500 Tubes
38025	1000 Tubes

Product Description

Falcon® Round-Bottom Polystyrene Tubes, 5 mL are sterile, high-quality tubes ideal for a wide variety of sample preparation, sample storage, and sample centrifugation tasks. They are required for EasySep™ cell isolation when using EasySep™ Magnet (Catalog #18000) or EasyEights™ EasySep™ Magnet (5 mL side; Catalog #18103), and can be used for flow cytometry and fluorescence-activated cell sorting (FACS). Tubes are RNase-/DNase-free.

Catalog #100-0088, 5 mL - 25 Tubes with caps; 25 Tubes/pack (Corning Catalog #352058)

Catalog #38007, 5 mL - 500 Tubes with caps; 25 Tubes/pack (Corning Catalog #352058)

Catalog #38025, 5 mL - 1000 Tubes with caps; 125 Tubes/pack (Corning Catalog #352054)

APPLICATIONS

- Cell separation
- Cell culture
- Sample storage
- Liquid handling

Properties

Storage: Store at 15 - 25°C.

Shelf Life: For Catalog #100-0088: Stable until expiry date on STEMCELL label.

For Catalog #38007 and #38025: Stable until expiry date on manufacturer label.

Features

MATERIALS

- Tube: clear and transparent polystyrene
- Cap: low-density polyethylene
- Dual-position snap cap with a tight seal or a loose position for sterile aerobic culture
- Tube and cap are biologically inert, RNase-/DNase-free, and non-pyrogenic

DIMENSIONS

- 12 x 75 mm (outer diameter x height)

NOMINAL VOLUME

- 5 mL
- No graduations

CENTRIFUGAL FORCE

- Up to 1400 x g

STERILITY

- Gamma irradiation

Handling / Directions For Use

Use as directed in protocol of choice.

Notes and Tips

For chemical resistance or thermoplastic property information, contact us at techsupport@stemcell.com.

Related Products

For a complete list of related products available from STEMCELL Technologies, visit www.stemcell.com or contact us at techsupport@stemcell.com.

STEMCELL TECHNOLOGIES INC.'S QUALITY MANAGEMENT SYSTEM IS CERTIFIED TO ISO 13485. PRODUCTS ARE FOR RESEARCH USE ONLY AND NOT INTENDED FOR HUMAN OR ANIMAL DIAGNOSTIC OR THERAPEUTIC USES UNLESS OTHERWISE STATED.

Copyright © 2019 by STEMCELL Technologies Inc. All rights reserved including graphics and images. STEMCELL Technologies & Design, STEMCELL Shield Design, Scientists Helping Scientists, EasyEights, and EasySep are trademarks of STEMCELL Technologies Canada Inc. Falcon is a registered trademark of Corning Incorporated. All other trademarks are the property of their respective holders. While STEMCELL has made all reasonable efforts to ensure that the information provided by STEMCELL and its suppliers is correct, it makes no warranties or representations as to the accuracy or completeness of such information.