

Falcon® Conical Tubes, 50 mL



Scientists Helping Scientists™ | WWW.STEMCELL.COM

TOLL FREE PHONE 1 800 667 0322 • PHONE +1 604 877 0713

INFO@STEMCELL.COM • TECHSUPPORT@STEMCELL.COM

FOR GLOBAL CONTACT DETAILS VISIT OUR WEBSITE

Sterile polypropylene conical tubes

Catalog # 100-0090	25 Tubes
38010	500 Tubes

Product Description

Falcon® Conical Tubes, 50 mL are recommended for use with EasySep™ cell separation when using the Easy50 EasySep™ Magnet (Catalog #18002). These tubes are also recommended for cell centrifugation and are suitable for a variety of other applications.

Catalog #100-0090, 50 mL - 25 Tubes; 25 Tubes/pack (Corning Catalog #352070)

Catalog #38010, 50 mL - 500 Tubes; 25 Tubes/pack (Corning Catalog #352070)

APPLICATIONS

- Cell centrifugation: separation by density gradients or pelleting
- Dilution of samples or solutions
- Sample storage
- Molecular biology: DNA purification

Properties

Storage: Store at 15 - 25°C.

Shelf Life: For Catalog #100-0090: Stable until expiry date (EXP) on STEMCELL label.

For Catalog #38010: Stable until expiry date on manufacturer label.

Features

MATERIALS

- Tube: translucent polypropylene
- Cap: high-density polyethylene
- Tube and cap are biologically inert, RNase/DNase-free, and non-pyrogenic

DIMENSIONS

- 30 x 115 mm (outer diameter x height)

NOMINAL VOLUME

- 50 mL
- Graduations every 5 mL

CENTRIFUGAL FORCE

- Up to 16,000 x g

STERILITY

- Gamma irradiation

Handling / Directions For Use

Use as directed in protocol of choice.

Notes and Tips

For information on chemical resistance or thermoplastic property, contact us at techsupport@stemcell.com.

Related Products

For a complete list of related products available from STEMCELL Technologies, visit www.stemcell.com or contact us at techsupport@stemcell.com.

Copyright © 2019 by STEMCELL Technologies Inc. All rights reserved including graphics and images. STEMCELL Technologies & Design, STEMCELL Shield Design, Scientists Helping Scientists, and EasySep are trademarks of STEMCELL Technologies Canada Inc. Falcon is a registered trademark of Corning Incorporated. All other trademarks are the property of their respective holders. While STEMCELL has made all reasonable efforts to ensure that the information provided by STEMCELL and its suppliers is correct, it makes no warranties or representations as to the accuracy or completeness of such information.

STEMCELL TECHNOLOGIES INC.'S QUALITY MANAGEMENT SYSTEM IS CERTIFIED TO ISO 13485. PRODUCTS ARE FOR RESEARCH USE ONLY AND NOT INTENDED FOR HUMAN OR ANIMAL DIAGNOSTIC OR THERAPEUTIC USES UNLESS OTHERWISE STATED.