

# Corning® 96-Well High-Binding Flat-Bottom Microplate

Clear polystyrene flat-bottom, high-binding microplate without lid

Catalog # 38019      6 Plates/Pack  
Catalog # 200-0480      100 Plates/Case



Scientists Helping Scientists™ | [WWW.STEMCELL.COM](http://WWW.STEMCELL.COM)

TOLL FREE PHONE 1 800 667 0322 • PHONE +1 604 877 0713  
[INFO@STEMCELL.COM](mailto:INFO@STEMCELL.COM) • [TECHSUPPORT@STEMCELL.COM](mailto:TECHSUPPORT@STEMCELL.COM)  
FOR GLOBAL CONTACT DETAILS VISIT OUR WEBSITE

## Product Description

The Corning® 96-Well High-Binding Flat-Bottom Microplate is recommended for use in radioimmunoassays (RIAs) and enzyme immunoassays (EIAs).

Catalog #38019, 6 Plates/pack, individually wrapped (Corning Catalog #3590)

Catalog #200-0480, 100 Plates/case, individually wrapped (Corning Catalog #3590)

### APPLICATIONS

- Immunoassays

## Properties

**Storage:** Store at 15 - 25°C.

**Shelf Life:** Stable until expiry date on outer package label.

## Features

### MATERIALS

- Plate: clear and transparent polystyrene
- Binding capacity is approximately 500 ng IgG/cm<sup>2</sup>

### DIMENSIONS

- Conforms to standard microplate footprint and dimensions

### CAPACITY

- 0.36 mL per well
- Recommended working volume of 0.075 - 0.2 mL

## Handling / Directions For Use

Use as directed in immunoassay of choice.

## Notes and Tips

For information on chemical resistance or thermoplastic property, contact us at [techsupport@stemcell.com](mailto:techsupport@stemcell.com).

## Related Products

For a complete list of related products available from STEMCELL Technologies, visit [www.stemcell.com](http://www.stemcell.com) or contact us at [techsupport@stemcell.com](mailto:techsupport@stemcell.com).

PRODUCTS ARE FOR RESEARCH USE ONLY AND NOT INTENDED FOR HUMAN OR ANIMAL DIAGNOSTIC OR THERAPEUTIC USES UNLESS OTHERWISE STATED.

Copyright © 2021 by STEMCELL Technologies Inc. All rights reserved including graphics and images. STEMCELL Technologies & Design, STEMCELL Shield Design, and Scientists Helping Scientists are trademarks of STEMCELL Technologies Canada Inc. Corning is a registered trademark of Corning Incorporated. All other trademarks are the property of their respective holders. While STEMCELL has made all reasonable efforts to ensure that the information provided by STEMCELL and its suppliers is correct, it makes no warranties or representations as to the accuracy or completeness of such information.