

Corning® 245 mm Square Dish, Non-Treated



Scientists Helping Scientists™ | WWW.STEMCELL.COM

TOLL FREE PHONE 1 800 667 0322 • PHONE +1 604 877 0713

INFO@STEMCELL.COM • TECHSUPPORT@STEMCELL.COM

FOR GLOBAL CONTACT DETAILS VISIT OUR WEBSITE

Square dish with handles for use as a moisture chamber

Catalog #38020

16 Dishes

Product Description

The Corning® 245 mm Square Dish, Non-Treated is optically clear, sterile, and stackable. It can be used in any non-adherent application; for example, it can be used as a moisture chamber for incubating plates with medium such as MyeloCult™ H5100 (Catalog #05150).

Catalog #38020, 16 Dishes; 4 Dishes/bag (Corning Catalog #431111)

APPLICATIONS

- Cell culture
- Sample storage
- Screening libraries or isolating clones
- Bacterial culture
- Suitable for agar diffusion assay

Properties

Storage: Store at 15 - 25°C.

Shelf Life: Stable until expiry date on manufacturer label.

Features

MATERIALS

- Plate and lid: clear and transparent polystyrene
- Non-pyrogenic

DIMENSIONS

- 245 x 245 mm

CAPACITY

- 150 mL

STERILITY

- Gamma irradiation

Directions for Use

Use as directed in protocol of choice.

Notes and Tips

For chemical resistance or thermoplastic property information, contact us at techsupport@stemcell.com.

Related Products

For a complete list of related products available from STEMCELL Technologies, visit www.stemcell.com, or contact us at techsupport@stemcell.com.

PRODUCTS ARE FOR RESEARCH USE ONLY AND NOT INTENDED FOR HUMAN OR ANIMAL DIAGNOSTIC OR THERAPEUTIC USES UNLESS OTHERWISE STATED.

Copyright © 2023 by STEMCELL Technologies Inc. All rights reserved including graphics and images. STEMCELL Technologies & Design, STEMCELL Shield Design, Scientists Helping Scientists, and MyeloCult are trademarks of STEMCELL Technologies Canada Inc. Corning is a registered trademark of Corning Incorporated. All other trademarks are the property of their respective holders. While STEMCELL has made all reasonable efforts to ensure that the information provided by STEMCELL and its suppliers is correct, it makes no warranties or representations as to the accuracy or completeness of such information.